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TABLE OF CONTENTS

Important notes	20	5. Product operation	26
Safety warning	20	6. Troubleshooting and maintenance	35
1. Matters needing attention	21	Appendix: Comparison table of tube sizes ...	37
2. Unboxing	22		
3. Product description	22		
4. Product installation	24		

Important notes

- Please read this operation manual carefully before using the product.
- The manufacturer reserves the right to change the product (design or specification) without prior notice.

Safety warning

- Before any cleaning or maintenance work, be sure to cut off the power supply.
- The tube may have cracks due to wear and tear, causing liquid to overflow from the tube, which may cause harm to the human body and equipment, so check it frequently and replace the tube in time!
- Please connect the power cord directly to the wall outlet and avoid using extension cords
- If the power cord is damaged, unplug it.
- If the following situations occur, turn off the power and unplug the power plug:
 1. Fluid is spilled on the machine
 2. You think this machine needs maintenance or repair
- The power supply must have a reliable grounding.
- The foot switch or other external control plug must be installed and unloaded when the power is off to prevent the external control interface from being burned.

Appendix: Comparison table of tube sizes

Tube	13#	14#	19#	16#	25#	17#	18#
Tube section (1:1)							
Wall thickness (mm)	1,6						
Inner diameter (mm)	0,8	1,6	2,4	3,1	4,8	6,4	7,9
Pressure (MPa)	Continuous	0,17				0,14	0,07
	Interval	0,27				0,14	0,1


Tube	15#	24#	35#	36#
Tube section (1:1)				
Wall thickness (mm)	2,4			
Inner diameter (mm)	4,8	6,4	7,9	9,6
Pressure (MPa)	Continuous	0,17		0,14
	Interval	0,27		0,24

6.2. Product maintenance

Warning: Before attempting any maintenance, be sure to cut off the power to the pump.


Replace the fuse

1. Place the power switch in the "OFF" position ("|" ON, "O" OFF)
2. Disconnect the AC power input cord from the outlet.
3. Take out the built-in fuse from the power socket of the pump.
4. Replace the original fuse.

Replace the pump head coupling

1. Place the power switch in the "OFF" position ("|" ON, "O" OFF)
2. Disconnect the AC power input cord from the outlet.
3. Remove the tube, the pump head and take out the coupling.
4. Install a new coupling, install the pump head and pump tube.

Basic maintenance and cleaning

Basic maintenance

1. Open the pump head when it is not working to avoid tube deformation caused by prolonged extrusion.
2. Keep the pump head rollers clean and dry to prevent surface damage and reduce tube wear; if there is splashing liquid, please wipe it dry as soon as possible.
3. Check the wear of the tube regularly and replace it in time to prevent leakage.
4. The pump head rollers do not need to add lubricating oil; an improper operation may cause the tube to shift or damage.
5. Not used to deal with chemical substances incompatible with the pump head or tube.
6. The pump head is not resistant to organic solvents and strong corrosive liquids. Please deal with it in time if there is effusion.
7. Built-in fuse, pump head shaft and other replaceable accessories, need to be installed under the guidance of professionals.
8. It is recommended that the working environment temperature be between 0-40°C.

Cleaning

Warning: Before attempting any cleaning, be sure to cut off the power to the pump.
When there are stubborn stains on the pump housing, please use a mild detergent to scrub the surface. Do not immerse the pump in liquid or use too much liquid to clean it.

1. Matters needing attention

Safety:

- The staff responsible for the installation or maintenance of this equipment should have the experience and ability to carry out related work.
- This product is not applicable to the ATEX explosion-proof directive and cannot be used in flammable and explosive environments.
- When pumping dangerous liquids, please follow safety precautions.
- Please determine whether you need to wear personal protective equipment when operating the pump in accordance with the nature of the transfer fluid and industry specifications.
- Non-professionals should not install this pump with other equipment to reduce safety risks.
- For hazardous fluids, a dedicated operation process must be specified to prevent personal injury.
- The power plug can disconnect the power supply and drive in an emergency. Do not place the pump in a workplace where it is difficult to cut off the power supply, which is conducive to emergency stop operations.

Tube:

- In the event of a tube failure, ensure that the fluid in the pump tube of the pump head can be discharged to a suitable container or drain.
- A ruptured tube may cause fluid to splash. Please take appropriate protective measures.
- When disassembling the tube, it is necessary to drain the medium and cut off the power supply to ensure that the pipeline is pressure-free.
- Ensure that the chemicals to be handled are compatible with the pump head, tubes and Accessories.

Rollers:

- Do not touch the rollers while the pump is running.
- Keep the rollers clean and dry to reduce tube wear.
- Do not lubricate the pump head rollers by yourself. Improper operation may cause the tube to run out or the pump head shell to corrode.

Drive:

- There are no user-serviceable parts in the pump.
- The power socket on the back of the driver is equipped with a user-replaceable built-in fuse. Only products of the same category can be used to replace the fuse.
- The surface of the driver and the pump head are not resistant to organic solvents and strong corrosive fluids. If the liquid is splashed or accumulated, please remove and clean it in time.

2. Unboxing

2.1. Unpacking inspection


Confirm that the pump is packaged in good condition. Please check the packing list, when unpacking, check the product model and the number of accessories, and check whether the parts are damaged during transportation. If you have any questions, please contact us immediately.

2.2. Product storage

This product can be stored for a long time, but before putting it into operation, please confirm that the drive, pump head, tubes and other accessories can be used normally.

3. Product description

3.1. Principle of Peristaltic Pump Operation

The peristaltic pump uses the rotor to alternately squeeze and release the tube to transfer fluid, just like squeezing a tube full of fluid with a finger. As the finger slides forward, when a negative pressure is formed in the tube, the liquid flows with it.

3.2. Product features

1. LCD 3.5-inch color touch screen.
2. Flow calibration: With flow calibration function, it is suitable for transferring liquids with higher precision.
3. Compatible with multiple pump heads: YZ and TX series pump heads can be installed, and the pump heads can be cascaded.
4. Multiple working modes: transmission mode, timing mode, quantitative mode, timing quantitative mode.
5. Data storage: 9 sets of data can be stored for users.
6. Support external control: support RS485/MODBUS communication protocol, analog control, foot switch control.

6. Troubleshooting and maintenance

6.1. Troubleshooting

No response at boot	> If a circuit protection device is installed, confirm that the circuit has not tripped > Confirm that the power plug is inserted into a working socket > Check whether the power cord plug is firmly inserted into the pump socket > Check whether the fuse at the power interface is blown
The fan and display screen are normal, but pump cannot be started	> Check if the device is in external control mode > Check if the keys are working
The pump is turned on and the pump head cannot run	> Check if the coupling is damaged > After cutting off the power, manually check whether the pump head is rotating normally
Low or no flow when the pump is running	> Check whether the material supply is normal > Check if the pipe is entangled or blocked > Check that all valves are open > Check if the tube is in the middle of the roller > Check whether the tube is cracked or damaged > Check the running direction > Check whether the pump head rollers can rotate flexibly
Pump cannot be controlled in external control mode	> Check the upper right corner of the LCD, the external control icon  is displayed as  > Check whether the external control settings are correctly connected > Check if the signal source is normal

5.5. External control operation

Foot switch operating instructions

The pedal can be operated by foot or hand. This device only controls the start and stop status of the pump; the running speed and direction are set through the operation panel. The foot switch is connected to the DB-15 external control interface of the peristaltic pump. Under external control, the Start/Stop key and the Full Speed key are invalid. The HOME key, direction key, and touch screen settings are used normally.

Operation setting:

Main interface → Click  → Click **System** → Click **ExInput**:

Ctrl → Analog → Ok

Item → Start Stop → Ok

Signal → Level/Pulse → Ok

After setting, exit to the main interface

The main interface  is displayed as , the setting is successful, and the external control mode is entered.

*Liquid volume transfer

[Level signal] press the foot switch for a long time, release it to continue running, press it again to stop.

[Pulse signal] press to start, press again to stop.


Accessories:

The following accessories are included with the pump

Accessories	Image	Features
Silicone tube #18		Inner diameter: 7,9 mm Wall thickness: 1,6 mm
Filling nozzle		Precise filling, # 18
Foot switch		With 15-pin RS485 connector Control of peristaltic pump start and stop
Check valve		Prevents backflow of transferred liquid, 10 mm
Straight connector		Connects 2 tubes of the system, 7.9 mm

3.3. Product structure

PUMP

A: Drive

B: Pump head

C: Power switch

D: External interface

E: Fuse

F: Power outlet

PUMP HEAD/TUBE/REFERENCE FLOW RATE

Max. Speed		13# (0.8*1.6)	14# (1.6*1.6)	19# (2.4*1.6)	16# (3.1*1.6)	25# (4.8*1.6)	17# (6.4*1.6)	18# (7.9*1.6)
 YZ1515x	600 rpm	42	162	306	492	1020	1740	2280
 TX315	400 rpm	28	106	227	399	830	1400	1930
Max. Speed		15# (4.8*2.4)			24# (6.4*2.4)			
 YZ2515x	600 rpm	1020			1740			
 TX325	400 rpm	760			1230			

Flow rate until: (ml/min)

Tube model: (DI*Wall thickness) mm

Note: The flow rate values were obtained with water at normal temperature and pressure and are for reference only. Obtain the flow rate value under your actual conditions.

3.4. Technical parameters

Max Speed	600 rpm (reversible)
Speed resolution	0.1 rpm
Max flow	2280 ml/min
Motor type	Stepper motor
Screen	LCD 3.5-inch color touch screen
Power supply	220VAC, 50/60Hz
Suction angle	10°-720°(0° is no return suction)
Suction speed	10-300 rpm
Power	<50W
External control interface	DB-15
External control method	Start control/direction control/speed control (0-5V, 0-10V, 4-20mA optional) RS485 serial communication
Ambient temperature	0°C-40°C
Protection level	IP31 (indoor use, avoid long-term exposure to UV rays)

4. Product installation

4.1. Pump head/Pump tube installation

Be sure to disconnect the pump from the main power supply.

PUMP HEAD INSTALLATION DIAGRAM

YZ1515x (YZ2515x)

TX315 (TX325)


Turn clockwise to buckle on the connecting plate


Mode selection instructions

CT series peristaltic pump has four modes to choose from: transmission mode, timing mode, quantitative mode, timing quantitative mode


* Main interface → Click → Click **Mode** → Select **Transfer/Timing/Ration/Filling** (Or click **Exit** to return to the previous menu)


Transmission mode (Transfer)

Quantitative mode (Ration)

Timing mode (Timing)

Timing quantitative mode (Filling)

Data operation instructions

CT series peristaltic pump can realize user data storage function: save 9 groups of filling data.

Click the data shortcut shortcut on the main interface to enter the data operation interface, or Main interface → Click → Click **Data** enter the data operation interface.

Num	Head	Tube	Flow	VOL	Time
1	YZ1515x	25#	60ml/m	50ml	50s
2					
3					
4					
5					
6					
7					
8					
9					
Del Save Use Exit					

- Num: 1 is the first set of saved data
- Head: : the pump head model currently in use according to the settings
- Tube: the model of the pump tubing that is currently in use
- Flow: display of liquid flow rate
- VOL: display of the volume of fluid transferred
- Time: transmission time display

Del: Select a group of data → Click **Del** → **Yes** → completed data deletion

Save: Select a location → Click **Save** → **Yes** → Completed the storage of the current data (this group of data is stored in the selected location; if there are data in this location, the original data will be overwritten).


Use: Select a set of data → Click **Use** → Enter the main interface automatically

Exit: Click **Exit** → Automatically return to the main interface

Calibration Guide

This peristaltic pump has a flow calibration function, which enables high-precision pumping. When a new tube model is used, external conditions change, or other conditions cause a deviation in the flow value, the flow needs to be calibrated to meet the user's requirements for flow accuracy.

1. Main Interface → Click  → Click **Calibration** → Keyboard input calibration value → Click **Next** to enter the calibration interface (or click **Exit** to return to the previous menu).


- **Tube:** Automatically display the tube model currently in use
- **Flowrate:** Automatically display the flow rate generated when filling volume and filling time (ml/min)
- **Filling volume:** The initial default value is the filling volume displayed on the main interface (Can be input and adjusted according to actual needs)


Note: If the amount of liquid to be calibrated is very large and the required time is long, the set value can be reduced according to the filling situation.

2. Calibration interface:

1. Press  to start transmission and display the remaining time. When the remaining time is 0, it will automatically stop and return to the calibration interface (If you need to stop, click **Cancel**, return to the calibration interface, the value of this transfer fluid is invalid)
2. Repeat the calibration 3 times → Enter the actual measured value in turn on the keypad → Click **Next** → Click **Ok**, Save the results to complete the calibration (Or click **Back**, return to the calibration interface, this time the result is invalid)


Countdown interface


Actual filling volume interface

PUMP TUBE INSTALLATION DIAGRAM

YZ1515x (YZ2515x)

TX315 (TX325)

4.2. Installation suggestions and precautions

SUGGESTIONS

- > Application accessories such as foot switch, countersunk head, check valve, filling nozzle, connector, etc. can be selected according to actual conditions.
- > For the size and selection of the tube, please refer to 3.3 Pump head/Tube/Reference flow.
- > For pump head models, please refer to 3.3 Pump Head/Tube/Reference flow.

1. Before cleaning, maintaining and installing the equipment, be sure to disconnect the control power supply;
2. The driver should be placed on a flat and rigid surface;
3. The ambient temperature of the pump should not exceed 104°F (40°C), and air circulation should be ensured to ensure the heat dissipation of the pump;
4. The start-stop key (shortcut key) on the operation panel can quickly change the direction and control the start-stop, but it is recommended to install an emergency stop device on the main circuit of the power supply to ensure higher safety;
5. Make sure that the inner wall of the tube is clean and free of foreign matter before use. The shorter the pipeline, the better, and the suction and lift should not be too long;
6. Determine the running direction of the pump (forward and reverse) according to the specific location of the fluid placement and supporting machinery on site, which is conducive to later operation;
7. In order to meet the requirements of flow and flow rate, a peristaltic pump tube with matching diameter is required;
8. The pump itself has self-priming characteristics, which can effectively prevent the backflow of liquid. Generally, there is no need to install valves at the outlet and inlet of the tube. You can also install a one-way valve in the pipeline according to actual needs to avoid fluid leakage when the pump head and tube fail.

PRECAUTIONS

- The diameter of the pipeline at the inlet is not less than the inner diameter of the pump tube, and a delivery pipe with a diameter $\geq$ the inner diameter of the pump tube should be selected.
- When transferring viscous liquids, it is necessary to maintain a low speed operation to improve the filling efficiency. It is recommended to connect a flexible tube no less than 1 meter between the inlet and the outlet to reduce the pulse and reduce the pulse loss.
- Try to put the pump at the same level or a lower level of the liquid to be transferred to improve the transfer efficiency of the pump.
- To replace a new tube or liquid, re-calibrate the liquid volume to ensure the accuracy of liquid transmission.
- When the peristaltic pump is running, all valves in the pipeline must be opened normally.
- Control wires and power wires are not allowed to have sharp bends, and it is not recommended to bundle them together.
- This product cannot be used for the transmission of any chemical substances incompatible with the pump head and tube.


5. Product operation

5.1. Line connection


Power connection:

- 220VAC $\pm$ 10%, 50/60Hz.
- Ensure that all power supplies are matched to equipment power and are well grounded.
- The position of the pump should ensure that it is convenient to disconnect the power supply when operating the equipment

Power supply wiring diagram

External control wiring diagram

5.2. Power-on

Power-on inspection

- Check whether the pump pipe has been installed correctly, and whether the tube inlet pipe and outlet pipe have been correctly connected.
- Check whether it is connected to a matching power supply.
- Check whether the peristaltic pump has been installed according to “4.2 Installation suggestions and precautions”.

ExOutput Select External control output parameters

Main interface $\rightarrow$ Click $\rightarrow$ **system** $\rightarrow$ Click **ExOutput** $\rightarrow$ External control output option interface.

Item $\rightarrow$ Start Stop/Direction/Speed $\rightarrow$ Click **Enter**, save and exit; Click **Cancel**, exit without saving.

Signal $\rightarrow$ Level/Pulse $\rightarrow$ Click **Enter** save and exit; Click **Cancel**, exit without saving.

Speed $\rightarrow$ 0-5V/0-10V/4-20mA/0-10kHz $\rightarrow$ Click **Enter**, save and exit; Click **Cancel**, exit without saving.

Uart Select the pump communication parameters

Main interface $\rightarrow$ Click $\rightarrow$ **system** $\rightarrow$ Click **Uart** $\rightarrow$ Enter the address of the pump on the keypad (default 1) / Baud rate (1200 /2400 /4800/ 9600) $\rightarrow$ Click **Enter**, save and exit; Click **Cancel**, exit without saving.

- CT pump also has 485 serial communication bus interfaces to be connected to the host computer (computer, PLC, SCM).
- PC can simultaneously connect up to 30 devices with 485 serial communication bus interface functions; when using multiple devices to communicate with the host computer, it must know the machine number for each device, this device ID is the only number that should be connected together, and each machine number of the devices is not the same.
- The factory default value for each pump is “1”; it could be reset by remote controlling through the host computer or by manually changing from the setting menu.
- Change it with the following steps: Set $\rightarrow$ System $\rightarrow$ Uart $\rightarrow$ Address (Default 1) $\rightarrow$ small keypad input.

Suck Back Select Suck-back parameters

Main interface $\rightarrow$ Click $\rightarrow$ **system** $\rightarrow$ Click **Suck Back** $\rightarrow$ Small keyboard input suck back angle (10-720deg)

$\rightarrow$ Small keyboard input suck back speed (10-300 rpm)

$\rightarrow$ Click **Enter**, save and exit; Click **Cancel**, exit without saving.

Note: In filling mode, pause time should be greater than the suck back time.

Language Choose Chinese/English

Main interface $\rightarrow$ Click $\rightarrow$ **system** $\rightarrow$ Click **Language** $\rightarrow$ Chinese/English; Click **Enter**, save and exit; Click **Cancel**, exit without saving.

*After completing the required system settings, return to the previous menu or exit to the main interface to perform other operations.

System setting instructions

*Main interface → Click  → Click System → Enter the system setting options interface (as follows):

Head	Tube	ExInput	ExOutput	Uart	Suck Back	Languaje	Exit
------	------	---------	----------	------	-----------	----------	------

System setting interface

System:

Head	Tube
ExInput	ExOutput
Uart	SuckBack
Language	Exit

Head Select the current use of the pump head models.

Main interface → Click  → **system** → Click **Head** → Select pump head model → Click **Enter**, save and exit; Click **Cancel**, exit without saving.

Tube Select the current use of the pump tubing models.

Main interface → Click  → **system** → Click **Tube** → Select pump tube model → Click **Enter**, save and exit; Click **Cancel**, exit without saving.

Note: If you only replace the tube, you do not need to select the pump head model again.

Choose a different type of pump head, and the matching tube model will automatically change to provide options.

Enter Cancel

Tube :

0. 5*0. 8
1*1
2*1
2. 4*0. 8

▲
▼

Enter Cancel

ExInput Select External control input parameters

Main interface → Click  → **system** → Click **ExInput** → External control input option interface.

Ctrl → Usrc/Analog → Click **Enter**, save and exit; Click **Cancel**, exit without saving.

Item → Start Stop/Direction/Speed → Pinchar **Enter**, save and exit; Click **Cancel**, exit without saving.


Signal → Level/Pulse → Click **Enter**, save and exit; Click **Cancel**, exit without saving.

Speed → 0-5V/0-10V/4-20mA/0-10kHz → Click **Enter**, save and exit; Click **Cancel**, exit without saving.

After the pump is turned on, the LCD display is powered on to enter the main interface, and you can start specific settings and operations.

The default factory settings for the first boot

- Factory setting: CT Series peristaltic pump, the pump head model/tube model has been set before leaving the factory, in accordance with the customer's request. If multiple types of tubes are needed, the model set is the tube with the largest diameter.
- After the pump is turned on, it starts to run with the default set values, and all operating parameters can be changed by touch screen operation (see 5.4 "Menu function operation").
Note: It supports communication in advance, and the factory settings can be adjusted according to user needs.

5.3. Operation panel and display**(1) Shortcut key area**

-  **Full Speed function:** Press this button for fast filling/emptying/cleaning of the tube (at maximum speed); press it again to return to the original state.
-  **HOME:** In other operation interfaces, press and quickly return to the main interface.
-  **Direction:** Change the direction of drive rotation, means counterclockwise rotation and means clockwise rotation.
-  **Start & Stop:** Each time it is pressed, the start/stop state is changed.

(2) Graphic menu area

Set

Data

Lock

Voice

Remote

Set: Function operation main menu, peristaltic pump can be set accordingly

Data: Save new data quickly and apply saved data


Lock: Lock screen, unable to touch screen operation / Unlock, touch screen operation

Voice: Voice mode / Silent mode

Remote: Pump operation / Remote operation

(3) Flow parameter input area


Mode 1: Transmission (default mode; on the display is shown as Transfer mode)

Flow rate value (ml/min)

Liquid delivery volume (OFF) Running time (OFF)


Mode 2: Quantitative (on the display is shown as Ration mode)

Flow rate value (ml/min)

Liquid delivery volume (ml) Running time (OFF)


Mode 3: Timing

Flow rate value (ml/min)

Liquid delivery volume (OFF) Running time (sec/minutes)


Mode 4: Timing quantitative (on the display is shown as Filling mode)

Flow rate value (ml/min)

Liquid delivery volume (ml) Running time (sec/minutes)

(4) Status parameter display area

Pump head model / Cumulative liquid volume (ml) / Running time

Pump Tube Model / Speed (rpm) / Mode Display

Note: A single start and stop is a complete timing measurement, stop and then restart, and re-time measurement.

(5) Image indication

1. Running animation

Display the running status of the pump: display rotation when running; display no rotation when stopped.

2. Direction/full speed indication, ► Clockwise, ◄ Counterclockwise, ►► Full speed-clockwise, ◄◄ Full speed-counterclockwise.

5.4. Menu function operation

Flow parameter setting instructions

Delivery volume → target demand for a single delivery

Transfer time → the time required to complete the target transfer volume

Stop time → the stop time between transfers

Fluid volume setting

Fluid volume:			
Max:9999L		Min:0.001ml	
124.5 ml			
1	2	3	C
4	5	6	U
7	8	9	Esc
0	.	OK	

1. Main interface → click on the value of the delivered liquid volume, display the keypad → input the value (as shown in the figure)
2. Max/Min: - Displays the upper and lower setting limits of the current pump head and pump tube, which is the input range of the liquid delivery volume.
3. U-unit setting: ml/liter
4. C-Clear Esc-Exit OK-Confirm | **Note:** The pump does not support setting when it is running

Flow rate setting

Flow range:			
Max:9999L		Min:0.001ml	
124.5 ml			
1	2	3	C
4	5	6	U
7	8	9	Esc
0	.	OK	

1. Main interface → click the flow rate value to display the keypad → input the value (as shown in the figure)
2. Max/Min: - Displays the upper and lower set limits of the current pump head and pump tube, which is the input range of the flow rate in ml/min. According to different pump heads, pump tubes, the maximum and minimum values displayed by the keypad are also different.
3. U: not shown in this keypad
4. C-Clear Esc-Exit OK-Confirm | **Note:** The pump does not support setting when it is running

Running time setting

Running time:			
Max:9999min		Min:0.1sec	
0 min			
1	2	3	C
4	5	6	U
7	8	9	Esc
0	.	OK	

1. Main interface → click on the running time value, display the keypad → enter the value (as shown in the figure)
2. Max/Min: Displays the current running time range 0.1s-9999m
3. U-unit: s, m (seconds, minutes)
4. The value sets the running time of the filling. You can click to set when the filling stops, but cannot be set during the filling cycle. | **Note:** The pump does not support setting when it is running