# LABORGLASKATALOG

# LABORATORY GLASSWARE CATALOGUE

**favs** S.N.C.

di GNUDI GIOVANNI & C.

40127 Bologna - Via Tartini, 14/A  
Tel. 051501153 - Fax 0516336182  
www.favs.it - info@favs.it

**DURAN GROUP**  
magic of precision

A man with short dark hair and a light beard, wearing a white lab coat over a white shirt and a blue tie. He has his arms crossed and is smiling slightly at the camera. The background is a blurred cityscape seen through a window.

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“I'D LIKE TO SAY THANKS FOR THE CONSTANT  
SERVICE, FOR THE RELIABILITY, FOR **THE MANY  
DIFFERENT OPTIONS** ... IN SHORT, CAN  
ANYBODY TELL ME HOW I CAN SAY THANKS  
TO SOME GLASSWARE?“

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
  
Michael Merz  
General Manager Marketing & Sales

  
Tobias A. Thiele  
Vice President Marketing

# WELCOME TO THE DURAN GROUP.

We are delighted with your interest in our chemical and heat-resistant DURAN® laboratory glassware.  
Based on the know-how of more than 650 employees at sites in Mainz and Wertheim in Germany  
and Pula in Croatia and with more than 100 years experience in manufacturing borosilicate glass 3.3,  
you will find we are the right contact for standard laboratory glassware as well as exacting,  
customized, special articles.

With our comprehensive range and our wide range of production and finishing options, we are sure  
we can find, in cooperation with you, the optimum laboratory glass for your special application.

This product catalogue provides an overview of our products and services within the field of DURAN®  
laboratory glassware. Our product managers and sales staff will be happy to discuss your individual  
needs as well as the properties of our speciality material, DURAN®.

You can find more information and details for your contact partner on our website  
[www.duran-group.com](http://www.duran-group.com).

Please do not hesitate to contact us.

  
Michael Merz  
General Manager Marketing & Sales

  
Tobias A. Thiele  
Vice President Marketing

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**DURAN GROUP. Magic of precision.**

## THERE ARE GOOD REASONS FOR GETTING TO KNOW US.

*The DURAN GROUP is one of the world's leading manufacturers of borosilicate glass, a special glass invented by Otto Schott in 1887 and registered in 1938 under the trade name DURAN®.*

*With sites in Mainz (DURAN Produktions GmbH & Co. KG), Wertheim (DURAN Group GmbH) and Pula (DURAN d.d.), the DURAN GROUP has a complete value-added chain from the glass melting process right through to glass forming and processing.*

*Over the last 70 years, DURAN® glass has proven itself both in laboratories as well as in industrial and domestic applications. At the same time our special industrial glass is used in numerous industries, amongst others, mechanical engineering, the electrical industry and medical technology. Well-known companies spread across all continents trust in the quality of our products and the creativity of our employees. Based on close dialogue with our customers, we achieve solutions, which meet individual requirements and permit highly precise applications in widely varying sectors.*

*After its spin-off from SCHOTT AG, the DURAN GROUP currently combines the strengths of its long history with the constantly present visionary spirit of the inventor Otto Schott.*

*The unique precision involved in the manufacturing and processing of DURAN® glass is the foundation of our corporate philosophy and guideline for our day-to-day operations. We offer our customers the flexibility of a medium-sized company with comprehensive know-how built up from many years of experience.*

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QUALITY WITHOUT COMPROMISE.

*The requirement of our customers for our company to develop and manufacture reliable and safe products in accordance with the highest possible quality standards is at the very centre of our quality policy.*

*In close cooperation with all our staff and with the active involvement of our customers and suppliers, the DURAN GROUP has established a quality management system that conforms to DIN EN ISO 9001:2000 and which is integrated into daily practice.*

*This quality management system is the foundation for all the processes that our products have to go through: from the customer's initial enquiry, through to order processing and delivery and up to customer feedback.*

We value the success of DURAN® products as a sign of confidence from our customers in our quality system, in our logistics and in our service.

Moreover, we are constantly working on improving our performance by assessment of the status quo and from it deriving new, ambitious goals, which we want to attain through process-orientated thinking, planning and action.

*In everything we do, customer and staff satisfaction, as well as adherence to the principles of environmental sustainability are the important criteria that we constantly measure and aim to improve upon.*

*Safety and environmental protection as corporate goals.*

Human safety and protection of the environment are important goals for our company. To achieve this, we are constantly working on the development of more environmentally sustainable products and manufacturing processes while ensuring we use resources in a sparing and responsible manner and set high standards for the safe operation of our plant.

*The DURAN GROUP is committed to responsible and fair behaviour towards society and the environment. Legislative guidelines and official ordinance are considered a minimum requirement. In our daily operations we go well beyond this minimum.*

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## DURAN® LABORATORY GLASSWARE WITH BATCH IDENTIFICATION AND QUALITY CERTIFICATE.

*Increasing requirements and standards in terms of quality assurance as well as the traceability of primary packaging and auxiliary materials (EN-ISO 9000/8402, GMP, EU 178/2002) are becoming more and more significant.*

*To meet these requirements, all DURAN® laboratory bottles, beakers, Erlenmeyer flasks and premium caps are provided with a **Retrace Code**.*

*This eight-digit batch code (two digits for premium caps) allows the products to be retraced to the point of production and the matching batch. The retrace code is a contribution to a continuous documentation within the user's quality management system and it is therefore particularly important in the areas of medicine, the pharmaceutical and food industries.*

*By using this "Retrace Code" and the corresponding DURAN® item number (directly or via a selection menu), users can readily retrieve a batch and quality certificate via the internet. Besides the production date and the signature of the Quality Manager, the certificate also provides information on conformity with standards and USP/EP/JP.*

*Please see [www.duran-group.com](http://www.duran-group.com) for more information about the Retrace Code.*

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“MY BEST ASSISTANT ... **CHEMICALLY RESISTANT,**  
**PHYSICALLY RESISTANT ... OR BETTER PUT:**  
UNBEATABLE.”

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I | LABORATORY GLASS BOTTLES AND ACCESSORIES

Due to their high chemical and thermal-shock resistance, DURAN® laboratory bottles are indispensable for daily work in the laboratory. As neutral glass with high hydrolytic resistance, DURAN® belongs to glass type I in accordance with EP, JP and USP.

The extensive range of original accessories includes caps for the widest possible range of applications. DURAN® laboratory bottles are completed by suitable pouring rings from different plastics, which enable drip-free working.

Alongside clear glass laboratory bottles, amber (absorbs light up to approx. 500 nm) and plastic-coated (absorbs light up to approx. 380 nm) bottles are also available to protect light-sensitive substances. The PU plastic coating also provides protection against scratching and, in the event of breakage, retains contents and glass splinters.

Usage tips:

- DURAN® laboratory bottles are, with the exception of the pressure-resistant DURAN® pressure plus bottles, not designed for use under pressure or vacuum conditions.
- When sterilising / autoclaving, the screw cap should

only be loosely applied since with a closed bottle pressure equalisation does not take place. An ideal complement is the membrane cap (see pages 18, 30).

- When cleaning in a dishwasher, load such that the glass bottles, especially the threads, are not in contact with each other.
- The printed scale on many items of DURAN® laboratory glassware is indicative with an accuracy of  $\pm 10\%$ . Therefore the items are not suitable for use as volumetric glassware.
- DURAN® protect bottles can be autoclaved. The maximum permissible temperature is  $+135^\circ\text{C}$  (max. 30 Min.).
- When freezing media, store the bottles in an inclined position and do not fill more than three quarters full. Bear in mind that plastic caps and pouring rings have restricted service temperature ranges.

DURAN® laboratory glass and reagent bottles are provided with a retrace code. Using the eight-character code and the corresponding article number, a batch and quality certificate for the DURAN® laboratory glass bottles can be obtained at [www.duran-group.com](http://www.duran-group.com).

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
With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable white ceramic. Complete with blue quick release closure (PP<sup>1</sup>, integral lip seal) and pouring ring (PP<sup>1</sup>) for drip-free pouring and clean, safe working. Service temperature limit of closure and pouring ring:  $+140^\circ\text{C}$ . Special thread means opening takes less than a turn. The 80 mm wide outer diameter of the bottle neck permits easy filling and pouring out of powders and viscous substances.

Typical applications: storage, transport, safekeeping and sampling of substances, easy to use with granulated material, powders and viscous media, sampling of hot media.

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DURAN® GLS 80 LABORATORY BOTTLE, WIDE NECK  
with GLS 80 thread

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DURAN® GLS 80 LABORATORY  
BOTTLE, WIDE NECK, AMBER  
with GLS 80 thread

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
With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable white ceramic. Complete with blue quick release closure (PP<sup>1</sup>, integral lip seal) and pouring ring (PP<sup>1</sup>) for drip-free pouring and clean, safe working. Service temperature limit of closure and pouring ring: +140 °C. Alongside easy handling, UV protection up to 500 nm. Unchanged DURAN® properties within the bottle, as colouration is only to the outer surface. Very uniform, durable and chemically resistant amber colour due to use of innovative technology.

Typical applications: storage, transport and safekeeping of light-sensitive substances, easy to use with granulated material, powders and viscous media.

<sup>1</sup> For chemical and thermal resistance, see page 208.

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DURAN® GLS 80 PROTECT  
LABORATORY BOTTLE  
with GLS 80 thread, plastic coated

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With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable white ceramic. Complete with blue quick release closure (PP<sup>1</sup>, integral lip seal) and pouring ring (PP<sup>1</sup>) for drip-free pouring and clean, safe working. Service temperature limit of closure and pouring ring: +140 °C. Service temperature limit of the PU<sup>1</sup> plastic coating: -30 °C to +135 °C. The coating provides scratch, leak and splinter protection and is ideally suited to both the transport and storage of hazardous media or valuable samples. UV protection up to approx. 380 nm wavelength. High transparency. Suitable for microwaving.

Typical applications: storage, transport and safe handling of hazardous substances.

<sup>1</sup> For chemical and thermal resistance, see page 208.

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
Permits opening and closing of the DURAN® GLS 80 bottle with only a three-quarter turn. A matching PP<sup>1</sup> pouring ring is also available, permitting clean, drop-free use.

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<sup>1</sup> For chemical and thermal resistance, see page 208.


QUICK RELEASE CLOSURE FOR  
DURAN® GLS 80 LABORATORY BOTTLE  
from PP<sup>1</sup>, blue, with lip seal

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QUICK RELEASE CLOSURE FOR  
DURAN® GLS 80 LABORATORY BOTTLE  
with sealing disc, (PSU compound<sup>1</sup>)

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
The material used is a special compound based on a polyarylsulphone I. Consequently the chemical, thermal and mechanical properties of the material are noticeably improved and matched to laboratory requirements. Thanks to the thread, the DURAN® GLS 80 bottle can be opened and closed with only a three-quarter turn. The seal, which is coated on both sides with PTFE<sup>1</sup>, ensures the bottle can be tightly closed. A matching PTFE<sup>1</sup> pouring ring is also available, permitting clean, drop-free use.

<sup>1</sup> For chemical and thermal resistance, see page 208.

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MEMBRANE SCREW CAP FOR  
DURAN® GLS 80 LABORATORY BOTTLE  
from PP<sup>1</sup>, blue, with welded-in PTFE<sup>1</sup> membrane  
for pressure equalisation

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For GLS 80 thread. Ideal for autoclaving processes because the membrane permits pressure equalisation and tight sealing. Hence the risk of contamination is greatly reduced. Ingress of liquids or solids is prevented and the bottle contents remain sterile.

Typical applications: storage or transport of gas generating media, autoclaving of media.

<sup>1</sup> For chemical and thermal resistance, see page 208.

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Materials used: PP<sup>1</sup> and PTFE<sup>1</sup>. Flexible modular system. Five different hose diameters (3.2 mm; 6.0 mm; 8.0 mm; 10.0 mm and 12.0 mm) can be connected. Sterile pressure equalisation is possible by using the membrane filter. Unused ports can be provided with a blind cap.

Typical applications: safe transfer of liquid media within a closed and sterile system (evaporation is reduced).

<sup>1</sup> For chemical and thermal resistance, see page 208.

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
GLS 80 CONNECTION SYSTEM  
screw cap GLS 80 with four ports (GL 18 thread)

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Schematic diagram of GLS 80 connection system

- ① z. B. Best.-Nr. 11 601 76    ① e.g. Cat. No. 11 601 76  
② z. B. Best.-Nr. 11 601 69    ② e.g. Cat. No. 11 601 69  
③ z. B. Best.-Nr. 11 601 74    ③ e.g. Cat. No. 11 601 74

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
GLS 80 STIRRED REACTOR

materials used: PP<sup>1</sup>/PTFE<sup>1</sup>/PEEK<sup>1</sup>/stainless steel

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The GLS 80 stirred reactor is suitable for a wide range of laboratory mixing processes. The connections (2 x GL 14, 2 x GL 18) provided permit addition or removal of media from the bottle during the mixing process. The whole assembly can be autoclaved and is therefore suitable for use in the biological sector. By using components from the GLS 80 connection system, an additional media bottle (OD hose: 1.6 - 12.0 mm) can be connected or a sterile pressure equalizer attached. Drive for the stirrer is provided by a standard commercial magnetic stirrer. The variable stirrer shaft can be used in DURAN® GLS 80 laboratory glass bottles (1 000 and 2 000 ml) and provides notably improved mixing in comparison with standard magnetic stir bars. The agitating element is exchangeable and use up to 500 rpm is possible.

Typical applications: mixing of liquids, dissolving of solids, simple fermentation processes.

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DURAN® LABORATORY BOTTLE  
with DIN thread, GL 45

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*With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable white ceramic. With proven DURAN® properties. Complete with blue screw cap (PP<sup>1</sup>, integral lip seal) and pouring ring (PP<sup>1</sup>) for drip-free pouring and clean, safe working. Service temperature level of screw cap and pouring ring: + 140 °C. For screw cap assortment, see pages 29–36.*

*Typical applications: storage, sample preparation, transport.*


Best.-Nr. Cat. No.	Inhalt Capacity ml	DIN-Gewinde DIN Thread GL	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
Mit Schraubverschluss und Ausgießring aus PP (blau), klar With screw-cap and pouring ring from PP (blue), clear					
21 801 14 5 <sup>2</sup>	25	25	36	74	10
21 801 17 5	50	32	46	91	10
21 801 24 5	100 <sup>3</sup>	45	56	105	10
21 801 36 5	250 <sup>3</sup>	45	70	143	10
21 801 44 5	500 <sup>3</sup>	45	86	181	10
21 801 54 5	1 000 <sup>3</sup>	45	101	230	10
21 801 63 5	2 000 <sup>3</sup>	45	136	265	10
21 801 73 5	5 000	45	182	335	1
21 801 86 5	10 000	45	227	415	1
21 801 88 5	15 000	45	268	450	1
21 801 91 5	20 000	45	288	510	1
Ohne Verschluss und Ausgießring, klar Without cap and pouring ring, clear					
21 801 14 <sup>2</sup>	25	25	36	70	10
21 801 17	50	32	46	87	10
21 801 24	100 <sup>3</sup>	45	56	100	10
21 801 36	250 <sup>3</sup>	45	70	138	10
21 801 44	500 <sup>3</sup>	45	86	176	10
21 801 54	1 000 <sup>3</sup>	45	101	225	10
21 801 63	2 000 <sup>3</sup>	45	136	260	10
21 801 73	5 000	45	182	330	1
21 801 86	10 000	45	227	410	1
21 801 88	15 000	45	268	445	1
21 801 91	20 000	45	288	505	1

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<sup>1</sup> For chemical and thermal resistance, see page 208.  
<sup>2</sup> With specially shaped glass edge for improved pouring out (so that an additional plastic pouring ring is not required).  
<sup>3</sup> Bottles of 100 – 2 000 ml capacity have a shoulder reinforcement ring corresponding to the bottles nominal capacity.

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DURAN® LABORATORY BOTTLE, AMBER  
with DIN thread, GL 45

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*With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable white ceramic. UV protection up to approx. 500 nm wavelength. Unchanged DURAN® properties within the bottle, as colouration is only to the outer surface. Very uniform, durable and chemically resistant amber colour due to use of innovative technology. For screw cap assortment, see pages 29–36.*

*Typical applications: storage and transport of light-sensitive substances.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	DIN-Gewinde DIN Thread GL	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
Braun   Amber					
21 806 14 <sup>1</sup>	25	25	36	70	10
21 806 17	50	32	46	87	10
21 806 24	100 <sup>2</sup>	45	56	100	10
21 806 36	250 <sup>2</sup>	45	70	138	10
21 806 44	500 <sup>2</sup>	45	86	176	10
21 806 54	1 000 <sup>2</sup>	45	101	225	10
21 806 63	2 000 <sup>2</sup>	45	136	260	10
21 806 73	5 000	45	182	330	1
21 806 86	10 000	45	227	410	1
21 806 88	15 000	45	268	445	1
21 806 91	20 000	45	288	505	1

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
<sup>1</sup> With specially shaped glass edge for improved pouring out (so that an additional plastic pouring ring is not required).  
<sup>2</sup> Bottles of 100 – 2 000 ml capacity have a shoulder reinforcement ring corresponding to the bottles nominal capacity.

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DURAN® PROTECT  
LABORATORY BOTTLE  
with DIN thread, GL 45, plastic coated

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With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable white ceramic. Service temperature limit of the PU<sup>1</sup> plastic coating: –30 °C to +135 °C (instructions for autoclaving, see page 217). The coating provides scratch, leak and splinter protection and is ideally suited to both the transport and storage of hazardous media or valuable samples. UV protection up to approx. 380 nm wavelength. High transparency. Suitable for microwaving. Caps and pouring rings see pages 29–36.

Typical applications: storage, transport and safe handling of hazardous substances.


Best.-Nr. Cat. No.	Inhalt Capacity ml	DIN-Gewinde DIN Thread GL	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
Ohne Verschluss und Ausgießring, klar Without cap and pouring ring, clear					
10 926 76 <sup>2</sup>	25	25	36	70	10
10 926 77	50	32	46	87	10
21 805 24	100 <sup>3</sup>	45	56	100	10
21 805 36	250 <sup>3</sup>	45	70	138	10
21 805 44	500 <sup>3</sup>	45	86	176	10
21 805 54	1 000 <sup>3</sup>	45	101	225	10
21 805 63	2 000 <sup>3</sup>	45	136	260	10
21 805 73	5 000	45	182	330	1

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DURAN® PRESSURE PLUS  
LABORATORY BOTTLE<sup>1</sup>  
with DIN thread, GL 45

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With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable ceramic. Pressure resistance conforms to DIN ISO 1595, confirmed by GS marking (TÜV ID: 0000020716). Vacuum and/or pressure resistant from –1 to +1.5 bar due to modified geometry (based on ISO 4796-1). When pressure loaded the following apply: thermal shock resistance 30 K and maximum usage temperature +140 °C. Blue scale for visual differentiation from the standard laboratory bottle. Also available in amber. See pages 29-36 for spare caps and pouring rings.

Typical applications: safe working under pressure or vacuum, sampling under pressure, storage of gas generating media.


Best.-Nr. Cat. No.	Inhalt Capacity ml	DIN-Gewinde DIN Thread GL	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
Ohne Verschluss und Ausgießring, klar Without cap and pouring ring, clear					
10 922 34	250	45	70	138	10
10 922 35	500	45	86	176	10
21 810 54	1 000	45	101	225	10
Braun   Amber					
10 943 67	250	45	70	138	10
10 943 68	500	45	86	176	10
21 816 24	1 000	45	101	225	10

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<sup>1</sup> Bottle with plastic coating is available on demand.

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DURAN® HPLC BOTTLE  
with DIN thread, GL 45

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With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable ceramic. Complete system comprising: DURAN® pressure plus laboratory bottle with 4-port screw cap (PP<sup>1</sup>); four screw connections (black, M8 thread) and silicone seals. Connection of different hose diameters (1.6 mm and 3.2 mm) as well as sterile pressure equalisation sets (membrane filter 0.2 µm) is possible. Unused ports can be sealed with silicone blanking seals.

Typical applications: safe transfer of liquid media within a closed and sterile system (evaporation is reduced).


Best.-Nr. Cat. No.	Inhalt Capacity ml	DIN-Gewinde DIN Thread GL	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
I I 298 21 <sup>3</sup>	500	45	86	176 <sup>2</sup>	2
I I 298 20 <sup>3</sup>	1 000	45	101	225 <sup>2</sup>	2
Zubehör   Accessories					
I I 298 12	Schraubverschluss HPLC, GL 45, 4 Ports, komplett <sup>4</sup> Screw cap HPLC, GL 45, 4 ports, complete <sup>4</sup>				2
I I 298 13	Ersatzset <sup>5</sup> für HPLC-Schraubverschluss Spare part set <sup>5</sup> for HPLC screw cap				1
I I 378 01	Druckausgleichset für 4-Port-Verschluss (0,2 µm inkl. Membranfilter) Pressure equalization set 4-port cap (0.2 µm incl. Membrane filter)				1
I I 298 19	Ersatz-Membranfilter für Druckausgleichset, 0,2 µm Spare membrane filter for pressure equalization set, 0.2 µm				2

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<sup>1</sup> For chemical and thermal resistance, see page 208.  
<sup>2</sup> Bottle only.  
<sup>3</sup> Replacement bottle, see page 25.  
<sup>4</sup> Including silicone seals.  
<sup>5</sup> GL 45 screw cap, 4 x M8 screw cap, 12 x silicone seal.

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DURAN® LABORATORY BOTTLE,  
SQUARE  
with DIN thread, GL 45

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With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable white ceramic. Complete with blue screw cap (PP<sup>1</sup>, integral lip seal) and pouring ring (PP<sup>1</sup>) for drip-free pouring and clean, safe working. Service temperature limit of screw cap and pouring ring: +140 °C. Ergonomic handling due to angular shape, highly stable, good stackability. Alongside proven DURAN® properties, a space saving of 44% in comparison with standard laboratory bottles (example applies to 100 ml bottles). Screw caps are also available in the following colours: green, yellow and grey (see pages 29–36).

Typical applications: space-saving storage, space-saving transport.

Best.-Nr. Cat. No.	Inhalt Capacity ml	DIN-Gewinde DIN Thread GL	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
21 820 24 5	100	32	50	109	10
21 820 36 5	250	45	64	143	10
21 820 44 5	500	45	78	181	10
21 820 54 5	1 000	45	94	222	10

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*With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable white ceramic. Complete with pouring ring and cap from TpCh260<sup>1</sup> (similar to PFA). The premium cap with its PTFE<sup>1</sup> coated silicone seal is colourless and temperature resistant from –196 °C to +260 °C. Together with proven DURAN® properties, TÜV tested thermal shock resistance of 160 K, confirmed by GS-marking (TÜV ID: 0000020715). USP/FDA conformity for the entire system comprised of bottle, screw cap and pouring ring. DMF available upon request (DMF no. 19757). Accurate scale: ± 5 %. Additional graduations as well as additional opposing scale simplify reading off.*

*Typical applications: due to its properties, ideal for applications in the pharmaceutical industry, handling of aggressive media, sterilisation processes (hot air and dry sterilisation) and depyrogenation.*

<sup>1</sup> For chemical and thermal resistance, see page 208.

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*Compound formulated without colour dye to avoid risk of leaching. Service temperature limit from –196 °C to +260 °C. Very high chemical resistance. High leak tightness through use of PTFE<sup>1</sup> coated silicone seal. DMF (no. 19757, in conjunction with premium bottle) available upon request. A matching TpCh260<sup>1</sup> pouring ring is also available, permitting clean, drop-free use.*

*Typical applications: due to its properties, ideal for applications in the pharmaceutical industry, handling of aggressive media, demanding sterilisation processes such as hot air or dry sterilisation and depyrogenation.*

<sup>1</sup> For chemical and thermal resistance, see page 208.

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*The tamper-evident screw cap is available with lip seal or PTFE<sup>1</sup> coated silicone seal. The integral plastic ring tears when the cap is first opened and is retained on the bottle neck. Therefore it can be clearly seen whether the bottle retains its original seal. After first untwisting, can be used as “normal” screw cap.*

*Typical applications: secure storage and transport/shipping of valuable media.*

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

PREMIUM CAP  
from TpCh260<sup>1</sup> (similar to PFA) with  
PTFE<sup>1</sup> coated silicone seal

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
TAMPER-EVIDENT SCREW CAP  
from PP<sup>1</sup>, for laboratory bottles with DIN thread

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MEMBRANE SCREW CAP  
from PP<sup>1</sup>, blue with welded-in PTFE<sup>1</sup>  
membrane for pressure equalisation

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
For all GL 45, GL 32 and GL 25 threads.  
Ideal for autoclaving processes because  
the membrane permits pressure equalisation  
and tight sealing, greatly reducing the risk of  
contamination. Ingress of liquids or solids is  
prevented and the bottle contents remain  
sterile.

Typical applications: storage or transport  
of gas generating media, autoclaving of  
media.

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
Available in the colours blue, green, yellow  
and grey with matching colour pouring  
rings. Distinguishing media types is simplified  
and interchanging of screw caps and carry  
over of substances is effectively excluded.

Typical applications: visual aid to the safe  
identification of different media.

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High leak tightness through use of  
PTFE<sup>1</sup> coated silicone seal. More chemically  
resistant than PP<sup>1</sup> screw cap. A matching  
ETFE<sup>1</sup> pouring ring is also available,  
permitting clean, drop-free use.

Best.-Nr. Cat. No.	DIN-Gewinde DIN Thread GL	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
Verschluss   Cap <sup>2</sup>				
29 240 13	25	33	23	10
29 240 19	32	42	26	10
29 240 28	45	54	28	10
Ausgießring   Pouring ring				
29 244 19	32		4	10
29 244 28	45		4	10

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<sup>1</sup> For chemical and thermal resistance, see page 208.

<sup>2</sup> For screw connection caps with aperture, silicone seals  
with aperture and septa, see pages 167–169.

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CONNECTION SYSTEM  
screw cap GL 45 with two or three ports (GL 14 thread)

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Materials used: PP<sup>1</sup> and PTFE<sup>1</sup>. Flexible modular system. Four different hose diameters (1.6 mm; 3.0 mm; 3.2 mm and 6.0 mm) can be connected. Sterile pressure equalisation is possible through use of membrane filter. Unused ports can be provided with a blind cap.

Typical applications: safe transfer of liquid media within a closed and sterile system (evaporation is reduced).

<sup>1</sup> For chemical and thermal resistance, see page 208.

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
Ordering example

You would like to equip two bottles with a 3-port connection system. You would like to work with two different tube diameters (1.6 mm and 3 mm) and pressure equalization. Consequently, you need the following individual components:

- Screw cap GL 45 3-port (Art. No. 11 297 51) 1 x [2 per pack]
- Screw cap GL 14 (Art. No. 11 298 14) 2 x [2 per pack]
- Insert for screw cap GL14 1.6 mm ID (Art. No. 11 298 15) 2 x [1 per pack]
- Insert for screw cap GL14 3.0 mm ID, (Art. No. 11 298 16) 2 x [1 per pack]
- Pressure equalization set (complete) (Art. No. 11 377 99) 2 x [1 per pack].

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Schematic diagram of GL 45 connection system

- ① z. B. Best.-Nr. 11 297 51      ① e.g. Cat. No. 11 297 51  
② z. B. Best.-Nr. 11 298 15      ② e.g. Cat. No. 11 298 15  
③ z. B. Best.-Nr. 11 298 14      ③ e.g. Cat. No. 11 298 14

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DG SAFETY CAPS

screw cap GL 45 with one, two or three ports

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High-quality and flexible connection system, specifically designed for handling aggressive, gas-generating media. Materials used: PP<sup>1</sup>, PBT<sup>1</sup> and PTFE<sup>1</sup>. All parts coming into contact with media are from PTFE<sup>1</sup>. Three different hose diameters (1.6 mm; 2.3 mm; 3.2 mm) can be connected. Pressure equalisation is possible through use of membrane filter. Outgassing of media is reduced by use of a one-way valve that decreases media loss and protects the user from vapours. A free moving PTFE<sup>1</sup> core prevents rotation of the hoses when screwing onto the DURAN® bottle. Unused ports can be provided with a blind cap. Optionally available with stopcock to prevent running dry of the hoses or air intakes during bottle changeover. DG safety caps are also available for bottles with standard ground joint.

Typical applications: safe transfer of aggressive media within a closed system (evaporation is reduced).

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The activated carbon filter prevents the escaping of hazardous gases and thus increases safety in the laboratory. Materials used: PBT<sup>1</sup> and PTFE<sup>1</sup>. Optional special connection for installations with a leakage system. A free moving PTFE<sup>1</sup> core prevents rotation of the hoses when screwing onto the DURAN® bottle. Unused ports can be provided with a blind cap.


Typical applications: safe storage and collection of liquid laboratory waste.

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DG SAFETY WASTE CAPS

additional connection option for activated carbon filter and leak hose, GL 45

A  
121 °C

Tmax.  
140 °C

<sup>1</sup> For chemical and thermal resistance, see page 208.  
<sup>2</sup> Complete with GL 45 screw cap, connection screws and conical nipple (ID 3.2 mm, PTFE).  
<sup>3</sup> Complete with GL 45 screw cap, connection screws, conical nipple (ID 3.2 mm, PTFE) and connection for leak hose.

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Without pouring ring and dust cover cap<sup>1</sup>.  
A maximum usage temperature of +100 °C  
is recommended. Thermal shock resistance  
30 K. Hydrolytic class 3.

Best.-Nr. Cat. No.	Inhalt Capacity ml	Gewinde Thread	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
Klarglas   Clear glass					
23 810 24 5	100	32 <sup>2</sup>	49	119	10
23 810 36 5	250	32 <sup>2</sup>	64	155	10
23 810 44 5	500	32 <sup>2</sup>	77	186	10
23 810 54 5	1 000	45 <sup>2</sup>	97	223	10
Braunglas   Amber glass					
23 816 24 5	100	32 <sup>2</sup>	49	119	10
23 816 36 5	250	32 <sup>2</sup>	64	155	10
23 816 44 5	500	32 <sup>2</sup>	77	186	10
23 816 54 5	1 000	45 <sup>2</sup>	97	223	10

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


SCREW CAP BOTTLE, ROUND, AMBER  
FROM SODA-LIME GLASS  
with thread

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Complete with pouring ring and dust cover  
cap<sup>1</sup>. A maximum usage temperature of  
+100 °C is recommended. Thermal shock  
resistance 30 K. Hydrolytic class 3.


Best.-Nr. Cat. No.	Inhalt Capacity ml	Gewinde Thread	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
Braunglas   Amber glass					
23 835 66 5	2 500	45 <sup>2</sup>	139	283	1

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TAMPER-EVIDENT SCREW CAP  
from PP<sup>1</sup>, blue, for soda-lime square screw  
cap bottle, narrow neck

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
<sup>1</sup> For chemical and thermal resistance, see page 208.  
<sup>2</sup> High form.

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SECURITY SCREW CAP  
from PP<sup>1</sup> with valve<sup>2</sup>, red, for soda-lime  
square screw cap bottles, narrow neck

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<sup>1</sup> For chemical and thermal resistance, see page 208.  
<sup>2</sup> Valve cap opens at 0.1 bar pressure.  
<sup>3</sup> High form.


Best.-Nr. Cat. No.	Gewinde Thread	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
29 303 19	32 <sup>2</sup>	44	23	10
29 303 28	45 <sup>2</sup>	58	27	10
29 303 32	54 <sup>2</sup>	69	29	10
29 303 35	60 <sup>2</sup>	78	29	10

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<sup>1</sup> For chemical and thermal resistance, see page 208.  
<sup>2</sup> Short form.

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TAMPER-EVIDENT SCREW CAP  
from PP<sup>1</sup>, blue, for soda-lime  
square screw cap bottles, wide neck

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DURAN® REAGENT BOTTLE,  
WIDE NECK  
neck with standard ground joint

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With precision ground neck. Only glass components, therefore also suitable for storage of aggressive media, which could attack plastic parts.

Typical application: storage of powders.

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DURAN® REAGENT BOTTLE,  
WIDE NECK  
neck with standard ground joint,  
with standard ground glass flat-head stopper

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With precision ground neck. Only glass components, therefore also suitable for storage of aggressive media, which could attack plastic parts. Also available in amber. Unchanged DURAN® properties within the bottle, as colouration is only to the outer surface. Very uniform amber colour thanks to innovative technology, moreover long durability and high chemical resistance of the amber colouring.

Typical applications: storage of powders.

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With precision ground neck. Only glass components, therefore also suitable for storage of aggressive media, which could attack plastic parts.

Typical application: storage of liquids.


<sup>1</sup> Narrow neck bottles with standard ground joint, 1 000 ml, also suitable for use as burette bottles, see page 113.

With precision ground neck. Only glass components, therefore also suitable for storage of aggressive media, which could attack plastic parts. Also available in amber. Unchanged DURAN® properties within the bottle, as colouration is only to the outer surface. Very uniform, durable and chemically resistant amber colour due to use of innovative technology.

Typical application: storage of liquids.

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
DURAN® REAGENT BOTTLE,  
NARROW NECK  
neck with standard ground joint

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
DURAN® REAGENT BOTTLE,  
NARROW NECK  
neck with standard ground joint,  
with standard ground glass flat-head stopper

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REAGENT BOTTLE, WIDE NECK  
FROM SODA-LIME GLASS  
neck with standard ground joint

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
Without stopper<sup>1</sup>. A maximum usage  
temperature of +100 °C is recommended.  
Thermal shock resistance 30 K. Hydrolytic  
class 3.

<sup>1</sup> Stopper see pages 44–45.  
<sup>2</sup> Neck NS 45/27

**Note** on using DURAN® glass stoppers  
with soda-lime glass bottles: If the bottle  
and the stopper have a temperature  
difference greater than 30 °C, the glass  
stoppers can become stuck!

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REAGENT BOTTLE, WIDE NECK  
FROM SODA-LIME GLASS  
neck with standard ground joint,  
with standard ground glass flat-head stopper<sup>1</sup>

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
REAGENT BOTTLE, NARROW NECK  
FROM SODA-LIME GLASS  
neck with standard ground joint

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A maximum usage temperature of +100 °C  
is recommended. Thermal shock resistance  
30 K. Hydrolytic class 3.

<sup>1</sup> Narrow neck bottles with standard ground joint, 1 000 ml,  
also suitable for use as burette bottles, see page 113.

**Note** on using DURAN® glass stoppers  
with soda-lime glass bottles: If the bottle  
and the stopper have a temperature  
difference greater than 30 °C, the glass  
stoppers can become stuck!

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
A maximum usage temperature of  
+100 °C is recommended. Thermal shock  
resistance 30 K. Hydrolytic class 3.

**Note** on using DURAN® glass stoppers  
with soda-lime glass bottles: If the bottle  
and the stopper have a temperature  
difference greater than 30 °C, the glass  
stoppers can become stuck!

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REAGENT BOTTLE, NARROW NECK  
FROM SODA-LIME GLASS  
neck with standard ground joint,  
with standard ground glass flat-head stopper

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
GLASS STOPPER, AMBER  
with standard ground joint, octagonal

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PLASTIC STOPPER  
from polyethylene<sup>1</sup>

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
GLASS STOPPER  
with standard ground joint, octagonal

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

GLASS STOPPER  
with short ground joint, octagonal

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DURAN® GLASS STOPPER  
ground conical, for reagent bottles  
(oxygen bottles according to Winkler)

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**Note** on using DURAN® glass stoppers with soda-lime glass bottles: If the bottle and the stopper have a temperature difference greater than 30 °C, the glass stoppers can become stuck!

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DURAN® LEVELLING BOTTLE  
with stopper-bed and olive

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
*Dosing of liquids is possible via an outlet.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
24 708 36	250	73	130	10
24 708 44	500	89	164	10
24 708 54	1 000	111	200	1

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DURAN® ACID BOTTLE<sup>1</sup>  
with standard ground “pennyhead” stopper,  
conical shoulders, interchangeable glass cap

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*In addition to the ground stopper, a glass cap with ground joint is supplied. This provides an improved seal and protection against acid vapours.*


Best.-Nr. Cat. No.	Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	h mm	NS	Verpackungseinheit/ Stück Pack/Quantity
klar   clear	braun   amber					
21 275 24	21 275 24 6	100	55	145	19/17	10
21 275 36	21 275 36 6	250	75	180	19/17	10
21 275 44	21 275 44 6	500	82	220	24/20	10
21 275 54	21 275 54 6	1 000	109	260	29/32	10

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DROPPING BOTTLE  
FROM SODA-LIME GLASS  
with interchangeable clear glass pipette<sup>1</sup> standard  
ground joint, complete with rubber teats<sup>2</sup>

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

DURAN® ASPIRATOR BOTTLE  
tubulated with standard ground joint,  
without stoppers, neck unground

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DURAN® ASPIRATOR BOTTLE  
tubulated with standard ground joint,  
complete with standard ground stopcock  
and standard ground stopper

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STOPCOCK WITH STANDARD  
GROUND JOINT FOR ASPIRATOR  
BOTTLE

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DURAN® ASPIRATOR BOTTLE  
neck with DIN thread GL 45,  
tabulator with GL 32

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*Dosing of liquids is possible via a stopcock.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	h mm	Bodentubus Tubulature NS	Hals Neck NS	Orig. Pack./Stück Stand. pack/quantity
24 702 44	500	86	164	19/26	24/29	10
24 702 54	1 000	107	200	19/26	29/32	10
24 702 63	2 000	134	249	19/26	29/32	10
24 702 73	5 000	182	320	29/32	45/40	1
24 702 86	10 000	228	398	29/32	60/46	1
24 702 91	20 000	289	492	29/32	60/46	1

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*Complete with screw connection cap<sup>1</sup>,  
silicone seal and stopcock with PTFE<sup>2</sup>  
spindle. Dosing of liquids is possible via a  
stopcock.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	DIN-Gewinde DIN Thread GL	d mm	h mm	Verpackungseinheit/ Stück Pack./Quantity
24 703 54	1 000	45	101	225	1
24 703 63	2 000	45	136	260	1
24 703 73	5 000	45	182	330	1
24 703 86	10 000	45	230	410	1

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
<sup>2</sup> For chemical and thermal resistance, see page 208.

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STOPCOCK FOR ASPIRATOR BOTTLE  
with PTFE<sup>1</sup> spindle, for GL 32 screw thread

Best.-Nr. Cat. No.	Bohrung Hole mm	für Inhalt for Capacity ml	Orig. Pack./Stück Stand. pack/quantity
24 147 03	6	1 000 + 2 000	1
24 147 04	8	5 000 + 10 000	1

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
PTFE-SILICONE<sup>1</sup> SEAL  
for GL 32 screw thread

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“LONDON, BOSTON, SHANGHAI ... IN REALITY,  
MASTERMINDS ARE VALUED ANYWHERE IN THE  
WORLD ... AS IS **GLASSWARE MADE IN GERMANY.**”

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2 | BOILING FLASKS AND GENERAL LABORATORY GLASSWARE

DURAN® laboratory glassware, including heating vessels, has very good thermal-shock resistance ( $\Delta T=100\text{ K}$ ) and a high operating temperature ( $+500\text{ °C}$ ). Not only the glass type but also its uniform wall thickness distribution are critical in preventing uneven expansion and stressing of the glass which could result in failure. For this reason, wall thickness distribution is, as a vital quality characteristic, continuously checked during the production process.

The beakers are primarily used as heating vessels. The tall shape is particularly suited to heating in liquid baths where the beaker contents are protected against the surrounding medium.

Erlenmeyer flasks are well suitable for mixing, because of their conical shape.

Weighing bottles are used when accurately weighing out substances. Close fitting lids with moulded grips are used to prevent the substances from being lost, e.g. during transport within the laboratory.

Watch glass dishes can be used both for covering beakers and Erlenmeyer flasks as well as for weighing small quantities of substances.

Our product range also includes a wide range of test tubes. In addition to DURAN® glass, other glass types are available (FIOLAX®, soda-lime). The characteristics of each glass type may be found on pages 203–207, 219.

Usage tips:

- Due to the uniform wall thickness distribution suitable for very high temperature changes.
- The printed scale on many items of DURAN® laboratory glassware is indicative with an accuracy of  $\pm 10\%$ . Therefore the items are not suitable for use as volumetric glassware.
- The products are not designed for use under differential pressure or vacuum conditions.

DURAN® beakers and Erlenmeyer flasks are provided with a retrace code. Using the eight-character code and the corresponding article number, a batch and quality certificate can be obtained at [www.duran-group.com](http://www.duran-group.com).

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
With easy-to-read scale and large labelling field for easy marking in fired-on, highly durable, white ceramic. Spout for clean pouring. Uniform wall thickness distribution makes these beakers ideal for heating applications.

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

DURAN® BEAKER  
high form, without spout

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
With easy-to-read scale and large labelling field for easy marking in fired-on, highly durable, white ceramic. Uniform wall thickness distribution makes these beakers ideal for heating applications.

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
21 117 17	50	38	70	10
21 117 24	100	48	80	10
21 117 29	150	54	95	10
21 117 36	250	60	120	10
21 117 41	400	70	130	10
21 117 48	600	80	150	10
21 117 54	1 000	95	180	10

DURAN® BECHER  
dickwandig (Filterierbecher)

DURAN® BEAKER  
heavy-wall (filtering beaker)

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With easy-to-read scale and large labelling field for easy marking in fired-on, highly durable, white ceramic. Has, due to the increased wall thickness, better mechanical properties than the standard beaker. Thermal shock resistance, however, is reduced so only limited application for heating. Spout for clean pouring.


Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
21 131 24	100	52	85	10
21 131 29	150	54	93	10
21 131 36	250	70	94	10
21 131 44	500	89	124	10
21 131 54	1 000	105	160	10
21 131 63	2 000	135	195	10
21 131 68	3 000	157	205	4
21 131 73	5 000	182	256	1
21 131 86 <sup>1</sup>	10 000	225	340	1
21 131 88 <sup>1</sup>	15 000	260	390	1
21 131 91 <sup>1</sup>	20 000	285	430	1

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

DURAN® BERZELIUS BEAKER  
without spout

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
With easy-to-read scale and large labelling field for easy marking in fired-on, highly durable, white ceramic. Due to conical form, suited to the mixing of liquids. Uniform wall thickness distribution makes these flasks ideal for heating applications.

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	d <sub>i</sub> mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
21 216 14	25	42	22	75	10
21 216 17	50	51	22	90	10
21 216 24	100	64	22	105	10
21 990 27 <sup>1</sup>	150	74	28	118	10
21 216 32 <sup>1</sup>	200	79	34	131	10
21 216 36	250	85	34	145	10
21 216 39 <sup>1</sup>	300	87	34	156	10
21 216 44	500	105	34	180	10
21 216 54	1 000	131	42	220	10
21 216 63	2 000	166	50	280	10
21 216 68	3 000	187	52	310	2
21 216 73	5 000	220	52	365	1

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With easy-to-read scale and large labelling field for easy marking in fired-on, highly durable, white ceramic. Due to conical form, suited to the mixing of liquids. Uniform wall thickness distribution makes these flasks ideal for heating applications. The wide neck enables easy filling and cleaning.

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	d <sub>i</sub> mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
21 226 14 <sup>1</sup>	25	43	31	70	10
21 226 17	50	51	34	85	10
21 226 24	100	64	34	105	10
21 226 32 <sup>1</sup>	200	79	50	131	10
21 226 36	250	85	50	140	10
21 226 39 <sup>1</sup>	300	87	50	156	10
21 226 44	500	105	50	175	10
21 226 54	1 000	131	50	220	10
21 226 63 <sup>1</sup>	2 000	153	72	276	10

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With easy-to-read scale and large labelling field for easy marking in fired-on, highly durable, white ceramic. The flask can be closed with a PBT cap<sup>1</sup> or membrane cap<sup>2</sup> (permits gas exchange).

Typical applications: The flask is suitable for storage, media preparation and cultivation.

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


DURAN® IODINE FLASK  
Erlenmeyer shape, with standard ground joint  
and glass stopper

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
With easy-to-read scale and large labelling field for easy marking in fired-on, highly durable, white ceramic. The flask can be closed with a glass stopper.

Typical applications: the iodine flask is suitable for determining the iodine number, i.e. the content of unsaturated fatty acids in oils and fats.

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DURAN® CONICAL FLASK  
Erlenmeyer shape, wide neck

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	d <sub>1</sub> mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
21 227 68 <sup>1</sup>	3 000	190	106	285	1
21 227 73 <sup>1</sup>	5 000	220	108	322	1
21 227 86 <sup>1</sup>	10 000	285	147	420	1

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<sup>2</sup> Nach ASTM E 1403.

Uniform wall thickness distribution makes these flasks ideal for heating applications. The geometry permits very uniform heating. Flasks with a neck diameter of 65 mm or more have a reinforced rim.

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
<sup>1</sup> Non-DIN EN ISO size.

Uniform wall thickness distribution makes these flasks ideal for heating applications. The geometry permits very uniform heating. The wide neck permits easy filling and removal of flask contents. Flasks with a neck diameter of 76 mm or more have a reinforced rim.

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
DURAN® ROUND BOTTOM FLASK  
narrow neck, with beaded rim

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
DURAN® ROUND BOTTOM FLASK  
wide neck, with beaded rim

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DURAN® FLAT BOTTOM FLASK  
narrow neck, with beaded rim

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Uniform wall thickness distribution makes these flasks ideal for heating applications. Flat base means flasks can be set down without a supporting ring. Flasks with a neck diameter of 65 mm have a reinforced rim.


Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	d <sub>1</sub> mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
21 711 17	50	51	26	90	10
21 711 24	100	64	26	105	10
21 711 36	250	85	34	138	10
21 711 44	500	105	34	163	10
21 711 54	1 000	131	42	190	10
21 711 64 <sup>1</sup>	2 000	166	42	250	10
21 711 68 <sup>1</sup>	3 000	185	50	250	1
21 711 71	4 000	207	50	275	1
21 711 73 <sup>1</sup>	5 000	223	50	290	1
21 711 76	6 000	237	65	315	1
21 711 86	10 000	280	65	360	1

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DURAN® FLAT BOTTOM FLASK  
wide neck, with beaded rim

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Uniform wall thickness distribution makes these flasks ideal for heating applications. Flat base means flasks can be set down without a supporting ring. The wide neck permits easy filling and removal of flask contents. Flasks with a neck diameter of 76 mm have a reinforced rim.


Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	d <sub>1</sub> mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
21 731 17	50	51	34	90	10
21 731 24	100	64	34	105	10
21 731 36	250	85	50	138	10
21 731 44	500	105	50	163	10
21 731 54	1 000	131	50	190	10
21 731 63 <sup>1</sup>	2 000	166	76	230	10
21 731 64	2 000	166	50	230	10

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DURAN® ENGLER FLASK  
with side outlet

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Uniform wall thickness distribution makes these flasks ideal for heating applications and distillations.

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	d <sub>1</sub> mm	h mm	Seitenrohr Side arm d <sub>2</sub> mm l mm	Verpackungseinheit/ Stück Pack/Quantity
21 653 24	100	66	20	215	6 100	10
21 653 28	125	68	22	215	7 100	10
21 653 29	150	73	20	223	6 100	10

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
DURAN® ENGLER FLASK  
with standard ground joint 19/26, side outlet

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

GLASKERAMIK LABOR-  
SCHUTZPLATTE

GLASS CERAMIC LABORATORY  
PROTECTION PLATE

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SQUARE QUADRUPOD  
for glass ceramic laboratory protection plate

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
Plate holder for glass ceramic plates. Made from heat-resistant chrome-nickel steel with four legs for extra stability.

Best.-Nr. Cat. No.	Für Platte For plates mm	Abmessungen der Platte Plate dimensions a x b mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
29 077 53	23 821 53	135 x 135	210	5
29 077 57	23 821 57	155 x 155	210	5
29 077 58	23 821 58	175 x 175	220	5

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PLATE HOLDER  
for glass ceramic laboratory protection plate

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
Plate holder for glass ceramic plates. Made from heat-resistant chrome-nickel steel.

Best.-Nr. Cat. No.	Für Platte For plates mm	Abmessungen der Platte Plate dimensions a x b mm	Verpackungseinheit/ Stück Pack/Quantity
29 078 53	23 821 53	135 x 135	5
29 078 57	23 821 57	155 x 155	5
29 078 58	23 821 58	175 x 175	5

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CLEANING SCRAPER  
for glass ceramic laboratory protection plate

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
Ideal for cleaning glass ceramic plates.

Best.-Nr. Cat. No.	Verpackungseinheit/Stück Pack/Quantity
29 079 01	10

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SPIRIT LAMP  
FROM SODA-LIME GLASS  
without filling tubulature, with ground over-cap

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
<sup>1</sup> Without print.  
<sup>2</sup> Without labelling field.

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
DURAN® EVAPORATING DISHES  
with spout

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DURAN® CRYSTALLIZING DISH  
with and without spout

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
DURAN® WEIGHING BOTTLE  
with ground lid

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
Close-fitting lid prevents any sample loss during transport after weighing. Available in low and high forms.

Best.-Nr. Cat. No.	Inhalt Capacity ml	d x h mm	Verpackungseinheit/ Stück Pack/Quantity
Niedrige Form Flat shape			
24 210 13	5	28 x 25	10
24 210 23	15	38 x 30	10
24 210 32	30	54 x 30	10
24 210 41	80	85 x 30	10
Hohe Form Tall shape			
24 211 13	10	28 x 40	10
24 211 18	20	32 x 50	10
24 211 23	45	38 x 70	10
24 211 23	70	44 x 80	10

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
Available in DURAN® and also in soda-lime glass (for properties, see pages 203–207, 219.

Best.-Nr. Cat. No.	d mm	Verpackungseinheit/Stück Pack/Quantity
DURAN®		
21 321 24	40	10
21 321 32	50	10
21 321 34	60	10
21 321 41	80	10
21 321 46	100	10
21 321 52	125	10
21 321 57	150	10
21 321 61	200	10
21 321 66	250	1
Kalk-Soda-Glas Soda-lime glass		
23 321 24	40	10
23 321 32	50	10
23 321 34	60	10
23 321 38	70	10
23 321 41	80	10
23 321 43	90	10
23 321 46	100	10
23 321 51	120	10
23 321 52	125	10
23 321 57	150	10
23 321 61	200	10
23 321 66	250	10

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
24 204 23	75	50	70	10
24 204 24	100	54	75	10

DURAN® ORGANGLAS  
ohne Stopfen


DURAN® ORGAN STORAGE JAR  
without stopper

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

DURAN® CYLINDER  
with knobbed lid, polished rim

Best.-Nr. Cat. No.	d x h mm	Volumen ca. Volume approx. ml	Verpackungseinheit/ Stück Pack/Quantity
24 205 01	80 x 80	250	10
24 205 03	100 x 100	500	10
24 205 05	120 x 120	1 000	1
24 205 10	150 x 150	2 000	1
24 205 21	210 x 210	6 000	1
24 205 32	260 x 260	12 000	1

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
The precise grinding of the knobbed lid and base vessel enables a very tight seal.

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DURAN® SPECIMEN JAR  
with ground, knobbed lid


Best.-Nr. Cat. No.	d x h mm	Volumen ca. Volume approx. ml	Verpackungseinheit/ Stück Pack/Quantity
24 209 02	65 x 63	125	10
24 209 09	65 x 103	250	10
24 209 11	110 x 103	800	10
24 209 16	85 x 123	600	10
24 209 17	127 x 123	1 400	1
24 209 24	85 x 153	750	1
24 209 26	110 x 153	1 250	1
24 209 28	157 x 153	2 500	1
24 209 38	110 x 203	1 750	1
24 209 39	157 x 203	4 000	1
24 209 49	110 x 253	2 000	1
24 209 50	157 x 253	4 500	1
24 209 57	120 x 300	2 250	1
24 209 59	252 x 303	13 000	1

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The precise grinding of the glass plate and base vessel enables a very tight seal.

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
DURAN® MUSEUM JAR  
with ground glass plate

Best.-Nr. Cat. No.	l x b x h mm	Verpackungseinheit/Stück Pack/Quantity
21 363 05	60 x 50 x 100	10
21 363 11	100 x 50 x 120	1
21 363 13	130 x 50 x 130	1
21 363 19	150 x 50 x 150	1
21 363 28	120 x 60 x 180	1
21 363 47	210 x 100 x 210	1
21 363 58	250 x 140 x 250	1

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
DURAN® MULTI-PURPOSE CYLINDER  
with round base

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
DURAN® STANDING CYLINDER  
with round base

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
Wall thickness and geometry designed  
to suit vacuum applications.

Best.-Nr. Cat. No.	h + h <sub>1</sub> mm	h x d mm	Verpackungseinheit/ Stück Pack/Quantity
24 460 59	300	250 x 185	1
24 460 66	305	255 x 260	1
24 460 69	350	300 x 315	1

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DURAN® BELL JAR  
with aperture in neck

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Wall thickness and geometry designed to  
suit vacuum applications. Neck aperture,  
standard ground joint 34/35.

Best.-Nr. Cat. No.	h + h <sub>1</sub> mm	h x d mm	NS	Verpackungseinheit/ Stück Pack/Quantity
24 465 59	300	250 x 185	34/35	1
24 465 61	350	300 x 215	34/35	1
24 465 69	550	500 x 315	34/35	1

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
For scale divisions and accuracy limits,  
see table:

Best.-Nr. Cat. No.	Inhalt Capacity ml	Größter d Max. d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
21 401 54	1 000	120	470 max.	10

Graduierung Scale ml	Teilung Division ml	Toleranz +/- Tolerance +/- ml
0-2	0,1	0,1
2-10	0,5	0,5
10-40	1	1
40-100	2	2
1 000	Ringmarke / Circular marking	10

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DURAN® SEDIMENTATION CONE  
Imhoff type, graduated

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
DURAN® TEST TUBE  
with beaded rim or straight rim

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
The test tubes are thick-walled and therefore mechanically very resistant, yet still retain good thermal shock resistance.

Best.-Nr. Bördelrand Cat. No. Beaded rim	Best.-Nr. Rand gerade Cat. No. Straight rim	d x h mm	ca. Volumen approx. vol. ml	Wanddicke Wall thickness mm	Verpackungseinheit/ Stück Pack/Quantity
26 130 01	26 131 01	8 x 70	2	0,8-1,0	100
26 130 03	26 131 03	10 x 75	4	0,8-1,0	100
26 130 06	26 131 06	10 x 100	5	0,8-1,0	100
26 130 08	26 131 08	12 x 75	6	0,8-1,0	100
26 130 11	26 131 11	12 x 100	8	0,8-1,0	100
26 130 13	26 131 13	14 x 130	16	0,8-1,0	100
26 130 16	26 131 16	16 x 130	17	1,0-1,2	100
26 130 21	26 131 21	16 x 160	21	1,0-1,2	100
26 130 23	26 131 23	18 x 180	32	1,0-1,2	100
26 130 26	26 131 26	20 x 150	34	1,0-1,2	100
26 130 28	26 131 28	20 x 180	40	1,0-1,2	100
26 130 33	26 131 33	25 x 150	55	1,0-1,2	50
26 130 36	26 131 36	25 x 200	70	1,0-1,2	50
26 130 38	26 131 38	30 x 200	100	1,0-1,4	50

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FIOLAX® BOROSILICATE TEST TUBE  
with beaded rim

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Thin-walled test tubes suited to rapid temperature changes or localized heating.

Best.-Nr. Cat. No.	d x h mm	ca. Volumen approx. vol. ml	Wanddicke Wall thickness mm	Verpackungseinheit/ Stück Pack/Quantity
26 110 01	8 x 70	2	0,4-0,5	100
26 110 03	10 x 75	4	0,4-0,5	100
26 110 06	10 x 100	6	0,4-0,5	100
26 110 08	12 x 75	6,5	0,4-0,5	100
26 110 11	12 x 100	9	0,4-0,5	100
26 110 13	14 x 130	16	0,4-0,5	100
26 110 16	16 x 130	20	0,5-0,6	100
26 110 21	16 x 160	25	0,5-0,6	100
26 110 23	18 x 180	35	0,5-0,6	100
26 110 26	20 x 150	39	0,5-0,6	100
26 110 28	20 x 180	45	0,5-0,6	100
26 110 33	25 x 150	60	0,6-0,7	50
26 110 36	25 x 200	80	0,6-0,7	50
26 110 38	30 x 200	110	0,7-0,8	50

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Centrifuge tubes are thick-walled and very resistant to mechanical loading (For further details needed for calculation of a maximum speed in rpm, see page 219). The higher density fraction collects in the bottom. Consequently solids can be collected and separated.

Best.-Nr. Cat. No.	d x h mm	Nenninhalt Nominal capacity ml	Verpackungseinheit/ Stück Pack/Quantity
21 601 10	12 x 100	6	50
21 601 11	16 x 100	12	50
21 601 14	24 x 100	25	10
21 601 17	34 x 100	50	10
21 601 24 <sup>1</sup>	40 x 115	80	10
21 601 26	44 x 100	80	10
21 601 36 <sup>1</sup>	56 x 147	250	10

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Centrifuge tubes are thick-walled and very resistant to mechanical loading (For further details needed for calculation of a maximum speed in rpm, see page 219). The higher density fraction collects in the pointed centre of the bottom. Consequently even small amounts of solids can be collected and separated.

Best.-Nr. Cat. No.	d x h mm	Nenninhalt Nominal capacity ml	Verpackungseinheit/ Stück Pack/Quantity
24 263 09	16 x 100	12	50

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




The straight rim permits the use of Kapsenberg caps; tubes are therefore well suited to the culture of micro-organisms.

Typical applications: growth and storage of sterile cultures.

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DURAN® CULTURE TUBE  
with DIN thread and screw cap  
from PBT<sup>1</sup>

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The DIN thread permits the use of PBT<sup>1</sup> screw caps; tubes are therefore well suited to the culture of micro-organisms. The contents only come into contact with the glass and PTFE<sup>1</sup> seal. For details of the PBT<sup>1</sup> screw cap, see pages 167–169.

Typical applications: growing and storage of sterile cultures.

<sup>1</sup> For chemical and thermal properties, see page 208.

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
DISPOSABLE CULTURE TUBE  
FROM SODA-LIME GLASS  
with DIN thread and PP<sup>1</sup> screw cap

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

SCREW CAP FOR CULTURE TUBES  
FROM SODA-LIME GLASS  
from PP<sup>1</sup> with seal

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DISPOSABLE CULTURE TUBE  
FROM SODA-LIME GLASS  
straight rim


Best.-Nr. Cat. No.	d x h mm	ca. Volumen approx. vol. ml	Wanddicke Wall thickness mm	Verpackungseinheit/ Stück Pack/Quantity
23 172 018	9,75 x 75	4	0,80	814
23 172 036	10,00 x 75	4	0,60	766
23 172 058	11,75 x 75	5	0,80	550
23 172 086	11,75 x 75	6	0,55	550
23 172 096	12,25 x 75	7	0,55	500
23 172 099	12,25 x 75	6	0,80	500
23 172 078	11,75 x 100	8	0,80	550
23 172 119	12,25 x 100	9	0,80	500
23 172 148	15,50 x 100	14	0,80	310
23 172 109	15,75 x 100	15	0,90	310
23 172 128	12,25 x 120	18	0,80	500
23 172 188	10,00 x 150	8	0,80	766
23 172 198	15,50 x 150	19	0,80	310
23 172 219	15,50 x 160	22	0,80	310

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<sup>1</sup> Other dimensions available upon request.

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NMR tubes are available, according to requirement, in three accuracy classes. The correct tube can be selected depending on magnetic field and spin. These tubes are noteworthy for their close tolerances and accuracy, especially to their straightness, wall thickness and wall thickness distribution. Consequently, quick and accurate test results are achievable.

Best.-Nr. Cat. No.	h mm	AD OD mm	ID ID mm	Gradheit Camber mm	MHZ mm	Verpackungseinheit/ Stück Pack/Quantity
Economic						
23 170 0117 <sup>1</sup>	178	4,95+/- 0,05	4,20 +/- 0,05	0,07	300	250
Professional						
23 170 0211	178	4,97+/- 0,025	4,20 +/- 0,025	0,03	400	250
Scientific						
23 170 0314	178	4,97+/- 0,013	4,20 +/- 0,025	0,013	500	5

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*“IS EVERYTHING MADE CORRECTLY? ... TO START  
WITH, WE HAD SO MANY QUESTIONS ... BUT AT  
LEAST ONE THING IS CLEAR: IT'S GOOD TO KNOW  
WE CAN RELY ON SOMEONE WHO HAS BEEN  
**AROUND FOR 100 YEARS.**”*

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All individual parts and a wide range of accessories such as lids, stopcocks, bases, etc. are compatible and can be interchanged as required. Always ensure the individual parts have the same DN (nominal diameter). For desiccators, the DN is based on the diameter of the sieve plate; this can be measured directly. For lids, measure the outside diameter of the flange and cross-reference with the tables on the product pages.

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An overview on page 80 indicates which individual parts are required to assemble the desired desiccator.

Usage tips:

- Designed for use under full vacuum.
- Due to the high wall thickness and the reduced thermal shock resistance under pressure loading, the desiccators must not be heated on one side only or heated using a naked flame.
- Before evacuation, it is recommended that the glass surfaces of the desiccator be checked for damage such as scratches, cracks or nicks. Damaged desiccators must not be used for safety reasons.
- Never expose desiccators to abrupt pressure changes (do not suddenly ventilate evacuated vessels).

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
Vacuum-tight, made from DURAN® borosilicate glass 3.3. To accelerate drying, a vacuum can be applied via the stopcock. Spare parts such as lids, bases, stopcocks and caps can be interchanged (observe DN).

Typical applications: drying of moist samples and storage of moisture-sensitive substances.

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DURAN® VACUUM DESICCATOR with NOVUS standard ground joint (NS 24/29) junction tube in the lid, stopcock and flat flange

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
Made from DURAN® borosilicate glass 3.3. Components are vacuum tight (DIN 12491). Spare parts such as lids and bases can be interchanged (observe DN).

Typical applications: drying of moist products and storage of moisture-sensitive substances.

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DURAN® DESICCATOR with flat flange and knobbed lid, no connection

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DURAN® VACUUM DESICCATOR  
flange with groove and NOVUS standard  
ground joint (NS 24/29) junction tube in the  
lid, with O-ring and stopcock

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DURAN® VACUUM DESICCATOR  
base with NOVUS standard ground joint  
(NS 24/29) and stopcock, flat flange and  
knobbed lid

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DURAN® VACUUM DESICCATOR  
base with flat flange, lid with VERTEX  
outlet and safety stopper

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DURAN® VACUUM DESICCATOR  
with MOBILEX screw thread (GL 32) in  
lid and base, with stopcock

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
**Note:** The desiccators on this page are  
no longer available as sets, but can be  
assembled by selecting components from  
the accompanying table.

<sup>1</sup> The characters XX in the 6th and 7th position should  
be replaced, dependent on DN size: DN 100 = 46,  
DN 150 = 57, DN 200 = 61, DN 250 = 66,  
DN 300 = 69

Best.-Nr. Cat. No.	DN	h mm	h <sub>1</sub> mm	ID Flansch ID Flange mm	AD Flansch OD Flange mm	ca. Volumen approx. Vol. l	Verpackungseinheit/ Stück Pack/Quantity
Klar   Clear							
24 770 46	100	112	58	119	153 +/- 2	0,7	1
24 770 57	150	154	81	172	215 +/- 2	2,4	1
24 770 61	200	202	115	224	270 +/- 2	5,8	1
24 770 66	250	235	120	274	320 +/- 2	10,5	1
24 770 69	300	283	150	332	380 +/- 2	18,5	1

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
DURAN® DESICCATOR BASE  
with flat flange, no outlet,  
suitable for all types of lids

Best.-Nr. Cat. No.	DN	h mm	h <sub>1</sub> mm	ID Flansch ID Flange mm	AD Flansch OD Flange mm	ca. Volumen approx. Vol. l	Verpackungseinheit/ Stück Pack/Quantity
Klar   Clear							
24 773 61	200	202	115	224	270 +/- 2	5,8	1

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
DURAN® DESICCATOR BASE  
with ring-grooved flange,  
suitable for all types of lids

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
DURAN® DESICCATOR BASE  
with flat flange, standard ground outlet (24/29),  
type NOVUS, suitable for all types of lids

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DURAN® DESICCATOR BASE  
with flat flange, screw thread outlet,  
type MOBILEX (GL 32), suitable for all  
types of lids

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

DURAN® DESICCATOR LID  
with knob, suitable for all types of bases

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
DURAN® DESICCATOR LID  
for standard ground joint stopcocks (NS 24/29),  
type NOVUS, suitable for all types of bases

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
DURAN® DESICCATOR LID  
with special tube (NS 24/29 type WERTEX),  
with ring grooved flange, suitable for all types  
of bases

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DURAN® DESICCATOR BASE  
with threaded outlet, type MOBILEX (GL 32),  
suitable for all bases

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


STAINLESS-STEEL DESICCATOR  
PLATE  
Material: ST 02 Z 275 NA, rust-free

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
O-RING<sup>1</sup>  
suitable for desiccators, from silicone (VMQ)

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DURAN® STOPCOCK WITH  
PTFE<sup>1</sup> SPINDLE  
for desiccator base side outlets, type NOVUS  
(NS 24/29)

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Ordering advice: the O-ring is dimensioned  
according to the groove diameter.

Best.-Nr. Cat. No.	DN	d mm	s mm	Verpackungseinheit/ Stück Pack/Quantity
29 224 57	150	174	5,7	I
29 224 60	200	229	5,7	I
29 224 66	250	274	6,5	I
für Modelle ab 1996 for articles since 1996				
29 215 57	150	185	5,3	I
29 215 60	200	236	5,3	I
29 215 66	250	290	5,3	I

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
<sup>2</sup> With self-sealing valve stopcock and PP tightening  
torque limiter.

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DURAN® STOPCOCK WITH  
PTFE<sup>1</sup> SPINDLE  
for desiccator lid outlets, type NOVUS (NS 24/29)

Best.-Nr. Cat. No.	l mm	d mm	NS	Verpackungseinheit/ Stück Pack/Quantity
24 799 04 <sup>2</sup>	85	8	24/29	I

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

DURAN® STOPCOCK WITH  
PTFE<sup>1</sup> SPINDLE  
for threaded outlets, type MOBILEX (GL 32)

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
SCREW CAP WITH APERTURE GL 32  
from PBT<sup>1</sup>, red, with pierced aperture for  
MOBILEX type threaded outlets

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SCREW CAP GL 32  
from PBT<sup>1</sup>, red, with PTFE<sup>1</sup> coated seal for  
MOBILEX type threaded outlets

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
Plastic. Cap for blanking unused outlets  
(type MOBILEX).

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SAFETY STOPCOCK  
for safety outlets type Wertex, NS 24/29

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*"GROWTH ... DEPENDS ON MANY FACTORS,  
NOT JUST FOR MY CULTURES ... IT'S GOOD TO  
KNOW, THAT IN A CHANGING WORLD, DURAN®  
IS ONE **RELIABLE CONSTANT**."*

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4 | GLASSWARE FOR MICROBIOLOGY

Due to its high thermal-shock resistance, DURAN® microbiology glassware is ideal for autoclaving and sterilisation processes and shows, even after multiple use no signs of wear. Unlike plastic items, it is very resistant to mechanical wear even after repeated use and sterilisation cycles.

Due to the nearly inert behaviour, there are no interactions (e.g. ion exchange) between medium and glass and any spurious influence on experiments is thereby effectively excluded.

DURAN® products are completely transparent in visible light and unlike many plastic products are ideal for use under the microscope.

DUROPLAN® Petri dishes are outstanding due to their distortion-free transparency and high planarity. These

excellent geometrical properties enable uniform agar distribution and reproducible culture growth.

Alongside the Petri dishes, the DURAN® range includes a wide range of culture bottles, culture flasks, roller bottles and spot plates. In addition, there are various types of staining dishes.

Usage tips:

- Only autoclave products, which are free from damage such as scratches, cracks or nicks.
- The outstanding thermal properties (max. operating temperature of +500 °C, thermal shock resistance ΔT=100 K) enable high temperature processes, such as hot air sterilisation.

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These Petri dishes are made from DURAN® borosilicate glass 3.3 using a special manufacturing process, which permits the uniform distribution of agar and guarantee distortion-free viewing.

Typical applications: biological and medical work, cultivation of microorganism, microscopy.


Best.-Nr. Cat. No.	d <sub>1</sub> x h <sub>1</sub> mm	d <sub>2</sub> mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
21 755 41	60 x 20	54	22	10
21 755 43	80 x 20	74	22	10
21 755 46	100 x 15	94	17	10
21 755 48	100 x 20	94	22	10
21 755 51	120 x 20	114	22	10
21 755 53	150 x 27	143	32	10

Best.-Nr. Cat. No.	d x h mm	Verpackungseinheit/Stück Pack/Quantity
23 755 39	40 x 12	10
23 755 40	60 x 15	10
23 755 42	80 x 15	10
11 840 71	90 x 15	10
23 755 45	100 x 10	10
23 755 46	100 x 15	10
23 755 48	100 x 20	10
23 755 51	120 x 20	10
23 755 52	150 x 25	10
23 755 56	180 x 30	10
23 755 59	200 x 30	10
23 755 61	200 x 45	10

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
DUROPLAN® PETRI DISH  
base and lid are flat inside and out, and free  
from bubbles and streaks

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
STERIPLAN® PETRI DISH  
MADE OF SODA-LIME GLASS

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
SPOT PLATE, TYPE FEIGL  
MADE OF SODA-LIME GLASS

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Typical application: detection reactions.

Best.-Nr. Cat. No.	l mm	b mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
23 671 52	130	100	14	10

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Bottom of the base has three rests for  
stability.

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Typical applications: biological and  
medical work, preparation of agars,  
microscopy.

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Large, flat bottom surface allows  
uniform culture thickness.

Typical application: preparation  
of cultures in nutrient media.

Best.-Nr. Cat. No.	Inhalt Capacity ml	l mm	b mm	h mm	d mm	Verpackungseinheit/ Stück Pack/Quantity
21 551 71	4 000	370	200	90	50	1

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

DURAN®, CULTURE FLASK,  
ROUX TYPE  
conical neck, excentric

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

DURAN®, CULTURE BOTTLE  
straight rim for Kapsenberg caps

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
Conical geometry makes the flasks particularly suited for shaking experiments (e.g. media optimisation).

Typical application: preparation of cultures in nutrient media.

Best.-Nr. Cat. No.	für Hals d <sub>i</sub> for neck d <sub>i</sub> mm	passend zu suitable for	Verpackungseinheit/ Stück Pack/Quantity
29 010 09	16	Kulturröhrchen Culture tubes 26 132 21	10
29 010 11	18	Kulturröhrchen Culture tubes 26 132 23	10
		Kulturflaschen Culture bottle 21 421 XX	10
		Erlenmeyerkolben Erlenmeyer tubes 21 491 24	10

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

DURAN® CULTURE MEDIA BOTTLE  
straight rim, for glass caps

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

DURAN® SQUARE BOTTLE  
after Breed-Demeter

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DURAN® CULTURE MEDIA BOTTLE  
with beaded rim

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
DURAN® CULTURE MEDIA BOTTLE  
straight neck for metal caps

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Conical geometry makes the flasks  
particularly suitable for shaking experiments  
(e.g. media optimisation). Erlenmeyer flasks  
with screw thread are also available. These  
flasks can be closed with a PBT cap<sup>1</sup> or  
membrane cap<sup>2</sup> (permits gas exchange).

Typical application: preparation  
of cultures in nutrient media.

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Large, flat, bottom surface allow  
uniform culture thickness.


Typical application: preparation  
of cultures in nutrient media.

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
DURAN® CULTURE FLASK,  
FERNBACH TYPE  
conical shape, straight neck for metal caps

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DURAN® CULTURE FLASK,  
ERLENMEYER SHAPE  
straight neck for metal caps

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
Suitable for: culture flask no. 21 771 XX,  
21 774 62 and culture media bottle  
Cat. No. 21 773 XX.


Best.-Nr. Cat. No.	Material Material	Für Hals For neck Ø mm	Verpackungseinheit/ Stück Pack/Quantity
29 012 24	Edelstahl Stainless steel	38	10
29 013 24	Aluminium, blau eloxiert Aluminium, anodised blue	38	10

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
DURAN® ROLLER BOTTLE  
FOR CELL CULTURES  
with DIN thread, GL 45

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
DURAN® ROLLED FLANGE BOTTLE  
without closure

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DURAN® ROLLED FLANGE BOTTLE  
with clamp closure

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With screw cap and pouring ring  
(blue, PP<sup>I</sup>). Spare caps and pouring  
rings see pages 29–31.

Typical application: preparation  
of cultures in nutrient media.

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In addition, the following individual parts  
are available: porcelain clamp closure  
(Cat. No. 29 701 08), replacement rubber  
seal (Cat. No. 29 990 31) and replace-  
ment silicone seal, autoclavable (Cat. No.  
29 990 10).

Typical applications: sampling and cultivation.

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Best.-Nr. Cat. No.	AD OD mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
23 319 00	66	108	10

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
STAINING JAR, COPLIN TYPE  
FROM SODA-LIME GLASS<sup>1</sup>

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Best.-Nr. Cat. No.	l x b x h mm	Verpackungseinheit/Stück Pack/Quantity
23 314 00	90 x 40 x 90	10

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
STAINING DISH,  
HELLENDAAHL TYPE  
FROM SODA-LIME GLASS<sup>1</sup>

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
For 16 microscope slides 76 x 26 mm,  
with widening towards the top.

<sup>1</sup> Hydrolytic class 5: DIN-ISO 719,  
acid resistance class 3: DIN 12 116,  
alkali resistance class 2: DIN ISO 695.

**Note:** Do not clean staining dishes and  
staining jars at temperatures above 60 °C  
(glass corrosion is possible).

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
STAINING DISH,  
SCHIEFFERDECKER TYPE  
FROM SODA-LIME GLASS<sup>1</sup>

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
GLASS BOX  
FROM SODA-LIME GLASS<sup>1</sup>  
for staining tray 21 317 00

Best.-Nr. Cat. No.	l x b x h mm	Verpackungseinheit/Stück Pack/Quantity
23 318 00	108 x 90 x 70	10

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

For 10 microscope slides 76 x 26 mm  
or each width up to 52 mm.

Best.-Nr. Cat. No.	l x b x h mm	Verpackungseinheit/Stück Pack/Quantity
21 317 00	88 x 40 x 70	10

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ACCESSORIES:  
STAINLESS STEEL HANDLE  
for staining tray

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“WHERE WORK'S CONCERNED, I TAKE THINGS  
PRETTY SERIOUSLY ... **EXACT RESULTS**,  
EVERYTHING UNDER CONTROL ... BUT THERE'S  
AN EVEN GREATER PERFECTIONIST IN MY TEAM –  
MY MEASURING INSTRUMENT.”

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5 | VOLUMETRIC GLASSWARE

DURAN® volumetric products have closely calibrated scales that permit very accurate determination and measurement of volumes. DURAN® products are available in two accuracy classes: class A/AS and class B (see the corresponding product descriptions on page 226). The two classes differ in the accuracy of measurement with class A being the highest accuracy class and class B being approximately half that of class A. Class AS has the same tolerances as class A, but is designed to permit more rapid outflow.

Volumetric flasks and cylinders are calibrated to measure the amount of fluid contained (“In”). Thus, for example the desired concentration can be precisely set.

Pipettes and burettes are calibrated to measure the amount of fluid delivered (“Ex”). This calibration takes into account surface adhesion to the glass / capillary effects; specified waiting times must be observed.

Volumetric flasks, volumetric and mixing cylinders as well as burettes are manufactured from DURAN® borosilicate glass 3.3 and have excellent chemical and thermal resistance. Measurement and bulb pipettes are made from soda-lime glass (see page 219 for more information on soda-lime glass).

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Calibration is based on the poured in volume (“In”) at a +20 °C reference temperature. The volume content tolerances conform to accuracy class A, the accuracy limits of the German weights and measures regulations and DIN and ISO specifications.

Typical applications: precise measurement of specified liquid amounts, preparation and storage of standard solutions.

<sup>1</sup> For chemical and thermal resistance, see page 208.  
<sup>2</sup> Non-ISO size.  
<sup>3</sup> For suitable glass stoppers, see pages 44–45.

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
The large hexagonal base prevents the cylinder from rolling. The base is equipped with three knobs which increase its stability. The cylinders have uniform wall thickness over the entire measurement range, so wedge errors are avoided. Calibration is based on the poured in volume (“In”) at a +20 °C reference temperature. Mixing cylinder accuracy limits conform to DIN and ISO standards.

Typical applications: diluting solutions, mixing several components with specified proportions.

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
DURAN® VOLUMETRIC FLASK  
with circular graduation mark and plastic stopper  
from PE<sup>1</sup>, finely ground stopperbed, suitable for  
standard ground stoppers

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
DURAN® MIXING CYLINDERS  
WITH HEXAGONAL BASE  
with graduation, standard ground joint  
and plastic PE<sup>1</sup> stopper

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DURAN® MEASURING CYLINDERS  
WITH HEXAGONAL BASE  
with graduation

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
The large hexagonal base prevents the cylinder from rolling. The base is equipped with three knobs which increase its stability. The cylinders have uniform wall thickness over the entire measurement range, so wedge errors are avoided. Calibration is based on the poured in volume („In“) at a +20 °C reference temperature. Mixing cylinder accuracy limits conform to DIN and ISO standards (class B).

Typical applications: holding and simultaneous measurement of varying liquid amounts.

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DURAN® MEASURING CYLINDER,  
GRADUATED LOW FORM  
with hexagonal base

Best.-Nr. Cat. No.	Inhalt Capacity ml	Fehlergrenze Accuracy limits ± ml	Teilung Graduation divisions ml	h mm	d mm	Verpackungseinheit/ Stück Pack/Quantity
21 395 08	10	0,2	1	90	21	10
21 395 14	25	0,5	1	115	25	10
21 395 17	50	1	2	145	29	10
21 395 24	100	1	2	165	39	10
21 395 36	250	2	5	195	54	10
21 395 44	500	5	10	250	65	10
21 395 54	1 000	10	20	285	85	1
21 395 63	2 000	20	50	340	105	1

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Calibration is based on the poured out volume („Ex“) at a +20 °C reference temperature. Volume content tolerances conform to DIN and ISO. The Class B accuracy limit is roughly one and a half times wider than for Class AS. The tolerances are thus more strict than specified by DIN.

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
With Schellbach stripe and main graduations as circular divisions. Calibration is based on the poured out volume („Ex“) at a +20 °C reference temperature. Volume content tolerances conform to DIN.

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
DURAN® BURETTE  
with straight standard ground stopcock, class B

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
DURAN® BURETTE  
with straight standard ground stopcock,  
class AS, 30 seconds waiting time

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
DURAN® MICRO-BURETTE  
with straight standard ground stopcock,  
class AS, 30 seconds waiting time

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DURAN® AUTOMATIC BURETTE  
Pellet-type, side-positioned standard ground  
stopcock, class AS, 30 seconds waiting time

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With Schellbach stripe and main graduations  
as circular divisions. Calibration is based on  
the poured out volume (“Ex”) at a +20 °C  
reference temperature. Volume content  
tolerances conform to DIN.

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With Schellbach stripe and main graduations  
as circular divisions, reservoir bottle (2,000 ml)  
and blowball.

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Replacement bottle for automatic burettes.

Best.-Nr. Cat. No. klar   clear	Best.-Nr. Cat. No. braun   amber	Inhalt Capacity ml	h mm	d mm	Verpackungseinheit/ Stück Pack/Quantity
Hals ungeschliffen Neck unground					
21 150 63 <sup>1</sup>	21 150 63 6	2 000	200	160	1
Mit Normalschliff NS 29/32 With standard ground NS 29/32					
21 159 63	21 159 63 6	2 000	200	160	1

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Numbering from the top down. Calibration  
is based on the poured out volume (“Ex”) at  
a +20 °C reference temperature. Due  
to the scale, variable volumes can be held  
and then dispensed in the same or differing  
increments.

Typical applications: accurate measurement  
and decanting of liquids.


Best.-Nr. Cat. No.	Inhalt Capacity ml	Fehlergrenze Accuracy limits ± ml	Teilung Graduation divisions ml	Farbkenn- zeichnung Colour code DIN 12 621	Gesamtlänge Overall length mm	Verpackungs- einheit/Stück Pack/Quantity
24 343 01 <sup>1</sup>	0,1	0,02	0,001	3 x grün   green	360	12
24 343 03 <sup>1</sup>	0,2	0,0045	0,001	3 x blau   blue	360	12
24 343 06	0,5	0,008	0,01	3 x gelb   yellow	360	12
24 343 11	1	0,008	0,01	2 x gelb   yellow	360	12
24 343 16	2	0,015	0,01	3 x weiß   white	360	12
24 343 23	5	0,04	0,05	2 x rot   red	360	12
24 343 29	10	0,08	0,1	2 x orange   orange	360	12
24 343 34	25	0,15	0,1	2 x weiß   white	450	12

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<sup>1</sup> Non-ISO size, calibrated to contain (“In”).

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

MEASURING PIPETTE  
FOR PARTIAL OUTFLOW  
FROM SODA-LIME GLASS  
class B, graduated

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MEASURING PIPETTE  
FOR COMPLETE OUTFLOW  
FROM SODA-LIME GLASS  
class B, graduated

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
Numbering from the top down. Calibration is based on the poured out volume („Ex“) at a + 20 °C reference temperature. Due to the scale, variable volumes can be held and then dispensed in the same or differing increments.

Typical applications: accurate measurement and decanting of liquids.

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MEASURING PIPETTE FOR  
COMPLETE AND PARTIAL OUTFLOW  
FROM SODA-LIME GLASS  
class AS, main graduations as circular divisions,  
15 seconds waiting time

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Numbering from the top down. Due to the scale, variable volumes can be held and then dispensed in the same or differing increments.

Typical applications: accurate measurement and decanting of liquids.

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Calibration is based on the poured out volume („Ex“) at a + 20 °C reference temperature. Calibrated to measure and discharge a single volume.


Typical applications: accurate measurement and decanting of liquids.

Calibrated to measure and discharge a single volume („Ex“) at a + 20 °C reference temperature. Calibrated to measure and discharge a single volume.

Typical applications: accurate measurement and decanting of liquids.

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
BULB PIPETTE  
FROM SODA-LIME GLASS  
class B

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BULB PIPETTE  
FROM SODA-LIME GLASS  
class AS, 15 seconds waiting time

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“MY LONG-TERM PROJECT ... TO SEPARATE  
THE IMPORTANT FROM THE UNIMPORTANT ...  
IF ONLY EVERYTHING WAS AS **EASY** AS  
SEPARATING USING DURAN®.”

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6 | GLASS FILTRATION APPARATUS AND ACCESSORIES

Due to their high chemical and thermal shock resistance DURAN® filters and the corresponding filter plates are ideal for separations, e. g. with strong acids or alkalis. Thus they offer advantages in comparison with other materials such as plastic or paper. DURAN® filter products have a maximum operating temperature of +450 °C and are therefore far superior to other materials.

The corresponding filtration vessels are specially optimised to the matching filtration apparatus and are vacuum-tight due to their special geometry and high wall thickness. This characteristic has been approved by the German TÜV and marked with the "GS"-indication.

The glass filters are classified as being in porosity classes 0 to 5 according to their nominal maximum pore size. The following table shows the corresponding porosity range. The specified pore sizes always relate to the largest pore in the plate. This specification also characterises the minimum nominal size of particles which may be retained by the filtration.

Usage tips:

- The maximum permissible operating temperature is +450 °C.

- Uniform heating is recommended to avoid thermal stresses and resultant breakages.
- Heat glass filtration apparatus with disk diameters of more than 20 mm in initially cold ovens or sterilisers only.
- The heating or cooling rate should not exceed 8 °C/min.
- When filtering hot substances observe the thermal shock resistance and, if necessary, preheat the filtration apparatus in a drying cabinet.
- Wet filtration apparatus should be heated slowly up to +80 °C and dried for one hour before increasing the temperature further.

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
Particular suited to use with hot or aggressive media thanks to the thermal shock and chemical resistance of DURAN®. Also available in **soda-lime glass**.

Typical applications: decanting and filtration of substances.

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
From DURAN® glass with its good thermal shock and chemical resistance. Also available in **soda-lime glass**.

Typical applications: decanting of powdered substances and granulated material.

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
DURAN® POWDER FUNNEL  
with short, wide stem

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
DURAN® ANALYTICAL FUNNEL  
for quick filtration

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From DURAN® glass with its good thermal shock and chemical resistance. Also available in **soda-lime glass**.

Typical applications: filtering and decanting of liquids of different densities.

<sup>1</sup> Non-DIN ISO size.  
<sup>2</sup> Details apply only to DURAN® version.

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From DURAN® glass with its good thermal shock and chemical resistance.

Typical application: for rapid liquid filtration.

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
From DURAN® glass with its good thermal shock and chemical resistance. The ribbed form is ideal for filtering with round-paper.

Typical application: filtering of liquids.

Made from DURAN® glass with its good thermal shock and chemical resistance. Filter funnels mate to the filtering flask via a conical rubber seal (GUKO).


Typical application: filtering of liquids.

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

From DURAN® glass with its good thermal shock and chemical resistance. Filter funnels mate to the filtering flask via a conical rubber seal (GUKO).

Typical applications: qualitative inorganic analysis and preparative chemistry.

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
From DURAN® glass with its good thermal shock and chemical resistance. The Buechner funnel features a glass support for membrane and paper filters.

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<sup>1</sup> Substitute required porosity (see page 222) for X in the Cat No.


Best.-Nr. Cat. No.	d mm	passender Filtrivorst- stoß Best.-Nr. suitable filter adapter Cat. No.	Verpackungseinheit/Stück Pack/Quantity
29 201 14	26	24 316 16	10
29 201 21	33	24 316 22	10
29 201 26	41	24 316 26	10
29 201 31	49	24 316 32	10

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
RUBBER SLEEVE  
from EPDM<sup>1</sup> for filter crucibles

Best.-Nr. Cat. No.	d mm	d <sub>i</sub> mm	h mm	passende Gummi- manschette Best.-Nr. suitable rubber sleeve Cat. No.	Verpackungseinheit/ Stück Pack/Quantity
24 316 16	27	10	108	29 201 14	10
24 316 22	34	10	110	29 201 21	10
24 316 26	41	10	125	29 201 26	10
24 316 32	50	10	132	29 201 31	10

DURAN® FILTRIERVORSTOSS


DURAN® FILTER CRUCIBLE /  
FUNNEL ADAPTER

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DURAN® FILTER FUNNEL HEAD  
WITH PP<sup>1</sup> FUNNEL AND TWO  
FKM<sup>1</sup> SEALS

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
Interchangeable filter disks. Available in three filter diameters, each of which is available in four different porosities. Important: insert filter disc between two FKM seals. After filtration, disk can be removed to allow simple and safe removal of the filtrand. Long filter disk service life, as disks are not damaged when the filtrand is scraped off. Easy cleaning of both sides is possible. Cost-effective as components and disks can be ordered individually as required.

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
DURAN® FILTER DISK  
WITH GLASS RIM  
patented

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From DURAN® glass with its good thermal shock and chemical resistance. Fused glass rim.

<sup>1</sup> Substitute required porosity (see page 222) for X in the Cat. No.

Best.-Nr. Cat. No.	für Platte for disc d mm	Verpackungseinheit/ Stück Pack/Quantity
29 220 24	24	10
29 220 50	50	10
29 220 90	90	10

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


FUNNEL FOR FILTER FUNNEL HEAD  
from PP<sup>1</sup>

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Typical application: in-line filtration of gas lines to remove solid impurities (e.g. dust).


Best.-Nr. Cat. No.	Porosität Porosity	Platte Disc d mm	d <sub>i</sub> mm	Verpackungseinheit/Stück Pack/Quantity
25 855 0X <sup>1</sup>	I-4	30	14	I
25 855 1X <sup>1</sup>	I-4	60	16	I
25 855 2X <sup>1</sup>	I-4	90	16	I

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DURAN® IMMERSION FILTER  
for reverse filtration

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Typical application: extraction of clear filtrate (not the filtrand).


Best.-Nr. Cat. No.	Porosität Porosity	Platte Disc d mm	d <sub>i</sub> mm	h mm	Verpackungseinheit/Stück Pack/Quantity
25 855 6X <sup>1</sup>	I-4	35	10	210	I

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Typical application: reactions between gases and liquids.

Best.-Nr. Cat. No.	Porosität Porosity	d mm	d <sub>i</sub> mm	h mm	Verpackungseinheit/Stück Pack/Quantity
25 856 0X <sup>1</sup>	0-2	6	22	250	5
25 856 1X <sup>1</sup>	0-2	9	25	250	5
25 856 2X <sup>1</sup>	I	11	34	250	5

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With screw-connection system. The insertion height of the head is adjustable. Individual parts can also be ordered separately (see page 128).

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Typical applications: cleaning (“washing”) of gases with solvents.

Best.-Nr. Cat. No.	Inhalt Capacity ml	DIN-Gewinde DIN Thread GL	Olive Hose connection d mm	Verpackungseinheit/ Stück Pack/Quantity
24 713 00	500	45	9	I

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With screw-connection system. The insertion height of the head is adjustable. The distribution of the gas in the liquid is carried out through the filter disk and thus the absorption efficiency significantly increased. Gas washing bottles also work reliably at high flow rates. Individual parts can also be ordered separately (see page 128).

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

DURAN® GAS WASHING BOTTLE  
head Drechsel-type, without filter disk

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DURAN® GAS WASHING BOTTLE  
head Drechsel-type, with filter disk

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
DURAN® GAS WASHING BOTTLE  
head with filter disk, with standard ground joint

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
Typical application: cleaning (“washing”) of gases with solvents.

Best.-Nr. Cat. No.	Porosität Porosity	Inhalt Capacity ml	NS Joint	d mm	h mm	h <sub>1</sub> mm	Olive Hose connection d mm	Verpackungs- einheit/Stück Pack/Quantity
25 703 01	I	350	29/32	60	250	180	10	I

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DURAN® GAS WASHING BOTTLE  
with fused-in filter disk,  
with standard ground joint and cap

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
DURAN® MICRO IMMERSION FILTER  
for reverse filtration


Best.-Nr. Cat. No.	Porosität Porosity	d mm	h mm	Verpackungseinheit/Stück Pack/Quantity
25 857 2X <sup>1</sup>	I-4	13	25	10

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<sup>1</sup> Substitute required porosity (see page 222) for X in the Cat. No.

DURAN® MIKRO-FILTERKERZE  
ohne Rohr


DURAN® MICRO FILTER CANDLE  
without tube

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DURAN® MICRO FILTER CANDLE  
with narrow tube

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

DURAN® MICRO FILTER FUNNEL  
Pregl type

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

Heavy walled for vacuum use. These filtering flasks fulfil the regulations of the "equipment and product safety regulations". For suitable rubber gaskets, see page 135.

Typical applications: separations by vacuum filtration.

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	d <sub>l</sub> mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
2I 183 36	250	85	34	155	10
2I 183 44	500	105	34	185	10
2I 183 54	1 000	135	45	230	10
2I 183 63	2 000	166	60	255	1

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DURAN® FILTERING FLASK  
WITH SIDE-ARM SOCKET  
Erlenmeyer shape

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DURAN® FILTERING FLASK  
WITH SIDE-ARM SOCKET  
bottle shape

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*Heavy walled for vacuum use. These filtering flasks fulfil the regulations of the "equipment and product safety regulations". For suitable rubber gaskets, see page 135.*


*Typical application: separations by vacuum filtration.*

**Note:** These filtering flasks have a ground socket 17.5/26 for vacuum tube of 15 to 18 mm OD (e.g. 6 x 5 mm or 8 x 5 mm, DIN 12 865).

*Provision of filtration flasks with a socket has not only made work in preparation and analytical laboratories easier and simpler, but has also reduced the risk of accidents.*

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DURAN® FILTERING FLASK  
WITH GLASS HOSE CONNECTION  
Erlenmeyer shape

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
*Due to the heavy wall thickness the apparatus is vacuum-tight. For suitable rubber gaskets, see page 135.*

*Typical application: separations by vacuum filtration.*

<sup>1</sup> Does not conform to the "equipment and product safety regulations".

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DURAN® FILTERING FLASK  
WITH GLASS HOSE CONNECTION  
bottle shape

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*Heavy walled for vacuum use. For suitable rubber gaskets, see page 135.*

*Typical application: separations by vacuum filtration.*

<sup>1</sup> Does not conform to the "equipment and product safety regulations".

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
*Heavy walled for vacuum use. These filtering flasks fulfil the regulations of the "equipment and product safety regulations". The plastic hose connection is replaceable. For suitable rubber seals, see page 135.*

*Typical application: Separations by vacuum filtration.*

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

DURAN® FILTERING FLASK  
WITH KECK™ ASSEMBLY SET  
bottle shape

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Due to the heavy wall thickness the apparatus is for vacuum use. These filtering flasks fulfil the regulations of the "equipment and product safety regulations". The plastic hose connections can be replaced. For suitable rubber gaskets, see page 135.

Typical application: separations by vacuum filtration.

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
With removable plastic hose connection (PBT<sup>1</sup>), screw (PP<sup>1</sup>), seals (VMQ<sup>1</sup>, EPDM<sup>1</sup>). Suitable for filtering flasks 100 - 20 000 ml.

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EPDM<sup>1</sup> GUKO SET, FROM EPDM<sup>1</sup>  
conical rubber gaskets for filtering flasks

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
GUKO FROM EPDM<sup>1</sup>  
conical rubber gaskets for filtering flasks

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
DURAN® LID FOR FILTER  
APPARATUS, WITT TYPE  
with tubulature, standard ground joint 29/32

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<sup>1</sup> Comprises the following order numbers:  
21 204 445, 29 202 27, 24 316 26, 29 201 26,  
25 851 23 and 25 851 24.

Best.-Nr. Cat. No.	DN	Hülse Socket size NS	Verpackungseinheit/ Stück Pack/Quantity
24 398 46	100	29/32	I
24 398 57	150	29/32	I
24 398 61	200	29/32	I

Best.-Nr. Cat. No.	DN	h <sub>1</sub> mm	h <sub>2</sub> mm	Hülse Socket size NS	Verpackungseinheit/ Stück Pack/Quantity
24 730 46	100	160	76	29/32	I
24 730 57	150	200	102	29/32	I
24 730 61	200	300	126	29/32	I

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
DURAN® FILTER APPARATUS,  
WITT TYPE  
complete with interchangeable lid and KECK™  
assembly set<sup>1</sup>, standard ground joint 29/32

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
DURAN® FILTER APPARATUS,  
WITT TYPE  
complete with interchangeable lid and KECK™  
assembly set<sup>1</sup>, standard ground joint 45/40

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
DURAN® LID FOR FILTER APPARATUS,  
WITT TYPE  
wide rough-ground tubulature, standard  
ground joint 45/40

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DURAN® WATER JET PUMP  
with non-return valve

Best.-Nr. Cat. No.	DN	Hülse Socket size NS	Verpackungseinheit/Stück Pack/Quantity
24 450 46	100	45/40	1
24 450 57	150	45/40	1

Best.-Nr. Cat. No.	Wasserverbrauch mind. Water flow rate min. l/h	Wasserdruck mind. Water pressure min. bar	h mm	Verpackungseinheit/ Stück Pack/Quantity
24 362 99	300	1,2	275	1

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Throughput: 400 l/h at 3.5 bar water  
pressure und 12 °C water temperature.


Best.-Nr. Cat. No.	Olive passend für Schlauch ID Connection suitable for hose ID mm	Wasserverbrauch mind. Water flow rate min. l/h	Wasserverbrauch mind. Water flow rate min. bar	h mm	Verpackungs- einheit/Stück Pack/Quantity
29 250 01	9-12	170	1	235	1

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WATER JET PUMP  
from plastic (PP<sup>1</sup>), with non-return valve,  
hose connection and adapters for 1/2" and 3/4"

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*“FLEXIBLE, ACCURATE, HIGHLY VERSATILE ...  
IT LOOKS AS IF I SHOULD FOLLOW THE EXAMPLE  
OF MY LABORATORY GLASSWARE.”*

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7 | INTERCHANGEABLE GLASSWARE

DURAN® interchangeable glassware is indispensable for laboratory work. DURAN GROUP offers a wide range of bottles and flasks with standard ground necks, vessels with flat flanges, condensers and stirrer shafts.

The DURAN® flat flange reaction vessels are valued for their universal suitability for use in the laboratories of a wide range of specialisations.

Whether reaction, distillation, evaporation or desiccation, DURAN® offers a wide range of unfinished and finished parts which always provide the optimum solution for the particular application. Due to the pure glass-glass connections, reactions with highly corrosive or highly chemically reactive substances can be carried out without problems.

The vessels are notable due to a robust glass flange design with an optimum external flange angle of 45°. Due to the precisely ground joint, the vessels can be closed tightly when using a sealing ring.

Matching stainless-steel quick release clamps, with three retaining clips, ensure easy and safe handling.

All individual parts and a wide range of accessories such as lids, seals, quick-release clamps etc. are compatible and can be interchanged as required. Vessels and lids can be matched using their DN (nominal diameter) number.

Usage tips:

- All components are suitable for use under full vacuum and approved for operating over-pressures (see product related pages).
- Before use, it is recommended that the glass surfaces of the vessels be checked for damage such as scratches, cracks or nicks.
- Damaged glass vessels should not be used for safety reasons.
- Due to the high wall thickness and reduced thermal shock resistance under pressure loading, the flat flange vessels should be heated uniformly and gradually.

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Pressure and vacuum resistant due to the wall thickness and geometry. Groove for O-ring seal (see pages 146–147 for suitable sealing rings). For suitable KPG stirrer bearings and shafts, see pages 159–161.

Typical applications: reactions under pressure and/or high temperature.

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

DURAN® FLAT FLANGE  
REACTION VESSEL  
flange with groove

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DURAN® FLAT FLANGE ROUND  
BOTTOMED FLASK  
flange with groove, for vacuum use

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
Pressure and vacuum resistant due to the wall thickness and geometry. Groove for O-ring seal (see pages 146–147 for suitable sealing rings).

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
**Note:** only heat flat-flange beakers in water or oil baths. At the maximum usage temperature of 250 °C and the maximum operating pressure, the temperature difference in the glass wall of the flat flange reaction vessels must not exceed 30 °C.

<sup>1</sup> Suitable for Witt-type filter apparatus.

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
Pressure and vacuum resistant due to the wall thickness and geometry. For suitable KPG stirrer bearings and shafts, see pages 159–161.

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**Note:** At the maximum usage temperature of 250 °C and the maximum operating pressure, the temperature difference in the glass wall of the flat flange reaction vessels must not exceed 30 °C.

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DURAN® FLAT FLANGE LID  
with centre neck

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
Pressure and vacuum resistant due to the wall thickness and geometry.

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
Accessories for flat flange vessels, comprising an elastic, silicone core with seamless FEP<sup>1</sup> coating that encloses the ring. The combination of these high-quality materials achieves good elasticity in conjunction with outstanding chemical resistance.

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O-RING, TRANSPARENT  
aus Silikon (VMQ<sup>1</sup>)

O-RING, TRANSPARENT  
from silicone (VMQ<sup>1</sup>)

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Accessories for flat flange vessels. From silicone (VMQ<sup>1</sup>), so highly elastic. The chemical resistance is reduced in comparison with FEP<sup>1</sup> coated O-rings.

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Accessories for flat flange vessels.

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
For secure fitting of the lid or the reaction vessel using two clamping rods.

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
QUICK RELEASE CLAMP  
from stainless steel, with 3 retaining clips

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HOLDING DEVICE FOR  
REACTION VESSELS  
from chrome-nickel steel

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DURAN® WOULFF BOTTLE  
3 standard ground necks

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DURAN® WOULFF BOTTLE  
3 standard ground necks and bottom tubulature

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Vacuum resistant due to the wall thickness and geometry.

Best.-Nr. Cat. No.	Inhalt Capacity ml	Hälse Necks NS	d mm	Verpackungseinheit/ Stück Pack/Quantity
24 709 44	500	19/26	87	I
24 709 54	I 000	24/29	I I3	I
24 709 63	2 000	29/32	I35	I
24 709 73	5 000	34/35	I85	I

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Thanks to the uniform wall thickness, round bottom flasks are ideal as heating vessels. The geometry permits very uniform heating.

Typical applications: distillation, extraction.

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
The conical geometry makes them ideal for small-scale reactions.

Best.-Nr. Cat. No.	Inhalt Capacity ml	Hals Neck NS	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
24 I95 08	I0	I4/23	30	75	I0
24 I95 I4	25	I4/23	40	90	I0
24 I95 20	50	I4/23	50	I I0	I0
24 I95 25	I00	I4/23	62	I25	I0

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
DURAN® ROUND BOTTOM FLASK  
with standard ground joint


DIN ISO 4797  
A 121 °C  
USP Standard

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DURAN® PEAR SHAPE FLASK  
with standard ground joint

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
See pages 164–165 for the Keck™ clips for conical joints.

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DURAN® ERLLENMEYER FLASK  
with standard ground joint

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With easy-to-read scale and large labelling field for easy marking. The conical shape makes these flasks ideal for mixing liquids and, due to the even wall thickness, also suitable for use as heating glassware.


Best.-Nr. Cat. No.	Inhalt Capacity ml	Hals Neck NS	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
24 193 13	25	14/23	42	75	10
24 193 20	50	14/23	51	85	10
24 193 17	50	19/26	51	85	10
24 193 18 <sup>1</sup>	50	24/29	51	85	10
24 193 19	50	29/32	51	85	10
24 193 24	100	19/26	64	105	10
24 193 26 <sup>1</sup>	100	24/29	64	105	10
24 193 27	100	29/32	64	105	10
24 193 32 <sup>1</sup>	200	29/32	79	131	10
24 193 36	250	24/29	85	140	10
24 193 37	250	29/32	85	140	10
24 193 38 <sup>1</sup>	250	45/40	85	140	10
24 193 39 <sup>1</sup>	300	29/32	87	155	10
24 193 44	500	24/29	105	175	10
24 193 46	500	29/32	105	175	10
24 193 47 <sup>1</sup>	500	45/40	105	175	10
24 193 54	1 000	24/29	131	220	10
24 193 56	1 000	29/32	131	220	10
24 193 57 <sup>1</sup>	1 000	45/40	131	220	10

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<sup>1</sup> Non-DIN ISO size.

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DURAN® FLAT BOTTOM FLASK  
with standard ground joint

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Due to the flat bottom the flask can be set upon a bench without a support ring.

Best.-Nr. Cat. No.	Inhalt Capacity ml	Hals Neck NS	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
24 171 19	50	29/32	51	85	10
24 171 24	100	19/26	64	103	10
24 171 26	100	24/29	64	103	10
24 171 27	100	29/32	64	103	10
24 171 36 <sup>1</sup>	250	24/29	85	130	10
24 171 37	250	29/32	85	130	10
24 171 44 <sup>1</sup>	500	24/29	105	160	10
24 171 46	500	29/32	105	160	10
24 171 54 <sup>1</sup>	1 000	24/29	131	187	10
24 171 56	1 000	29/32	131	187	10
24 171 63 <sup>1</sup>	2 000	29/32	166	230	10

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

DURAN® EVAPORATING FLASK  
with standard ground joint, pear shape

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DURAN® TWIN-NECK, ROUND  
BOTTOM FLASK  
with standard ground joint, inclined side neck

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
Thanks to the uniform wall thickness, round bottom flasks are ideal as heating vessels. The geometry permits very uniform heating. Depending upon the application, accessories, columns, thermometers, dropping funnels, boiling capillaries, etc. can be fitted.

Typical applications: distillation, extraction.

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DURAN® TRIPLE-NECK, ROUND  
BOTTOM FLASK  
with standard ground joint, inclined side necks

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
DURAN® TRIPLE-NECK, ROUND  
BOTTOM FLASK  
with standard ground joint, parallel side necks

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
Typical application: distillation.

Best.-Nr. Cat. No.	Nutzlänge Effective length mm	Hülse Socket size NS	Kern Cone NS	Gesamtlänge Overall length mm	Verpackungs- einheit/Stück Pack/Quantity
24 240 71	300	24/29	24/29	450	1
24 240 72	300	29/32	29/32	450	1
24 240 87	500	24/29	24/29	650	1
24 240 88	500	29/32	29/32	650	1

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
DURAN® VIGREUX COLUMN  
with 2 standard ground joints, complete with  
slide-on glass jacket

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
DURAN® LIEBIG CONDENSER  
(WEST CONDENSER)  
with 2 standard ground joints and  
2 screw-on plastic hose connections

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DURAN® BULB CONDENSER  
(ALLIHN CONDENSER)  
with 2 standard ground joints and  
2 screw-on plastic hose connections

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Relatively small heat exchange surface  
and thus relatively low cooling capacity.  
  
Typical application: product condenser  
for distillate separations.

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Typical application: product condenser  
for distillate separations.

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
The Dimroth condenser comprises a coil  
condenser located within a tube. This  
condenser type has a larger heat exchange  
surface and thus a better cooling effect  
than the Liebig or Allihn condenser.

Typical applications: product and reflux  
condenser.

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
DURAN® COIL DISTILLATE  
CONDENSER  
with 2 standard ground joints and  
2 screw-on plastic hose connections

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
DURAN® DIMROTH CONDENSER  
with 2 standard ground joints and  
2 screw-on plastic hose connections

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DURAN® JACKETED COIL  
CONDENSER  
with 2 standard ground joints and  
2 screw-on plastic hose connections

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
*This type of condenser has a very large  
heat exchange surface due to its coil  
condenser and double jacket and is thus  
especially suited to working with low-  
boiling point media.*

*Typical application: use as a reflux  
condenser for condensation and feedback  
of the (solvent) vapour to the reaction  
mixture.*

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
*With standard ground stopcock and  
retaining device.*

*Typical applications: uniform and metered  
liquid supply to a reaction mixture. The  
rate of supply can be adjusted.*

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DURAN® DROPPING FUNNEL  
cylindrical, with scale and ground joint

DIN ISO  
4800

*With standard ground joint, standard  
ground stopcock and retaining device.*

*Typical application: uniform and metered  
liquid supply to a reaction mixture. The  
rate of supply can be adjusted.*

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
<sup>1</sup> Special size, non-DIN size.

*With pressure equalisation tube, standard  
ground cone, standard ground stopcock  
and retaining device.*

*Typical application: uniform and metered  
liquid supply to a reaction mixture. The  
rate of supply can be adjusted.*

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
DURAN® DROPPING FUNNEL  
cylindrical, with scale

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DURAN® SEPARATING FUNNEL  
conical shape

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
With standard ground stopcock, retaining device and plastic stopper. The conical shape makes it highly suited to phase separation (See pages 44–45 for replacement stoppers.).

Best.-Nr. Cat. No.	Inhalt Capacity ml	Stiel Stern d <sub>i</sub> mm	h mm	NS-Glasküken Nenngröße (DIN 12 541) Standard solid key nominal size (DIN 12 541)	Hülse Socket size NS	Verpackungs- einheit/Stück Pack/Quantity
Mit Glasküken With solid key						
24 294 17	50	9	190	3 NS	19/26	1
24 294 24	100	9	230	3 NS	19/26	1
24 294 36	250	10	280	4 NS	29/32	1
24 294 44	500	10	320	4 NS	29/32	1
24 294 54	1 000	13	380	6 NS	29/32	1
24 294 63	2 000	13	430	6 NS	29/32	1
Mit PTFE-Küken With PTFE key						
10 648 05	100	9	230	3 NS	19/26	1
10 648 06	250	10	280	4 NS	29/32	1
10 648 07	500	10	320	4 NS	29/32	1
10 648 09	1 000	13	380	6 NS	29/32	1

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DURAN® SEPARATING FUNNEL  
spherical, heavy-duty version

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With standard ground stopcock and standard ground stopper (See pages 44–45 for replacement stoppers.).

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Typical application: Phase separation.

Best.-Nr. Cat. No.	Inhalt Capacity ml	Stopfen NS	h mm	d mm	Verpackungs- einheit/Stück Pack/Quantity
24 291 36	250	24/20	235	90	1
24 291 44	500	24/29	276	115	1
24 291 54	1 000	29/32	295	132	1
24 291 66	2 500	45/40	370	182	1
24 291 73	5 000	45/40	425	222	1
24 291 86	10 000	45/40	490	286	1

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


Ground and polished bearing surface, with standard ground cone.

Best.-Nr. Cat. No.	Kennzeichnung Designation	d mm	h mm	Kern Cone size NS	Verpackungs- einheit/Stück Pack/Quantity
24 523 55	HP 10	10	75	24/29	1
24 528 56	HQ 10	10	75	29/32	1
24 535 49	HS 10	10	65	24/29	1
24 540 51	HT 10	10	65	29/32	1
24 541 54	HT 16	16	85	29/32	1

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

Ground and polished bearing surface.

Best.-Nr. Cat. No.	Kennzeichnung Designation	d mm	h mm	h <sub>1</sub> mm	Verpackungs- einheit/Stück Pack/Quantity
24 500 42	HA 10	10	80	65	1
24 505 51	HB 10	10	120	75	1
24 506 57	HB 16	16	150	90	1

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
Ground and polished bearing surface,  
with GL screw thread tube and standard  
ground cone.

Best.-Nr. Cat. No.	Kennzeichnung Designation	d mm	h mm	DIN-Gewinde DIN thread GL	Kern Cone size NS	Verpackungs- einheit/Stück Pack/Quantity
24 750 08	HB 10	10	75	32	24/29	1
24 750 09	HB 10	10	75	32	29/32	1

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


Bearing surface ground and polished.

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

DURAN® DRYING TUBE, BENT  
with standard ground cone

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DURAN® CONNECTION PIECE  
with standard ground cone, 90° angle

Best.-Nr. Cat. No.	Kern Cone size NS	Verpackungseinheit/Stück Pack/Quantity
24 300 06	14/23	1
24 300 08	24/29	1
24 300 09	29/32	1

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DURAN® RECEIVER ADAPTER  
bent, with ground socket

Best.-Nr. Cat. No.	Hülse Socket size NS	Verpackungseinheit/Stück Pack/Quantity
24 310 06	14/23	1
24 310 08	24/29	1
24 310 09	29/32	1

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See pages 168–170 for screw caps and  
hose connections.


Best.-Nr. Cat. No.	Hülse Socket size NS	Kern Cone size NS	DIN-Gewinde DIN thread GL	Olive Hose connection d mm	Verpackungs- einheit/Stück Pack/Quantity
24 130 21	14/23	14/23	14	8,6	1
24 130 34	24/29	24/29	14	8,6	1
24 130 46	29/32	29/32	14	8,6	1

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DURAN® VACUUM RECEIVER  
ADAPTER

bent, with 2 standard ground joints  
and screw-on plastic hose connection

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

DURAN® VACUUM RECEIVER  
ADAPTER

straight, with 2 standard ground joints  
and screw-on plastic hose connection

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
DURAN® ADAPTER  
with standard ground joint

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KECK™ CLIP  
for conical, joints, from POM<sup>1</sup>

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
KECK™ CLIP ASSORTMENT  
for conical, joints, from POM<sup>1</sup>

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
KECK™ CLIP ASSORTMENT  
for conical, joints from metal<sup>1</sup>

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

KECK™ CLIP ASSORTMENT  
for spherical joints, from POM<sup>1</sup>

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
DURAN® SCREWTHREAD TUBE  
with DIN thread

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For flexible connection of two glass screw-  
thread connections. With integral PTFE<sup>1</sup>  
faced silicone (VMQ<sup>1</sup>) seal.

Best.-Nr. Cat. No.	DIN-Gewinde DIN thread GL	Verpackungseinheit/Stück Pack/Quantity
29 226 05 5	14	1
29 226 06 5	18	1
29 226 09 5	25	1
29 226 08 5	32	1

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High leak tightness through use of PTFE<sup>1</sup>  
coated silicone seal. More chemically  
resistant than PP<sup>1</sup> cap.


Best.-Nr. Cat. No.	DIN-Gewinde DIN thread GL	Verpackungseinheit/Stück Pack/Quantity
29 240 08	14	10
29 240 11	18	10
29 240 13	25	10
29 240 19	32	10
29 240 28	45	10

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DURAN® SCREWTHREAD  
TO GROUND JOINT ADAPTER  
with DIN thread and standard ground cone

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


SCREW CAP WITH APERTURE  
from PBT<sup>1</sup>, red

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Suitable for PBT<sup>1</sup> screw caps, premium and tamper-evident screw caps<sup>2</sup>. Heat resistance: 130 °C (vapour) and 200 °C (dry heat). Good chemical resistance due to PTFE<sup>1</sup> coating.

<sup>1</sup> For chemical and thermal resistance, see page 208.  
<sup>2</sup> For caps, see pages 29–31.

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Suitable for silicone seal for piercing (Septa). More chemically resistant than PP<sup>1</sup> cap.

Typical applications: Injection or removal of media.

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Suitable for PBT<sup>1</sup> screw cap with aperture. Heat resistance: 130 °C (vapour) and 200 °C (dry heat).

Typical application: fixing of tubes.

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
Suitable for PBT<sup>1</sup> screw cap with aperture. Heat resistance: 130 °C (vapour) and 200 °C (dry heat).

Typical applications: injection or removal of media.

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SILICONE SEALING RING (VMQ<sup>1</sup>)  
VMQ<sup>1</sup> with bonded PTFE<sup>1</sup> face

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

A  
121 °C

KUNSTSTOFF-OLIVE  
gerade, aus PP<sup>1</sup>

PLASTIC HOSE CONNECTION  
straight, from PP<sup>1</sup>

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With silicone seal (VMQ<sup>1</sup>). Suitable for  
GL 14 screw cap (Cat. No. 29 227 05).


Best.-Nr. Cat. No.	d mm	d <sub>i</sub> mm	Verpackungseinheit/ Stück Pack/Quantity
29 255 06	8,6	5	10

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KUNSTSTOFF-OLIVE  
gebogen, aus PP<sup>1</sup>

PLASTIC HOSE CONNECTION  
bent, from PP<sup>1</sup>

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
For plastic hose connections (Cat. No.  
29 255 and 29 247 05). From silicone  
(VMQ<sup>1</sup>).

Best.-Nr. Cat. No.	Material Material	Verpackungseinheit/Stück Pack/Quantity
29 220 09	Silikon (VMQ) <sup>1</sup> Silicone (VMQ) <sup>1</sup>	10

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“I USED TO THINK: GLASS IS JUST GLASS...  
TODAY I KNOW THAT YOU CAN'T TELL QUALITY  
AT FIRST GLANCE... BUT IT'S IN THE NAME..”

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WHAT IS GLASS ?

Glass is an inorganic mixture fused at high temperature which solidifies on cooling but does not crystallize. Its basic components, network formers and modifiers, are present in the common glasses in the form of oxides.

Typical glass formers (network formers) are silicon dioxide ( $\text{SiO}_2$ ), boric acid ( $\text{B}_2\text{O}_3$ ), phosphoric acid ( $\text{P}_2\text{O}_5$ ) and aluminium oxide ( $\text{Al}_2\text{O}_3$ ). These substances are capable of absorbing (dissolving) metal oxides up to a certain proportion without losing their glassy character. This means that the incorporated oxides are not involved in the formation of the glass but modify certain physical properties of the structure of the glass as “network modifiers”.

A large number of chemical substances have the property that they solidify from the molten state into a glassy state. The formation of glass depends on its cooling rate and a necessary prerequisite is the existence of mixed types of bond (covalent bonds and ionic bonds) between the atoms or groups of atoms.

As a result, glass-forming products show a strong tendency whilst still in the molten state towards amorphous three-dimensional networking though polymerisation. Crystals are formed when the individual atoms form a regular three-dimensional arrangement in what is known as a “crystal lattice” as soon as the particular substance changes from the liquid to the solid state. Glass, however, forms a largely amorphous “network” when it cools down from the molten state. The components mainly involved in the formation of the glass are therefore described as “network formers”. The glass-forming molecules in this network can incorporate ions that open up the network at certain points, changing its structure and thus the properties of the glass. They are therefore called “network modifiers”.

WHAT IS DURAN® ?

The special features of DURAN®

Very high chemical resistance, nearly inert behaviour, a high usage temperature, minimal thermal expansion and the resultant high resistance to thermal shock are its most significant properties. This optimum physical and chemical performance makes DURAN® the ideal material for use in the laboratory and for the manufacture of chemical apparatus used in large-scale industrial plant. It is also widely used on an industrial scale in all other application areas in which extreme heat resistance, resistance to thermal shock, mechanical strength and exceptional chemical resistance are required.

Chemical composition of DURAN®

DURAN® has the following approximate composition:

81	% by weight	$\text{SiO}_2$
13	% by weight	$\text{B}_2\text{O}_3$
4	% by weight	$\text{Na}_2\text{O}/\text{K}_2\text{O}$
2	% by weight	$\text{Al}_2\text{O}_3$

DURAN® properties are specified in DIN ISO 3585. In contrast to other borosilicate 3.3 glasses, DURAN® is notable for its highly consistent, technically reproducible quality.


Chemical properties

The chemical resistance especially of DURAN® glass is more comprehensive than that of all other known materials. DURAN® borosilicate glass is highly resistant to water, acids, saline solutions, organic substances and also halogens such as chlorine and bromine. Its resistance to alkali is also relatively good. Only hydrofluoric acid, concentrated phosphoric acid and strong alkalis cause appreciable surface removal of the glass (glass corrosion) at elevated temperatures (>100 °C). Due to the nearly inert behaviour, there are no interactions (e.g. ion exchange) between medium and glass and any spurious influence on experiments is thereby effectively excluded.

Hydrolytic resistance


DURAN® corresponds to Class 1 of the glasses that are divided into a total of 5 hydrolytic resistance classes in accordance with ISO 719 (98 °C). The amount of Na<sub>2</sub>O/g glass grain leached out after 1 hour in water at 98 °C is measured. For DURAN® the quantity of Na<sub>2</sub>O leached out is less than 31 µg/g of glass grain. DURAN® also corresponds to Class 1 of the glasses divided into a total of 3 hydrolytic resistance classes in accordance with ISO 720: (121 °C). The quantity of Na<sub>2</sub>O leached out after 1 hour in water at 121 °C is less than 62 µg/g of glass grain. Due to its good hydrolytic resistance, DURAN® meets the requirements of the USP, JP and EP for a neutral glass that corresponds to glass type 1. It can therefore be used in an almost unrestricted way in pharmaceutical applications and in contact with foodstuffs.

Hydrolytic attack on DURAN® as a function of time (100 °C)

Acid resistance


DURAN® corresponds to Class 1 of the glasses divided into 4 acid classes in accordance with DIN 12116. As the surface removal after boiling for 6 hours in normal HCl is less than 0.7 mg/100 cm<sup>2</sup>, DURAN® is classed as acid-resistant borosilicate glass. The quantity of alkaline metal oxides leached out in accordance with ISO 1776 is less than 100 µg Na<sub>2</sub>O/100 cm<sup>2</sup>.

Acid attack on DURAN® as a function of acid concentration

Alkali resistance

DURAN® corresponds to Class 2 of the glasses divided into 3 alkali classes in accordance with DIN ISO 695. The surface erosion after 3 hours boiling in a mixture of equal volume fractions of sodium hydroxide solution (concentration 1 mol/l) and sodium carbonate solution (concentration 0.5 mol/l) is only 134 mg/100 cm<sup>2</sup>.

Alkali attack on DURAN® as a function of temperature

Overview of the chemical properties of technical glasses

Description	Chemical resistance class		
	Hydrolytic resistance DIN ISO 719	Acid resistance DIN 12 116	Alkali resistance ISO 695
DURAN®	1	1	2
FIOLAX®	1	1	2
Soda-lime glass	3	1	2
SBW	1	1	1

Physical properties

Temperature resistance when heated and thermal shock resistance

The maximum permissible operating temperature for DURAN® is 500 °C. Above a temperature of 525 °C the glass begins to soften and above a temperature of 860 °C it changes to the liquid state. As it has a very low coefficient of linear expansion ( $\alpha = 3.3 \times 10^{-6} \text{ K}^{-1}$ ), a feature of DURAN® is its high thermal shock resistance up to  $\Delta T = 100 \text{ K}$ . For a temperature change of 1 K, the glass changes by only  $3.3 \times 10^{-6}$  relative length units, resulting in low levels of mechanical strain were a thermal gradient exists. The thermal shock resistance is influenced wall thickness and product geometry.

Temperature resistance at low temperatures

DURAN® can be cooled down to the maximum possible negative temperature and is therefore suitable for use with liquid nitrogen (approx.  $-196^\circ\text{C}$ ). During such use/ freezing In general DURAN® products are recommended for use down to  $-70^\circ\text{C}$ .

When working at low temperatures, the effect of any expansion of a DURAN® vessel's contents must be borne in mind. During cooling and thawing ensure that the temperature difference does not exceed 100 K. In practice, therefore, stepwise cooling and heating are recommended. When freezing substances in such items as DURAN® bottles or DURAN® test tubes, the container should only be filled to a maximum of 3/4 of its capacity. Moreover, it should be frozen slanted at an angle of  $45^\circ$  (to enlarge the surface area). The minimum service temperature is dependant upon the properties of any screw caps or other components used. For the blue PP screw cap the minimum temperature is  $-40^\circ\text{C}$ .

Use in the microwave

DURAN® laboratory glassware is suitable for use in microwaves. This also applies to plastic coated DURAN® products .

Overview of the physical properties of technical glasses

Description	Linear expansion coefficient $\alpha$ (20 °C/300 °C) [10 <sup>-6</sup> K <sup>-1</sup> ]	Transformation temperature [°C]	Density [g/ cm <sup>3</sup> ]
DURAN®	3.3	525	2.23
FIOLAX®	4.9	565	2.34
Soda-lime glass	9.1	525	2.5
SBW	6.5	555	2.45


Optical properties

In the spectral range from about 310 to 2200 nm the absorption of DURAN® is negligibly low. It is clear and colourless. Fairly large layer thicknesses (axial view through pipes) appear slightly yellow/ greenish. Amber-coloured DURAN® products are suited to use with light-sensitive substances (see amber colouring of DURAN®). This results in strong absorption in the short-wave region up to approx. 500 nm. In photochemical processes the light transmission of DURAN® in the ultraviolet range is of particular importance. The degree of light transmission in the UV range indicates the ease with which photochemical reactions can be carried out, for example chlorinations and sulfochlorination. The chlorine molecule absorbs light in the range from 280 to 400 nm and thus serves as a transmitter of the radiation energy.

Amber colouring of DURAN® laboratory glassware

Amber colouring enables storage of light sensitive substances in DURAN® products. Light transmission in the wavelength range between 300 and 500 nm is, in comparison with DURAN® clear glass, < 10%. Accordingly amber DURAN® glass corresponds to USP/EP/JP specifications.

To colour an existing piece of DURAN® glassware, it is sprayed using an innovative process with a special medium-diffusion ink solely on the outside of the clear glass. The piece is heated to a high temperature which mobilises the diffusing ions into the out layer of glass. On cooling, the ambering is very uniform, resistant to chemicals and cleaning in a dishwasher. The proven DURAN® properties within the bottle remain unaffected; there is no contact or interaction between contents and amber coating. The uniformity of the amber colouring process and the resultant quality of the amber colour is ensured by continuous monitoring.

Transmission curves for DURAN®

CONFORMITY WITH STANDARDS AND GUIDELINES

Alongside the international standard DIN ISO 3585, in which the properties of borosilicate glass 3.3 are defined, DURAN® laboratory glassware corresponds to the current standards for glass laboratory apparatus. The relevant DIN/ISO standards are given on the product pages of this catalogue. If the standard is changed, e.g. in case of harmonisation to ISO, our dimensions are adjusted accordingly within an appropriate time interval.

DURAN® is a neutral glass of high hydrolytic resistance and thus belongs to glass type 1 in accordance with the European pharmacopeia (EP, chapter 3.2.1), the Japanese pharmacopeia (JP, chapter 7.01) and the United States pharmacopeia (USP, section: 660) and National Formulary.

LABORATORY GLASSWARE AND PLASTICS

Plastics used with laboratory glass<sup>1</sup>

To complement DURAN® laboratory glassware products, various plastic products such as screw caps are available. Their properties are listed in the following table.

		Temperature resistance range °C
EPDM	Ethylene/propylene-diene-rubber	-45 to +150
ETFE	Partially crystalline ethylene/tetrafluoroethylene copolymer	-100 to +180
FEP	Tetra-Fluor-Ethylen/Hexafluor-Propylene	-200 to +200
FKM	Fluorinated rubber	-20 to +200
PBT <sup>2</sup>	Polybutylenterephthalat	-45 to +180
PE	Polyethylene	-40 to +80
POM	Polyoxymethylene	-40 to +90
PP	Polypropylene	-40 to +140
PTFE	Polytetrafluorethylene	-200 to +260
PU <sup>3</sup>	Polyurethane	-30 to +135
TpCh260	Thermoplastic/duroplastic	-196 to +260
VMQ	Silicone rubber	-50 to +200
PSU Compound	Compound polyarylsulfone based	-45 to +180

<sup>1</sup> The details in the table relating to temperature resistance relate to standard autoclaving processes with steam and a duration of 20 minutes (see page 213, sterilisation).

<sup>2</sup> Changes in colour may occur if exposed to temperatures above 180 °C.

<sup>3</sup> All laboratory glasses coated with PU may only be cleaned while moist in order to avoid any electrostatic charge which may form.

Chemical resistance of plastics

Classes of substances +20 °C	PE	PP	PBT	PTFE/ FEP	TpCh260	ETFE	VMQ	EPDM	PU	FKM	POM	PSU Compound
Alcohols, aliphatic	+	+	++	++	++	++	+	+	++	–	+	++
Aldehydes	+	+	+	++	++	++	+		++		+	+
Alkaline solutions	++	++	+	++	++	++	–	++	++	–	+	++
Esters	+	+	+	++	++	++	–	++	+	–	–	+
Ethers	–	–	+	++	++	++	–	–	+	–	+	+
Hydrocarbons, aliphatic	–	++	+	++	++	++	–	++	++	++	+	+
Hydrocarbons, aromatic	–	+	+	++	++	++	–	+	++	++	+	–
Hydrocarbons, halogenated	–	+		++	++	++	–	+	–	++	+	–
Ketones	+	+	+	++	++	+	–	++	+	–	+	–
Acids, dilute or weak	+	++	++	++	++	++	–	++	++	++	–	++
Acids, conc. or strong	+	+	+	++	++	++	–	++	+	++	–	++
Acids, oxidising	–	+	–	++	++	+	–	–	+	+	–	+

++ = very good resistance  
+ = good to limited resistance  
– = low resistance

CLEANING LABORATORY GLASSWARE

Special glass laboratory apparatus can be washed by hand in a soaking bath or by machine in a lab washer. Laboratory dealers can supply a wide range of detergents and cleaner-disinfectants for both methods. As contamination during the delivery of our laboratory glassware cannot be totally ruled out, we recommend washing laboratory glassware before it is used for the first time. To care properly for laboratory glassware, it should be washed immediately after use at low temperature, on a short cycle and with low alkalinity. Laboratory apparatus that has come into contact with infectious substances or microorganisms should be treated in accordance with the current guidelines. Dependent on the substance, autoclaving (e.g. to kill microorganisms) may be necessary prior to cleaning, but it is generally recommended that cleaning or washing of glass products be carried out prior to autoclaving or hot-air sterilisation. This prevents dirt or impurities from adhering to the glassware surfaces and prevent damage caused by any possibly adhering chemicals.

Manual cleaning

The generally recognized method is to wipe and rub the glass with a cloth or sponge soaked in cleaning solution. Abrasive cleaners and abrasive sponges should not be used on laboratory glassware as these can damage the surface of the glass. Surface damage can affect the glass properties and limit further use of the product. When soaking glassware it should generally be left in the cleaning solution for 20 to 30 minutes at room temperature, then rinsed with tap water followed by distilled water. To clean the glass as gently as possible, and thus extend its service life, a prolonged soaking time and higher temperature should only be used for stubborn soiling. Laboratory glassware should not be soaked for long periods in strongly alkaline media at more than 70 °C since this can have an adverse effect on the ceramic printing and may cause glass corrosion. Also to be avoided is severe mechanical action, e.g. scraping using a metal spoon.

Washer-disinfectors for automatic laboratory glassware reprocessing

Washer-disinfectors for automatic laboratory glassware preparation are available in various sizes and performance classes. The product range extends from compact machines of 60 to 90 cm width up to powerful, large capacity machines. The large capacity machines are specially intended for central reprocessing of large quantities of laboratory glassware and are available as both 1-door and 2-door barrier machines for installation in a diaphragm wall.

60 cm wide compact machine  
Performance/load: e.g. 37 narrow neck glasses, 96 pipettes


115 cm wide large capacity machine  
Performance/load: e.g. 232 narrow neck glasses, 232 pipettes

Before purchasing a washer-disinfector, you must first be clear which laboratory glassware, and how much of it, requires reprocessing on a day-to-day basis within the laboratory. Once the machine size is specified, the appropriate accessories can be individually selected. The accessories include trolleys and inserts for secure support of the laboratory glassware. Inserts are primarily for holding wide-necked laboratory glassware. Special injector trolleys are offered to thoroughly clean laboratory glassware with a narrow internal diameter. These couple directly to the water supply of the machine and thus ensure that even internal cleaning of the laboratory glassware is correctly carried out. This system ensures that even difficult-to-access points are cleaned, which would be very difficult, or even impossible, to clean manually.

Phases of machine-based reprocessing

Machined-based reprocessing comprises cleaning, rinsing, disinfection (if necessary) and drying of the laboratory apparatus. The following figure shows a typical programme cycle for laboratory glassware reprocessing.

Example

Pre-rinse: cold water without process chemicals  
Cleaning: cold or hot water with alkaline cleaning agent  
Neutralisation: cold or hot water with acidic neutralisation agent

Example

Rinse I: cold water  
Rinse II: deionised or ultrapure water  
Flushing: deionised or ultrapure water at 75 °C

Cleaning

Cleaning removes dirt from the surfaces. At this stage, process chemicals (e.g. cleaning agents, surfactants, emulsifiers, neutralizers) are used. Cleaning may comprise several programme blocks, such as pre-rinse, cleaning, neutralisation.

Rinsing

During rinsing the dissolved dirt and the process chemicals used are rinsed off. Rinsing can comprise a number of individual programme blocks. The choice of water quality (e.g. tap water, deionised water, ultrapure water) depends on the application (e.g. organic/inorganic analysis, microbiology).

Disinfection

During disinfection, infectious contamination is killed/inactivated to such a degree that the laboratory glassware no longer represents an infection risk. On the one hand, disinfection serves to protect personnel within laboratories who work with infectious contamination. On the other hand disinfection prevents transfer of germs from samples and preparations in medical laboratories, hygiene institutes, pharmaceutical laboratories and the food and cosmetic industries. Thus hygienic, problem-free working is guaranteed.

Drying

The washer-disinfectors have, dependent on model and construction, an active hot-air dryer which permits not only drying of the external surfaces, but also drying of narrow diameter laboratory glassware. Also laboratory glassware of complex shape is reliably dried using hot-air drying. To effectively protect the laboratory glassware against dust particles and microorganisms, the drying air is passed through a HEPA filter.

Typical programme using a Miele washer-disinfector for reprocessing of laboratory glassware:

Miele washer-disinfectors for laboratory glassware reprocessing have up to 10 standard programmes. Numerous programme parameters can be adjusted to adapt the standard programmes for particular customer applications. Moreover, customer-specific programmes can be created for special applications.

Inorganic	To remove acid-soluble inorganic residues
Organic	To remove heavy organic residues such as oil, grease, wax, agar
Standard	Simple standard programme for slightly soiled glassware with a low final-rinse requirement
Universal	To remove organic residues (e.g. proteins, oils), for medium-level dirt and a medium final-rinse requirement
Intensive	To remove organic residues (e.g. proteins, cell and tissue cultures, oil), for heavy levels of dirt and a high final-rinse requirement
Plastic	For temperature-sensitive laboratory equipment (e.g. plastic bottles) with a low to medium level of dirt and a medium final-rinse requirement
Vario TD	For cleaning and heat disinfection at 93 °C with 5 minutes temperature-holding time, in accordance with EN ISO 15883-1, disinfection in the last rinse block
Special 93°C-10	For cleaning and heat disinfection at 93 °C with 10 minutes temperature-holding time, disinfection in the first rinse block, used in the case of an outbreak of a notifiable disease.

Analysis purity through conductivity measurement in the final rinse

The requirements for analysis purity depend largely on the application of the laboratory glassware. To ensure analysis purity, washer-disinfectors for laboratory glassware reprocessing can optionally be provided with a conductivity measurement module. An integrated conductivity measurement offers the following advantages:

- Detection of undesirable contents in the rinse water (dissolved salts of alkaline or acidic process chemicals)
- Definition of a customer-specific permissible conductivity level

Process reliability for reproducible results

Automatic preparation is a validatable preparation process that delivers reproducible results. This is one reason why automatic preparation should be favoured over manual processes. To guarantee the reproducibility of the results, the machines have the following safety installations:

- Temperature monitoring using two redundant temperature sensors
- Automatic liquid dosing including dosing volume control
- Spray arm rotation speed monitoring

Process documentation

In applications which require high standardisation and reproducibility, process documentation contributes significantly to quality control. Process documentation can take place via documentation software or a printer.

Economy

Nowadays, laboratory glassware preparation must constantly meet ever higher requirements in respect of performance and economy. Machine-based reprocessing is by comparison with manual cleaning, much more efficient: for example, the economy arises from lower time / personnel expenses, shorter process cycles as well as lower power and water consumption. In particular, the short process cycles mean the laboratory glassware is quickly ready for its next use. Minimal handling of contaminated laboratory glassware simultaneously reduces the potential risk to personnel (injury, chemical burns and risk of infection).

Value retention through gentle preparation

Automatic laboratory glassware preparation is gentler than manual cleaning. The glass surfaces only comes into contact with the alkalinity of the detergent for a short, defined time interval, so that glass corrosion is minimized. The accessories include special holders and locks so that the laboratory glassware is securely fastened and protected against breakage.

DURAN GROUP recommends Miele

To guarantee thorough, gentle and safe laboratory glassware preparation, DURAN GROUP recommends Miele washer disinfectors.

Miele "Made in Germany" quality is notable for its high reliability and efficiency in day-to-day use in the laboratory. Short operating times and reliable results ensure that high-value laboratory glassware is once again ready for use after only a short period. In addition, the gentle preparation also ensures a long service life for DURAN® laboratory glassware.

STERILISATION

When preparing laboratory glass for sterile applications or as part of the cleaning process, sterilisation is a well-established process. DURAN® laboratory glassware is suitable both for autoclaving as well as for hot air and plasma sterilisation ( $H_2O_2$ ). Laboratory apparatus that has come into contact with infectious substances or microorganisms must be cleaned in accordance with the appropriate guidelines for handling these materials. As the case may be, this may include sterilisation.

When carrying out sterilisation, especially of laboratory glassware, the following instructions should be observed: To avoid overpressure, all vessels should always be kept open. When sterilising media, the use of a membrane cap is recommended. Such a cap permits pressure equalisation through a PTFE membrane and hence the cap can be tightly closed. Consequently, the risk of contamination is greatly reduced.

Alongside the standardised procedures described above, individually modified methods are also applicable to all DURAN® products, for example using higher temperatures. However, you must ensure, especially with bottles (due to the screw caps) that the permissible highest temperatures for the plastic used in the accessories is not exceeded (see page 208).

WORKING UNDER PRESSURE

Only products whose design includes the appropriate geometry and wall thickness, and which are explicitly designated as such, are suitable for working under pressure and / or vacuum (e.g. filtering flasks, desiccators or flat flange vessels).

When used under positive or negative pressure, and especially when also working with differential temperatures, additional care measures must be taken. Glass apparatus that is under pressure or vacuum should only be subject to further stress (e.g. significant temperature change) with extreme caution, as the individual resulting strains are additive and could readily result in failure.

To guarantee optimum user safety, the following points should be borne in mind:

- To avoid stresses in the glass, evacuated vessels or vessels under pressure should not be heated on one side or heated with an open flame.
- When working under pressure the maximum figures indicated in the catalogue should not be exceeded.
- Before using glass equipment under vacuum or pressure it must always be visually inspected to check that it is in perfect condition (no serious scratches, micro-cracks, abrasions, etc.). Damaged glassware should not be used for work under pressure or vacuum for safety reasons.
- Never subject glassware to sudden pressure changes, e.g. always re-pressurise evacuated glass apparatus slowly.
- Laboratory glassware with a flat bottom (e.g. Erlenmeyer and flat bottom flasks) should not be used under pressure or vacuum.
- The plastic coating of laboratory bottles (DURAN® protect) has no influence on pressure resistance. These products are not designed for use under pressure. For pressure applications using laboratory bottles, the DURAN® pressure plus bottle should be used. The DURAN® pressure plus bottle is pressure resistant from -1 to 1.5 bar due to a modified geometry and increased wall thickness.

SAFETY INSTRUCTIONS

When used according to our specifications, DURAN® glassware is very safe to use. The appropriate guidelines applicable for the use of special glass in laboratories in the country in question should always be complied with.

The following points should, however, be observed in every case:

- For safety reasons, before DURAN® laboratory glassware is used it should be checked to ensure that it is suitable for the intended purpose and that it can be used without problem.
- Defective laboratory glassware represents a risk (e. g. risk of cuts, burns, infection) that should not be underestimated. If appropriate repairs to any item cannot be carried out or cannot be justified for economic reasons, it must be disposed of in the proper manner.
- Repairs must only be carried out by competent, specialist, glassworkers. Poorly, repaired glassware can fail without warning and represents a significant hazard.
- Only subject DURAN® glassware to sudden temperature changes within the recommended limit for thermal shock resistance ( $\Delta T = 100K$ ). This means that hot laboratory glassware should not be taken out of a drying cabinet and placed on a cold or even wet laboratory bench. This applies in particular to thick-walled glassware such as filtration flasks and desiccators.
- When assembling apparatus make sure that it stands firmly and is not subjected to stress by using appropriate stands.

DISPOSAL

DURAN® laboratory glass should under no circumstances be disposed of in the domestic glass recycling stream (e. g. bottle banks), since its high melting point and different chemistry make it incompatible with other glass cullet (soda-lime glass) for recycling. The correct way to dispose of it is, in principle, to include it with general household waste (residual waste) in accordance with the relevant guidelines, provided that the glass is quite free of any harmful contamination (Waste code no: 17 02 99 D1).

DURAN® LABORATORY GLASS IS ECO-FRIENDLY

DURAN® laboratory glass is made from natural, mineral raw materials. Unlike other materials, laboratory glass, when used properly, will give years of service and this means that it is vastly superior to other materials from an ecological viewpoint too. Depending on its use, DURAN® can be disposed of as household waste and does not need to be dealt with as special waste which may be environmentally harmful. Toxic substances cannot leach out because of the raw materials used.


Production processes in our factories have been consistently optimized over recent years to ensure that they are environmentally friendly during the actual manufacturing stage through the minimum usage of valuable resources. Electrical heating and advanced technology in our melting units ensure that no pollutants are released during manufacture in our ultramodern factories which could harm our workers or people living nearby. In addition energy demand is kept as low as possible. The latest waste gas purification equipment is used to avoid emissions which could pollute the environment. A significant investment has been made in an enclosed cooling water recirculation system to cut the amount of fresh water required to a minimum thus helping to conserve vital water resources. We use packaging made from environmentally harmless, recycled paper which can be returned after use to the resource cycle.

FURTHER PROCESSING

DURAN® items made of borosilicate glass 3.3 are suitable for further processing such as the addition of screw thread tubes, olives, tubulatures, necks and ground glass joints. Preferred for further processing are round, flat bottom and Erlenmeyer flasks. Certain sections of the temperature / viscosity range are of particular importance for glassworking. In the transformation range the elastic-brittle behaviour of the glass changes with increasing temperature into a markedly viscous one, so that consequently all its physical and chemical properties change significantly with temperature. The transformation temperature range thus plays an important part in stress relief during heating up and the introduction of stress when the glass is cooled. The position of the transformation range is identified by the transformation temperature "Tg" DIN 52 324.

Note:

DURAN GROUP cannot accept any product liability where items are subjected to further processing. In this case the entire responsibility for quality lies with the glassworker. The latter is therefore responsible for ensuring that the further processed item conforms to current directives and safety requirements.

Normal temperature dependence/viscosity curve of, for example, DURAN®; viscosity ranges of important processing techniques, position of fixed points of viscosity and various limiting temperatures.

BOTTLES

Laboratory bottles

DURAN® laboratory bottles are chemically resistant and stable. The extensive range of original accessories includes screw caps for the widest possible range of applications. Alongside the standard PP screw cap for everyday laboratory use, further caps made from various plastics and having special properties are available. DURAN® laboratory bottles are completed by suitable pouring rings from different plastics, which enable drip-free working. As almost all GL 45 bottles of 100 ml capacity and above use the same thread size, screw caps and pouring rings are fully interchangeable. The bottles, pouring rings and caps are autoclavable/sterilisable.

Properties

Light protection

- amber bottles are opaque up to 500 nm
- plastic coated bottles are opaque up to 380 nm
- Application: storage of light sensitive substances

High thermal shock resistance

Due to their temperature properties, the bottles are suitable for autoclaving and sterilising (see page 213, general section). Because of the bottom geometry and the wall thickness, direct heating with an unshielded flame is not recommended. When using an electronic heating plate or water bath laboratory bottles should be heated gradually.

Notes on use:

Pressure resistance

DURAN® laboratory bottles are, with the exception of the pressure-resistant DURAN® pressure plus bottles, in general not suitable for use under pressure or in a vacuum. DURAN® pressure plus bottles are pressure resistant from -1 to +1.5 bar (overpressure) due to a modified geometry and increased wall thickness.

Sterilisation

When sterilising or autoclaving contents, the screw cap must only be loosely fitted (max. one turn). The contents may expand or boil causing a large pressure difference in a closed vessel, which may well result in explosive failure. Alternatively, a DURAN® membrane cap may be used. Pressure equalisation takes place through the PTFE membrane, while at the same time the membrane cap can remain tightly closed, greatly reducing the risk of contamination. See also page 213, general part.

Cleaning

Cleaning should be carried out manually in a soaking bath or automatically in a dishwasher (see page 209, general part). When cleaning in a dishwasher, load so that there is no glass-to-glass contact (especially the threads) to avoid chips or abrasions.

Freezing substances

Recommendation: The bottle should be frozen slanted at an angle of 45°, filled to a maximum 3/4 (to enlarge the surface area) and dependent on the properties of any screw caps or other components used. For the blue PP screw cap the minimum temperature is -40 °C. Alternatively the Premium screw cap can be used (min. working temperature: -196 °C). See also page 206, general part.

Thawing frozen substances

Frozen contents can be thawed by immersing the bottle in a liquid bath while taking care that the

temperature difference between the contents and the bath does not exceed  $\Delta T = 100\text{K}$ . This will ensure that the frozen material is warmed uniformly from every side without damaging the bottle. The contents can, however, also be thawed slowly from above, so that the surface melts first, allowing the material to expand.

Laboratory bottles with plastic coating

The coating of DURAN® Protect bottles is a resistant and transparent plastic coating based on a cross-linked copolymer.

The coating adheres securely to the glass surface and fulfils the following functions:

- Protects the glass surface against mechanical damage (scratch protection).
- Holds the fragments together in the event of the glass breaking (splinter protection).
- Minimises liquid loss if the glass breaks (protects against contents escaping and splash).
- Absorbs UV rays up to a light wavelength of 380 nm (light protection).

Usage tips:

- The plastic coating does not increase the pressure resistance. These bottles are designed for pressure or vacuum applications.
- If the plastic coated bottle breaks during use, the contents and the plastic coating are likely to come into contact. A test for any interaction between plastic and contents should be carried out to ensure that the contents remain unchanged and can be further used.

Temperature resistance

Do not expose DURAN® protect bottles to open flames or direct heat, e.g. on a laboratory hotplate. The maximum operating temperature is +135 °C and thus the bottle is suitable for use in an autoclave. Long-term exposure to temperature (> 30 minutes) should be avoided. DURAN® protect bottles can be used for freezing to -30 °C and used in microwaves. Thermal and chemical stresses can result in coating discolouration.

Autoclaving:

The following procedure, bearing in mind the maximum temperature resistance, is recommended:

- Steam sterilisation at +121 °C or +134 °C.
- The cycle duration should not exceed 20 minutes.

(See also page 213, general section)


When sterilising, the screw cap should only be loosely applied (max. one turn - do not tighten), or use a membrane cap that allows pressure equalisation.

Gas washing bottles

By distributing the gas through the liquid by use of a filter disk, the gas surface is significantly increased and the interchange between gas and medium is improved. DURAN® gas washing bottles also work reliably at high flow velocities. The graph illustrates the effectiveness of gas washing bottles with and without a gas filter disk.

Absorption efficiency of two gas wash bottles:  
A without gas filter and B with gas filter plate

Filtering flasks with side-arm socket or plastic hose connection

DURAN® filtering flasks are vacuum tight in accordance with DIN 12 476, ISO 6556. Alongside the filtering flasks with glass hose connections, versions are also available with a side-arm socket or plastic hose connection. The ground side-arm socket with dimensions 17.5/26 is suitable for vacuum hoses from 15 to 18 mm OD (e.g. 6 x 5 mm or 8 x 5 mm, DIN 12 865). The plastic hose connections are suitable for hoses of approx. 9 mm internal diameter. The versions with side-arm socket or plastic hose connection offer improved safety for the user.

DESICCATORS

DURAN® desiccators are used for drying moist substances or as storage vessels for moisture-sensitive products. To accelerate the drying process, the desiccators can be used under vacuum. Due to the high wall-thickness of the vessels and the exact machining of the vacuum-tight ground joints on the lid and base, storage under vacuum is possible - even over extremely long periods.

All individual parts and a wide range of accessories such as lids, stopcocks, bases, etc. are compatible and can be interchanged as required. Always ensure the individual parts have the same DN (nominal diameter in millimetres).

For desiccators, the DN is based on the diameter of the sieve plate; this, or the lip it rests upon in the desiccator base, can be measured directly. For lids, measure the outside diameter of the flange and cross-reference with the tables on the product pages. An overview on page 80 indicates which individual parts are required to assemble the desired desiccator.


Usage tips:

- Designed for use under full vacuum (-1 bar)
- Due to the high wall thickness and the reduced thermal shock resistance under pressure loading, the desiccators must not be heated on one side only or heated using a naked flame.
- Before evacuation, it is recommended that the glass surfaces of the desiccator be checked for damage such as scratches, cracks or nicks.  
Damaged desiccators must not be used for safety reasons.
- Never expose desiccators to abrupt pressure changes (do not suddenly ventilate evacuated vessels).

CENTRIFUGE TUBES AND CULTURE TUBES

Centrifuge tubes

DURAN® centrifuge tubes are approved in accordance with DIN 58 970 (Part 2) up to a maximum relative centrifugal acceleration (RZB = 4000) and for filling up to their capacity with contents having a maximum density of 1.2 g/ml.

Calculation:  
 $RZA = 1,118 \times 10^{-5} \times r \times n^2$   
 $n = \sqrt{\frac{4000}{1,118 \times 10^{-5} \times r}}$

Example:  $r = 15 \text{ cm}$   
Example in the diagram:  
number of revolutions (n)  
 $= 4900 \text{ min}^{-1}$

Culture tubes

In addition to DURAN® culture tubes, our product range also includes soda-lime culture tubes. This is a glass belonging to the third water resistance class and is one of the soda-lime glasses with a high fraction of alkaline and alkaline earth oxides.

Properties of soda-lime glass:

Physical data		Chemical data							
Linear expansion coefficient		Hydrolytic class		(ISO 719)		3			
$\alpha_{20/300}$ to DIN 52328:	$9,1 \times 10^{-6} \text{ K}^{-1}$	Acid class		(DIN 12 116)		1			
Transformation temperature Tg:	525 °C	Alkali class		(ISO 695)		2			
Temperature fixed points at viscosity $\eta$ n dPa x s:									
$10^{13}$ upper annealing temperature	530 °C	Chemical composition							
$10^{7,6}$ softing temperature	720 °C	(main components in approx. weight %)							
$10^4$ working temperature	1040 °C	SiO <sub>2</sub>	B <sub>2</sub> O <sub>3</sub>	K <sub>2</sub> O	Al <sub>2</sub> O <sub>3</sub>	Na <sub>2</sub> O	BaO	CaO	MgO
Density $\rho$ :	2,50 g/cm <sup>3</sup>	69	1	3	4	13	2	5	3


FLAT FLANGE RANGE

The DURAN® flat flange reaction vessels are valued for their universal suitability for use in the laboratories of a wide range of specialisations. Whether for reaction, distillation, evaporation or desiccation, DURAN Group offers a wide range of unfinished and finished parts which always provide the optimum solution for the particular application. Due to the pure glass-glass connections, reactions with highly corrosive or highly chemically reactive substances can be carried out without problem.

The vessels are notable due to a robust glass flange design with an optimum flange angle of 45°. The proven flange design (flat ground) is available with groove, consequently O-Rings can be used. The corresponding stainless-steel quick release clamps with three flexible retaining clips ensure easy and safe handling. All individual parts and a wide range of accessories such as lids, O-Rings and quick-release clamps etc. are compatible and can be interchanged as required. In so doing however, you must always ensure the same DN (nominal diameter) of the individual parts applies.

Usage tips:

- All components are suitable for use under full vacuum (-1 bar). Many are rated for positive pressure operation (see product descriptions for details)
- Before use, it is recommended that the glass surfaces be checked for damage such as scratches, cracks or nicks.
- Damaged glassware should not be used for safety reasons.
- Due to the high wall thickness and reduced thermal shock resistance under pressure loading, the flat flange vessels should be heated uniformly and gradually.

Beaded lid for safer handling of the reaction vessel

Accessories

Flat flange reaction vessels can be sealed by:

- a) O-rings (see below) for use at positive and negative pressures up to max. 230 °C (O-ring dependent)

Advantages:

  - Easy to open
  - The lid does not stick, even after operation for long periods under vacuum and at high temperatures
  - Reduced need to grease contact surfaces
- b) Grease for use at positive and negative pressures above 230 °C

Advantages:

  - Flange groove acts as an ideal grease groove
  - No substance contamination

The stainless steel quick release clamps with three holding segments are optimally designed to provide even distribution of contact pressure. The chromium nickel steel support comprising two clamping rods is designed for secure fitting of the reaction vessels or the lids in support bar. For example, if there is a need to change the lid or the vessel, this can be done without dismantling the entire apparatus.

Shape retentive O-rings

FEP seamlessly coated elastomer O-rings with silicone core

Comprising an elastic, silicone core with a seamless FEP coating that encloses the ring. The combination of these high-quality materials ensures good elasticity in conjunction with outstanding chemical resistance. The chemical resistance of FEP (tetrafluoroethylene hexafluor-propylene copolymer) is equal to that of PTFE. Hence the material is resistant to almost all chemicals and is suitable for temperature from - 200 °C to + 200 °C.

Silicone (VMQ) O-rings

These O-rings are made solely from silicone (VMQ) and therefore are highly elastic. Their chemical resistance, however, is reduced in comparison with FEP coated O-rings. Temperature resistance extends from -50 °C to +230 °C.

	Elasticity/recovery	Temperature resistance	Chemical resistance	Solvent resistance	Physiologically harmless
O-Rings, red FEP coated	+	++	++	++	++
O-Rings, transparent made of silicone (VMQ)	++	++	+	+	++

++ = very good resistance  
+ = good resistance

FILTERS AND FILTRATION APPARATUS

Porosity

DURAN® filters and the corresponding filter plates are precision manufactured from DURAN® borosilicate glass 3.3 with its high chemical and thermal-shock resistance. They are entirely inorganic and inert in most circumstances. There are therefore no leachable organic or ionic species present that could otherwise contaminate the filtrates. They are ideal for separations, e.g. with strong acids or alkalis and can likewise be readily cleaned and reused. DURAN® filter products have a maximum operating temperature of +450 °C.

DURAN® filtration vessels are specially optimised to the matching filtration apparatus (eg funnels with guko adapters) and are vacuum-tight due to their special geometry and high wall thickness. Their designs have been approved by the TÜV accreditation body and marked with the “GS” indication were appropriate; see specific products for details.

Porosity measurement is by the Bechhold bubble pressure method, which is widely described in the literature<sup>1</sup>. In the interests of rapid filtration every effort is made to produce filter disks with as many open pores as possible without blockages or closed cavities. This is one of the areas where DURAN® glass filters stand out.

Prerequisite for the successful use of glass filters is selection of the correct porosity. In this respect, the following table lists details of six porosity ranges with indications of their main areas of application. A point to be borne in mind is that the filtration equipment should ideally be selected to ensure that the nominal size of the largest pore is somewhat smaller than the smallest particles to be filtered out. This will prevent infiltration of particles into of the pores.

For quantitative analysis applications, porosity 3 or porosity 4 glass filtration apparatus is used almost exclusively. Different working methods often contain different porosity indications here for the same materials. This is because different processes used in the production of precipitations for gravimetric analysis often result in different grain sizes.

Porosity	New identification ISO 4793	Nominal max. pore size µm	Areas of application, examples
0	P 250	160–250	<b>Gas distribution:</b> Gas distribution in liquids with low gas pressure. Filtration of coarsest precipitates.
1	P 160	100–160	<b>Coarse filtration, Filtration of coarsest precipitates.</b> <b>Gas distribution in liquids</b> Liquid distribution, coarse glass filters, extraction apparatus for coarse-grained material. Substrates for loose filter layers against gelatinous precipitates.
2	P 100	40–100	<b>Preparative fine filtration:</b> Preparative work with crystalline precipitates. Mercury filtration.
3	P 40	16–40	<b>Analytical filtration:</b> Analytical work with medium fine precipitates. Preparative work with fine precipitates. Filtration in cellulose chemistry, fine glass filters. Extraction apparatus for fine grained material.
4	P 16	10–16	<b>Analytical fine filtration:</b> Analytical work with very fine precipitates. (e.g. BaSO <sub>4</sub> , Cu <sub>2</sub> O). Preparative work with correspondingly fine precipitates. Non-return valves and check valves for mercury.
5	P 1,6	1,0–1,6	<b>Ultrafine filtration</b>

Flow rate

To determine the possible applications of glass filter disks and filtration apparatus, it is necessary to know not only the porosity, but also the flow rates of liquids and gases. These are given in Figures 9 and 10 for water and air. The data applies to 30 mm diameter filter disks.


The flow rates for other disk diameters can be calculated by multiplying the value read off by the conversion factor given in Table the following table:

Filter disk diam. mm	10	20	30	40	60	90	120	150	175
Conversion factor	0,13	0,55	1	1,5	2,5	4,3	6,8	9,7	15

Example:


Suction filtration of an aqueous solution under vacuum using a suction filter with a 60 mm disk diameter and porosity 4. Figure 9 gives a flow rate of 200 ml/min for a pressure differential of about 900 mbar. Table 8 gives a flow volume of 200 x 2.5 = 500 ml/min for a 60 mm disk diameter. As the flow rate is heavily dependent on the pore diameter (pore radius to the power of 4), deviations from the values indicated may occur. Flow can also be obstructed by the formation of a filter cake over the surface of the filter disk. Further changes to the flow rate occur if liquids are used whose viscosity differs from that of water. The resultant flow rate is then inversely proportional to the viscosity. Differences for gases result when using filter disks that are coated with water or other liquids (gas flow in washing processes). More detailed information can be found in the literature<sup>1</sup>.

Water flow rate

Water flow rate through filter discs of various porosities as a function of pressure differential.  
For filter discs with Ø 30 mm

Dry-air flow rate

Air flow rate through filter discs of various porosities as a function of pressure differential.  
For filter discs with Ø 30 mm

<sup>1</sup> Frank, W: GIT (1967) Iss.7 pp. 683-688

Care and cleaning of filtration apparatus

In addition to the information in the general section on page 209, please also note the following guidelines relating to thermal stresses, which apply specifically to filtration apparatus, in order to avoid glass breakage.

Temperature changes (thermal shock), drying and sterilisation

- The maximum permissible operating temperature is +450 °C.
- Uniform heating is recommended to avoid thermal stresses and resultant breakages.
- Heat glass filtration apparatus with disk diameters of more than 20 mm in initially cold ovens or sterilisers only.
- The heating or cooling rate should not exceed 8°C/min.
- When filtering hot substances avoid temperatures differences of more than 100 K; if necessary, preheat the filtration apparatus in a drying cabinet.
- Wet filtration apparatus should be heated slowly up to 80 °C and dried for one hour before increasing the temperature further.

Whenever possible, filtration apparatus should be stood on its rim (stem upwards) to allow air convection between the inside of the vessel and the oven chamber. If placing the filtration apparatus in the oven at an angle cannot be avoided (as in the case of pipeline filters), any support point close to the position of the filter weld must be protected against heating up prematurely by placing heat-insulating material under it.

Cleaning new glass filtration apparatus

Before using glass filtration apparatus for the first time, it should be rinsed with water (if applicable, acid), to remove any minor contamination that may be present.

Mechanical cleaning

In many cases, if no precipitate has infiltrated the pores, simple spraying of the surface (e.g. with a spray bottle) will suffice. Brushes or rubber wipers can also be used to clean the surface of the filter disk. If some precipitate has infiltrated into the pores, then back-flushing of the disk is required.

Important instructions:

- Glass filters should always be cleaned immediately after use.
- Do not use sharp objects to remove the filtrate to prevent damage to the filter surface.

Chemical cleaning

If some of the pores on the filter disk still remain clogged after mechanical cleaning or if it is desirable to make sure that no residue from previous work remains before filtering a new substance, then thorough chemical cleaning is necessary. The choice of solvent used depends on the nature of the contamination (see example in the following overview).

Barium sulfate	hot conc. sulfuric acid
Silver chloride	hot ammonia liquor
Red copper oxide	hot hydrochloric acid and potassium chlorate
Mercury residue	hot conc. nitric acid
Mercury sulfide	hot aqua regia
Albumen	hot ammonia liquor or hydrochloric acid
Grease, oil	acetone, isopropanol
Other organic substances	hot conc. sulfuric acid with addition of nitric acid, sodium nitrate or potassium dichromate

When chemical cleaning is completed, it should be followed by thorough rinsing with copious amounts of water. Use of hot concentrated phosphoric acid and hot alkali solutions is not recommended, as these may attack the glass surface.

Screwfilters with interchangeable filter disks

With 3 filter sizes, each having 4 filter disks of varying porosity, 12 different filter rates are available. DURAN® screwfilters have a range of benefits compared with conventional filter apparatus:

- Interchangeable filter disks
- Safe and simple removal of the filtered material
- Disks have longer service life, as no damage is caused by scraping off the filtered material
- Filter disks are easy to clean from both sides
- Slit sieve (Cat. No. 21 340 31) can be used in the medium sized screwfilter to support membrane and paper filters
- Space saving
- Cost-effective; filter disks and apparatus can be ordered individually, as required.

Note:

The filter disk should be located between 2 FKM gaskets.

VOLUMETRIC PRODUCTS

DURAN® volumetric products have closely calibrated scales that permit very accurate determination and measurement of volumes. They are available in two accuracy classes: class A/AS and class B. The two classes differ in the accuracy of measurement with class A being the highest accuracy, and class B is approximately half that of class A. Class AS has the same tolerances as class A, but is designed to permit more rapid outflow; it is applicable to burettes and pipettes.

Volumetric flasks

DURAN® volumetric flasks are manufactured from the chemically highly resistant borosilicate glass 3.3. Used for the accurate measurement of specific quantities of liquid they are, like virtually all volumetric glassware, volumetric analysis aids. They are mainly used for preparation and storage of standard solutions. Calibration is based on the amount of fluid contained ("In") at a +20 °C reference temperature, which means that when the circular graduation mark is reached, exactly the specified liquid amount is contained in the vessel. Thus the desired concentration can be precisely set. The volume content tolerances for volumetric flasks conform to accuracy class A, the accuracy limits of the German weights and measures regulation and to DIN and ISO guidelines.

Measuring and mixing cylinders

DURAN® measuring and mixing cylinders are manufactured from borosilicate glass 3.3 and therefore are very resistant to mechanical and thermal stresses. Measuring cylinders are for holding and simultaneously measuring different liquid amounts. Mixing cylinders are for diluting solutions and mixing several components in a given quantity ratio. Their large hexagonal base prevents the cylinder from rolling. The base is equipped with three knobs that increase its stability. The cylinders have uniform wall thickness over the entire measurement range, so wedge errors are avoided. Calibration is based on contained fluid ("In") at a +20 °C reference temperature, which means that when the circular graduation mark is reached, exactly the specified liquid amount is contained in the vessel. Thus the desired concentration can be precisely set. Volume content tolerances for measuring and mixing cylinders conform to DIN and ISO accuracy limits.

Burettes

DURAN® burettes are manufactured from chemically highly resistant borosilicate glass 3.3. They are primarily used for titration. The precise scale permits exact reading of the liquid quantity required for the titration. Calibration is based on the released volume ("Ex") at a +20 °C reference temperature. The fluid quantity released can be taken exactly from the scale, as the liquid adhesion to the glass is taken into account in the calibration. This only applies, however, if the specified waiting times for reading the scale are adhered to. Volume content tolerances for burettes conform to DIN and ISO accuracy limits. The DURAN® Class B burettes' accuracy limits are roughly one and a half times the Class AS accuracy limit. The tolerances are thus stricter than specified by DIN.

By the specification of a class "AS", the German weights and measures regulations have, within the scope of the 15th Amendment Regulations, acknowledged that the great majority of volumetric measurements, especially in clinical laboratories, are carried out with water or dilute aqueous solutions; thus apparatus with considerably shorter draining times than previously required but with the same accuracy limits is now admitted by the calibration regulations.

Capacity  ml	Accuracy limits class AS suitable for official calibration DIN 12 700  ± ml	Accuracy limits class B	
		DIN 12 700  ± ml	DURAN  ± ml
1	0.01	-	-
2	0.01	-	-
5	0.01	-	-
10	0.02	0.05	0.03
25	0.03	0.05	0.04
50	0.05	0.1	0.08
100 <sup>1</sup>	0.08	0.2	0.15

<sup>1</sup> Non-DIN size.

Pipettes

Measurement and bulb pipettes are made from soda-lime glass (see page 219, for more information on soda-lime glass). Pipettes are for precise measurement and filling of liquids. Measurement pipettes are graduated to permit the taking up of varying liquid quantities and then dispensing of the same or different amounts. Bulb pipettes are designed to repeatedly take up and discharge a fixed volume for each pipette size. Calibration is based on the released volume ("Ex") at a +20 °C reference temperature. The fluid quantity released can be taken exactly from the scale, as the liquid adhesion to the glass is taken into account in the calibration. This only applies, however, if the specified waiting times for reading the scale are adhered to. Volume content tolerances for calibrated pipettes conform to DIN and ISO accuracy limits. DURAN® Class B pipettes' accuracy limits are roughly one and a half times the Class AS accuracy limit. The tolerances are thus stricter than specified by DIN.

By the specification of a class "AS", the German weights and measures regulations have, within the scope of the 15th Amendment Regulations, acknowledged that the great majority of volumetric measurements, especially in clinical laboratories, are carried out with water or dilute aqueous solutions; thus apparatus with considerably shorter draining times than previously required but with the same accuracy limits is now admitted by the calibration regulations.

Capacity  ml	Accuracy limits class AS suitable for official calibration ISO 385  ± ml	Accuracy limits class B	
		ISO 385  ± ml	DURAN  ± ml
0.1 <sup>1</sup>	-	-	0.003
0.2 <sup>1</sup>	-	-	0.004
0.5	-	0.01	0.008
1	0.006	0.01	0.008
2	0.010	0.02	0.015
5	0.030	0.05	0.040
10	0.050	0.10	0.080
25	0.100	0.20	0.150

<sup>1</sup> Non-ISO size, graduated pipettes 0.1 and 0.2 ml are calibrated to contain ("In").

Usage tips:

- To ensure a long service life for your volumetric glassware and to exclude possible volume changes, these products should not be heated above +180 °C in drying cabinets or sterilisers.
- Never heat volumetric glassware on a hot plate.
- Always heat up and cool down volumetric glassware gradually, to avoid thermal stresses and thus any possible breakage of the glass.

## GLASS-CERAMIC LABORATORY PROTECTION PLATES

Due to low thermal expansion stresses, these glass ceramic plates are well suited to heating glassware with a Bunsen burner.

**Energy and time savings**

The high transparency to infrared radiation means heat energy is transferred to the material being heated with low losses that shortens heating time and results in energy savings of 20% or more. In addition, several vessels can be placed on the plate's square, stable surface.

**Chemically resistant**

When working in the laboratory it is impossible in practice to avoid aggressive media boiling over or spilling. The glass-ceramic laboratory protection plate is resistant even against highly corrosive media.

**Trouble-free cleaning**

The pore-free smooth surface of the glass-ceramic laboratory protection plate can be cleaned easily either manually or mechanically.

**High temperature resistance**

Service temperature from  $-200\text{ }^{\circ}\text{C}$  to  $+700\text{ }^{\circ}\text{C}$ . The glass-ceramic laboratory protection plate is may be used continuously at high temperatures. Durability at  $700\text{ }^{\circ}\text{C}$ : 6000 h; at  $750\text{ }^{\circ}\text{C}$ : 750 h. Even when a hot plate is quenched with cold water, there is no risk of breakage, since it is resistant to thermal shock even with a  $\Delta T > 650\text{ K}$ . To avoid overheating, care must be taken not to exceed the above-mentioned limits when working with a Bunsen burner. The glass-ceramic laboratory protection plate retains its shape, remains flat and does not age.

**Note:** Further information about DURAN® laboratory glassware is available upon request.

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A

Acid bottle	46
Adapter	
filter crucible	123
receiver	162–163
stopper type*	
vacuum receiver	163
Allihn	
condenser	154
filter tube	131
Analytical funnel	120
Aspirator bottle	47–48
Assembly set, for filtering flask, KECK™	134
Automatic burette	112

B

Beaker	
Berzelius	55
filtering, heavy-wall	54
flat flange	144
high form	53–54
low form	53
Philips	55
Bearing, KPG stirrer	160–161
Bell jar	68–69
Berzelius beaker	55
Bloom test vessel	56
Blowball	112
Bottle	
aspirator	47–48
culture	97
culture media	99–101
gas washing	127–129
GL 45 laboratory bottle, amber	23
GL 45 laboratory bottle, clear	22
GL 45 laboratory bottle, Premium	28
GL 45 laboratory bottle, pressure plus	25
GL 45 laboratory bottle, protect, PU-coated	24
GLS 80 laboratory bottle, amber	16
GLS 80 laboratory bottle, clear	15
GLS 80 laboratory bottle, protect, PU-coated	17
HPLC	26
reagent, narrow neck	41, 43
reagent, wide neck	40, 42
rolled flange	102
Roller	102
soda lime	37–38, 42–43
square	27, 37, 99
weighing	64
Woulff	148
Breed-Demeter, square bottle	99
Buechner funnel	122
Bulb pipette	115
Burette	111–112
Burette bottle	41, 43, 113

C

Calcium chloride cylinder	162
Candles, filter micro	129–130

Caps, see screw caps	
glass	99
Kapsenberg	98
metal	101
Centrifuge tube	71
Cone*	
Coplin type staining jar	103
Clamp	
closure for rolled bottle	102
quick release	147
Cleaning scraper	62
Clip, KECK™	164–165
Column, Vigreux	153
Condenser	
Allihn	154
bulb	154
coil distillate	155
Dimroth	155
jacketed coil	156
Liebig	154
West	154
Cone	
Connection pieces	162
Connection system	
GL 45	32–33
GLS 80	19
DG Safety Caps	34–36
DG Safety Waste Caps	35
Coupling, screwthread	166
Crucible, filter	123
Crystallizing dish	64
Culture bottle	97
Culture flask	
Erlenmeyer shape	100
Fernbach type	95
Kolle type	96
penicillin type	97
Roux type	96
Culture media bottle	99–101
Culture tube	72–74
Cylinder	
calcium chloride	162
measuring	110
mixing	109
multi purpose	68
standing	68

D

Desiccator	79
base	80–82
lid	83–84
O-ring	86
plate	85
screw cap	88
stopcock	86–87, 89
DG Safety Caps	34–36
DG Safety Waste Caps	35
Dimroth coil condenser	155
Dish	
crystallizing	64
evaporating	63

Petri	93–94
watch glass	65
Disposable culture tube	73–74
Drechsel head	128
Dropping bottle	47
Dropping funnel	156, 157
Drying tube, bent	161
DUROPLAN® Petri dish	93

E

Engler flask	61
Erlenmeyer flask	
for Kapsenberg caps	98
narrow neck	56
screwthread (DIN thread)	57
wide neck	57
with ground joint	58, 150
Evaporating dish	63
Evaporating flask	151

F

Feigl type spot plate	94
Fernbach type culture flask	95, 101
Filter apparatus, Witt type	136–138, 144
Filter candle, micro	129–130
Filter crucible / funnel adapter	123
Filter disk*	
with glass rim	125
Filter funnel	121–122, 125, 130
Filter head	124
Filter, micro	129–130
Filter, immersion	126, 129
Filter tube, Allihn type	131
Filtering apparatus	136–137
Filtering flask	131–135
Flask	
culture	95, 97, 100–101
Engler	61
Erlenmeyer	56–58, 98, 100, 150
evaporating	151
flat bottom	60, 151
flat flange	143
filtering	131–135
iodine determination	58
pear shape	149
penicillin	97
round bottom	59, 143, 149, 152–153
triple-neck round bottom	152–153
twin-neck round bottom	152
volumetric	109
Flat bottom flask	
narrow neck	60
wide neck	60
with ground joint	151
Flat flange*	
beaker	144
lid	145–146
quick release clamp	147
reaction vessel	143, 144
round bottomed flask	143
Funnel	

analytical	120
Buechner	122
Bunsen	120
dropping	156–157
filter	121–122, 125, 130
filter adapter	123
filter, micro	130
short stem	119
long stem	120
powder	119
separating	158

G

Gas distribution tube	126
Gas washing bottle	127–129
Glass box	104
Glass cap	99
Glass ceramic laboratory protection plate	61
Glass stopper, amber	44
GL 45, GL 32, GL 25, see bottles	
GLS 80	
bottle	15–17
connection system	19
membrane cap	18
pouring ring	17–18
quick release screw cap	17–18
stirred reactor	20–21
GUKO	135

H

Head, Drechsel	126, 129
Hellendahl type staining dish	103
Hexagon base measuring cylinder	110
Holding device for reaction vessel	147
Hose connection, plastic	134, 170
HPLC bottle	26

I

Immersion filter	126, 129
Immersion filter, micro	129
Iodine flask	58

J

Jacketed coil condenser	156
Jar	65–69, 103–104

K

Kapsenberg cap	98
KECK™	
assembly set	134
for conical joint*	164–165
for spherical joint*	165
Kolle type culture flask	96
KPG stirrer bearing	159–160
KPG stirrer shaft	160–161

L

Laboratory bottle, see bottles	
Laboratory protection plate, glass ceramic	61
Levelling bottle	46
Lid	

\* This article you will find in the DURAN® laboratory glassware for manipulators catalogue.

\* This article you will find in the DURAN® laboratory glassware for manipulators catalogue.

for desiccator	83–84
for filter apparatus, Witt type	136, 138
for flat flange vessel	145–146
Liebig condenser (West condenser)	154

M

Measuring cylinder*	110
Measuring pipette	113–114
Membrane filter	19, 26, 32, 34
Membrane screw cap	
GL 25, GL 32, GL45	30
GLS 80	18
Metal cap	101
Micro burette	112
Micro filter candle	129, 130
Micro immersion filter	129
Mixing cylinder	109, 110
Museum jar	67

N

Narrow neck reagent bottle	41, 43
NMR tube	74

O

Organ storage jar	65
O-ring	
for desiccator	86
for flat flange vessel	146–147

P

Pear shape flask	149
Penicillin flask	97
Petri dish	93–94
Philips beaker	55
Pipeline filter	126
Pipette	
bulb	115
for dropping bottle	47
measuring	113–114
Plastic hose connection	134, 170
Plastic stopper	44
Plate holder	62
Porcelain desiccator plate	85
Pouring ring	
GL 32,	30–31
GL 45	29–31
GLS 80	17–18
Powder funnel	119
Pregl, micro filter funnel	130
Premium bottle	28
Premium screw cap	29
Pressure equalization set	19, 26, 32
Pressure plus bottle	25
Protect bottle	17, 24
Pump, water jet	138–139

Q

Quick release clamp	147
Quick release closure, GLS 80	17–18

R

Reaction vessel, flat flange	143
Reagent bottle	
narrow neck	41, 43
wide neck	40, 42
Receiver adapter, bent	162, 163
Replacement seal	170
Reservoir bottle	112
RODAVISS®*	
Rolled flange bottle	102
Roller bottle	102
Round bottom flask	
narrow neck	59
triple-neck	153–154
twin-neck	152
wide neck	59
with ground joint	149, 152–153
Roux type culture flask	96–97
Rubber sleeve	123

S

Safety stopcock	89
Schiefferdecker type staining dish	104
Screw cap	
GLS 80 membrane screw cap	18
GLS 80 quick release screw cap	17, 18
membrane screw cap	18, 30
Premium screw cap	29
screw cap, PBT, PP	30–31
tamper evident screw cap	29
with aperture	128, 168
Screwthread coupling	166
Screwthread to ground joint adapters	167
Screwthread tube	166
Seal, silicone	49, 168–170
Security screw cap	39
Sedimentation cone	69
Separating funnel	157
Shaft, KPG stirrer	160, 161
Silicone rubber seal	168
Silicone seal (Septa)	169
Silicone sealing ring (VMQ)	86
Socket* for spirit lamps see spirit lamp	
Soda-lime	
disposable culture tubes	73–74
funnel	119–120
pipette	47, 113–115
reagent bottle, narrow neck	43
reagent bottle, wide neck	42
screw cap bottle, round	38
screw cap bottle, square	37
staining dish	103–104
test tube	70
Specimen jar	67
Spherical joints*	
Spirit lamp	63
Spot plate, Feigl	94
Square bottle	
Breed-Demeter	99
laboratory bottle, DURAN®	27
screw cap bottle, soda lime	37, 38

Square quadrupod	62
Staining jar, Coplin type	103
Staining tray	103, 104
Stainless-steel desiccator plate	85
Stainless steel handle	105
Standing cylinder	68
STERIPLAN® Petri dish	93
Stirred Reactor, GLS 80	20
Stirrer	
bearings, KPG	159–160
shafts, KPG	160–161
Stopcock	
single way*	
three way*	
two way*	
with PTFE-spindle for desiccator	86–87
with standard ground joint for aspirator	
bottle	48–49
Stopper	
for desiccator type WERTEX	80, 89
glass	44–45
plastic	44

T

Tamper-evident screw cap	29
Test tube	70
Threaded filter head	124
Triple-neck round bottomed flask	152–153
Tube	
centrifuge	71
culture	72–74
NMR	74
test	70
Tubing clamps, KECK™*	
Twin-neck round bottomed flask	152

V

Vacuum	
Desiccator	79–80
receiver adapter, bent	163
Valves	
GU*	
PRODURAN®*	
Vigreux column	153
Volumetric flask	109

W

Watch glass dish	65
Washing bottle, gas	127–129
Water jet pump	138–139
Weighing bottle	64
Wick for spirit lamp	63
Wide neck	
GLS 80 bottle	15–16
reagent bottle	40, 42
Witt, filtration apparatus	137
Woulff bottle	148

\* This article you will find in the DURAN® laboratory glassware for manipulators catalogue.

\* This article you will find in the DURAN® laboratory glassware for manipulators catalogue.

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GENERAL NOTE

The DURAN® laboratory glassware catalogue provides a basic information source for ordering our products. It does not represent a proposal for concluding a concrete agreement and will only serve as the basis for a contract upon explicit inclusion in a contractual relationship. We reserve the right to make changes to technical specifications, article numbers, packaging and design (e.g. due to changes of directives and standards). The contents of the catalogue have been created with the greatest possible care. However, we can accept no liability for the correctness, completeness and actuality of the contents. The presented replicated images provide an illustration of the article, details may however differ from the actual article.

REGISTERED TRADEMARKS

DURAN®, DUROPLAN®, FIOLAX®, KPG®, PRODURAN®, RODAVISS® and KECK™ are registered trademarks.

TERMS AND CONDITIONS FOR DELIVERY AND PAYMENT

The following terms and conditions for delivery and payment shall govern all deliveries and services. Any conflicting purchasing terms and conditions of the customer are hereby expressly rejected. Any such purchasing terms and conditions shall apply only if we expressly confirm them in writing. Acceptance of the delivered goods shall be deemed acknowledgement of our terms and conditions.

1. Prices/ Terms and Conditions of Payment

**1.1** We calculate the applicable prices on the date of delivery, such prices being in EURO (EUR) unless otherwise stipulated, plus an additional amount for VAT as applicable from time to time. Unless special terms are agreed, the prices should be understood to be prices ex works, with no deduction or discount being granted for immediate payment.

**1.2** If payment deadlines specified in our order confirmation are not met, this will automatically give rise to all of the statutory consequences of default, without any special reminder being required. In particular, we reserve the right to charge interest at the applicable rate charged by our bank for utilised credit if such interest exceeds the interest rate prescribed by statute (8 % above the reference interest rate). Furthermore, the entire balance shall become due and payable immediately, irrespective of any payment targets.

**1.3** In the case of custom-made products, we reserve the right to increase the price by a reasonable amount and to deviate to a reasonable extent from the agreed quantity. The customer must take delivery of surplus quantities. In respect of pressed pieces (Presslinge) and glass feedthrough headers (Glasdurchführungen), a variation of +/-10 % in relation to the ordered quantity shall be deemed as the agreed tolerance.

2. Delivery Dates

We will endeavour to adhere to stipulated delivery deadlines. However, due to the hazards and peculiar features of glass processing, delivery deadlines will not be binding unless expressly agreed otherwise. Our delivery times are subject to our supplier delivering the correct products to us on time.

3. Place of Performance and Passing of Risk

**3.1** The place of performance for the delivery is the

principal place of business of our respective supplying factory. The place of performance for payment is our principal place of business.

**3.2** When goods are transported, the risk (of accidental loss, destruction or deterioration) (the "Risk") shall pass to the customer as soon as we have delivered the goods to the carrier chosen by us.

4. Packaging

Unless otherwise agreed, we will accept the return of packaging to the extent that we are obliged to do so under the German Packaging Regulation (Verpackungsverordnung) or other legal requirements.

5. Warranties in Respect of Material Defects and Notification of Defects

**5.1** If, despite the greatest of care being taken, the goods give rise to complaints, then, in accordance with § 377 of the German Commercial Code (Handelsgesetzbuch, or "HGB"), obvious defects must be notified without delay, in any case no later than 14 calendar days after receipt of the goods, and hidden defects must be notified without delay after their discovery, otherwise the goods shall be deemed accepted. We shall not be liable for damage to deliveries through the breakage of glass during transit ("breakages") where the cause of the damage arises after the transfer of risk. In such cases any claim for breakages must be made against the carrier or the policy of transit insurance. We shall not pay compensation for breakages where the value of the relevant item is EUR 20.00 or less unless the breakage is due to intentional conduct or gross negligence on the part of ourselves or our servants or agents.

**5.2** Claims on the basis of defects as to quality ("material defects") shall become statute-barred 12 months after delivery of our goods to our customer. The foregoing provisions shall not apply to the

extent that longer limitation periods are prescribed by statute pursuant to § 438 (1) No. 2 of the German Civil Code (Bürgerliches Gesetzbuch, or "BGB" – Physical Structures and Physical Objects used for Physical Structures), § 479 (1) BGB (Recourse Claim), and § 634a (1) BGB (Construction Defects).

**5.3** Delivered goods which are returned to us ("returns") shall only be accepted if we are notified of the returns before their dispatch and the following conditions are satisfied:

5.3.1. The orderer shall receive a returns number when the return of items is notified to us; this returns number must appear on the documentation for the returned items;

5.3.2. Each return must be reported to our freight centre by carriage documents marked with the returns number.

**5.4** If, despite all care being taken at our premises, the delivered goods contain a defect that already existed at the time that the Risk passed, then we will, at our election and subject to receiving notification of the defect within the required time, repair the goods or deliver substitute goods. We must always be given the opportunity to render subsequent performance within a reasonable time.

**5.5** If subsequent performance cannot be rendered, the customer may – notwithstanding any claims for compensatory damages – rescind the agreement or reduce the contractual fee. The customer may not demand compensation for expenses incurred in vain.

**5.6** The following shall not give rise to any claims based on defects: merely immaterial deviations from the agreed condition of the goods, merely immaterial impairments to their utility, natural wear and tear, or loss or damage that arises after the Risk has passed as a result of incorrect or careless treatment, overuse, unsuitable operating resources, defective building work, unsuitable building foundations or special external influences that are not requirements under the contract. In addition, if the customer or a third party improperly (in a non-workmanlike manner) carries out maintenance work or modifications, then no claims based on defects may be made in respect of such work or

modifications or the resulting consequences.

**5.7** Claims on the part of the customer for expenses necessary to enable subsequent performance, particularly transport, infrastructure (e.g., tolls) and labour costs and the cost of materials, are excluded to the extent that such expenses are higher because the goods delivered by us were subsequently taken to a location other than the customer's business premises, unless such transportation is consistent with the authorised use of the goods.

**5.8** Any recourse claims on the part of the customer against us shall exist only to the extent that the customer has not entered into any agreements with its customers going beyond the mandatory statutory claims based on defects. Item 5.7 shall apply mutatis mutandis to the scope of the customer's recourse claim against the supplier.

**5.9** Claims based on material defects on the part of the customer against us or our vicarious agents that go beyond or are not included in the claims governed by Item 5 are excluded, irrespective of the stipulations of Item 7 concerning claims for compensation for damages and reimbursement of costs on the part of the customer.

**5.10** If, within the meaning of § 443 BGB, a defect is fraudulently concealed or a warranty is given with respect to the condition of the goods as at the time the Risk passes (seller's representation that the subject matter of the sale has a particular characteristic at the time the Risk passes and that the seller wishes to be held responsible for all consequences flowing from the fact that the characteristic is absent, regardless of fault), the customer's rights shall be exclusively governed by the statutory provisions.

**6. Industrial Property Rights and Copyright; Title Defects**

**6.1** Unless otherwise agreed, we have an obligation (although such obligation exists only in the country in which the place of delivery is located) to deliver the goods free from the industrial property rights and copyrights of third parties (hereinafter referred to as "Proprietary Rights"). In the event that a third party makes legitimate claims against the customer for infringement of Proprietary Rights based on the goods delivered by the supplier and used in accordance with the contract, we shall be liable to the customer within the period specified in Item 5.1 above as follows:

a) At our election and at our own expense, we will either secure a licence for the goods concerned, modify them so that the Proprietary Right is not infringed, or exchange them. If we are unable to do any of the above on reasonable terms, then the customer shall be entitled to the statutory rights of rescission and reduction of the contract price. The customer may not demand compensation for expenses incurred in vain.

b) The provisions of Item 7 shall apply to any claims for compensatory damages.

c) Our obligations as described above shall exist only on the condition that the customer notifies us in writing without delay of the claims asserted by the third

party, the customer does not admit to the infringement and leaves in our hands any defence of the claims and settlement negotiations. If the customer discontinues using the delivered goods in order to mitigate loss or for any other good reason, then the customer shall bring to the attention of the third party the fact that discontinuing use of the goods in no way constitutes an admission of an infringement of Proprietary Rights.

**6.2** Claims on the part of the customer are excluded if the customer is responsible for the infringement of Proprietary Rights.

**6.3** Claims on the part of the customer shall be further excluded if the infringement of Proprietary Rights is a result of special instructions issued by the customer, an application or use of the goods that was not foreseeable by us, or as a result of the customer modifying the goods or using them together with goods not delivered by us.

**6.4** In the event of an infringement of Proprietary Rights, the provisions set forth under Items 5.4 and 5.8 shall otherwise apply mutatis mutandis to the customer's claims governed by Item 6.1 a).

**6.5** If other title defects exist, then the provisions of Item 5 shall apply mutatis mutandis.

**6.6** Claims based on title defects on the part of the customer against us or our vicarious agents that go beyond or are not included in the claims governed by Item 6 are excluded.

**6.7** If, within the meaning of § 444 BGB, a defect is fraudulently concealed or a warranty is given with respect to the condition of the goods as at the time the Risk passes (seller's representation that the subject matter of the sale has a particular characteristic at the time the Risk passes and that the seller wishes to be held responsible for all consequences flowing from the fact that the characteristic is absent, regardless of fault), the customer's rights shall be exclusively governed by the statutory provisions.

**7. Other Claims for Compensatory Damages**

**7.1** In the event of a breach of a pre-contractual, contractual and/or non-contractual obligation, including unsatisfactory delivery, tortious conduct and manufacturer's liability, we shall be liable for compensatory damages and the reimbursement of costs – subject to further contractual or statutory liability requirements – only in the case of wilful conduct, gross negligence or breach of a material contractual duty ("condition") (contractual duty, the infringement of which jeopardises the ultimate purpose of the contract) where such breach was due to ordinary negligence. However, except in the case of wilful conduct, our liability shall be limited to typical contractual loss or damage that was foreseeable at the time the contract was entered into. The customer is not permitted to make a claim for expenses incurred in vain.

**7.2** For loss or damage caused by delay due to ordinary negligence, we shall be liable for up to only 5 % of the purchase price agreed with us.

**7.3** Except in the case of breach of a condition,

liability for ordinary negligence is excluded. This shall not affect Item 7.2.

**7.4** The exclusions and limitations of liability set forth under Items 7.1 to 7.3 shall not apply in the event that a warranty is given with respect to the condition of the goods within the meaning of § 443 BGB, a defect is fraudulently concealed, or in the event of injury to life, physical injury or injury to health, or strict liability under the German Product Liability Act (Produkthaftungsgesetz).

**8. Non-binding Nature of Drawings, Diagrams, Measurements and Weights**

Drawings, diagrams, measurements and weights are approximate only, unless they are expressly stipulated to be binding. The customer must guarantee that working drawings (construction diagrams) supplied by it do not infringe the Proprietary Rights of third parties. The customer must hold us harmless in the event that rights of recourse are asserted by third parties.

**9. Documents**

Documents supplied by us may not be copied or made available to third parties, or used for any purpose other than the agreed purpose.

**10. Reservation of Title**

**10.1** We shall retain title to the goods until all of our claims, including claims arising in the future, are fully paid. The customer may process and sell the goods in accordance with the following conditions: If the goods are further processed or remodelled by the customer, then we shall be deemed the manufacturer the within the meaning of § 950 BGB and shall acquire direct title to the intermediate or final products. The customer hereby relinquishes in our favour the ownership of the new goods created by further processing or remodeling on conclusion of the relevant purchase contract. The customer shall be merely the custodian of the goods. If the goods subject to the reservation of title ("reserved goods") are mixed or processed with other property not belonging to us, then we shall acquire a co-ownership interest in the new item proportionate to the value of the reserved goods to the other property.

**10.2** The goods may be sold only in the normal and ordinary course of business and only if claims deriving from their resale are not assigned to third parties beforehand. The customer's claims deriving from resale of the reserved goods shall, upon execution of the purchase agreement between the customer, and us be deemed assigned to us in advance to the extent that our goods are mixed or interprocessed with other property. In such a case, the assigned claims shall serve as our security only up to the value of the reserved goods sold in each case. We will not collect on the assigned claims for as long as the customer complies with its payment obligations. However, the customer has an obligation to disclose to us the identity of the third party debtor at our request and to notify such debtor of the assignment. The customer may collect on the claims resulting from sale of the reserved goods unless

and until it receives instructions to the contrary from us. The customer must immediately transfer the amounts collected by it to the extent that our claims are due.

**10.3** Pledges or the granting of security interests or relinquishment of the reserved goods or the assigned claims are not permitted. The customer must inform us immediately of any action by third parties affecting the goods delivered subject to a reservation of title or the assigned claims. We agree to release the assigned claims at our election if they exceed the value of our claims to be secured by more than 20 % and are derived from fully paid goods.

**10.4** In the event of breach of duty by the customer, particularly in the case of default on payment, we are entitled to rescind the agreement in whole or in part and recover the goods. The customer has an obligation to deliver up the goods. The declaration of recovery or reservation of title or assignment of the goods by us constitute a declaration of withdrawal from the agreement with respect to the reserved goods.

**10.5** If, in the case of non-domestic sales, the reservation of title agreed under Item 10 is not permitted with the same effect as under German law, then we shall retain title to the goods until payment of all of our claims arising out of the contractual relationship formed through the sale of the goods. If the foregoing reservation of title is not permitted with the same effect as under German law either, but it is permissible to reserve other rights in respect of the goods, then we are authorised to exercise all of these rights. The customer shall cooperate in all actions we may wish to take in order to protect our ownership interest or alternative right in the goods.

**11. Returns**

The following rules shall apply to returns unless the returns comprise defective goods (clause 5.3):

**11.1** The returned goods must have been purchased within 4 weeks in the case of deliveries within the Federal Republic of Germany or within 8 weeks in the case of deliveries to European customers or within 12 weeks in the case of deliveries to overseas customers.

**11.2** The provisions of clause 5.3 shall apply mutatis mutandis to the acceptance, notification and labelling of returns.

**11.3** Only unopened and undamaged goods without additional stickers or labelling shall be accepted as returns. We must be able to resell the goods.

**11.4** The returns shall be sent at the orderer's cost and risk.

**11.5** We shall also charge a handling fee equivalent to 20 % of the value of the item subject to a minimum charge of EUR 20.00.

**11.6** With the exception of paragraphs (2) and (3) the above provisions shall also apply to returned goods which the customer did not order.

**12. Relevant Information regarding Electronic Commerce**

In case we use electronic means for the conclusion of a contract with regard to the supply of goods or

rendering of services ("e-commerce") in the sense of Sec. 312e of the German Civil Code (Bürgerliches Gesetzbuch) the customer waives its rights regarding a) availability and explanation of a system to recognize and correct input errors before submitting an order and b) the supply of information regarding (i) the necessary steps to conclude a contract, (ii) electronic storage of the contract so concluded and its accessibility to the customer, (iii) the possible languages available for contract conclusion.

**13. Applicable Law and Judicial Forum**

**13.1** With the exception of choice of law rules under German private international law and the provisions of the UN Convention on Contracts for the International Sale of Goods ("UN-CISG"), the substantive law of the Federal Republic of Germany shall apply to all legal relationships and transactions established by this purchase agreement.

**13.2** Sole place of adjudication for both parties for all legal disputes arising out of the relevant purchase contracts or in connection with the supply relationship, including matters of currency exchange, is our head offices. If we appear as the plaintiff, we are also entitled to bring an action before the court responsible for the customer's head office.

**14. Moulds and tools**

Moulds and tools whether manufactured by us or sourced from third parties on behalf of the customer shall remain in our ownership and possession. At the start of the contract, the customer shall pay the agreed mould and tool contribution which grants the right to exclusively be supplied from these moulds. At the end of the contract, or any other discontinuation of the project, no assignment or transfer of the moulds and tools will take place; they will remain the property of, and in the possession of, us. In these cases, however, the customer shall be entitled to demand that we scrap the moulds and tools at our own expense and provide evidence of the scrapping to the customer. An obligation by us to store project-related moulds and tools shall end automatically at the end of the contract or project. If there is no written agreement to the contrary, a project shall be deemed to have ended after the expiry of a two-year period after the confirmation, by us, of the customer's last order. We shall ensure proper storage, handling and maintenance of the moulds and tools within the usual scope, during the term of the project. If the moulds or tools are destroyed or damaged due to improper storage, handling or maintenance by us then they shall be repaired or newly acquired at our expenses. The same applies to loss, destruction or damage as a result of force majeure. In the case that moulds and tools are used beyond their limit of wear and tear, the customer shall bear the costs of the new moulds and tools to be acquired by us, up to the amount of the originally agreed cost contribution for the worn part. Should the limit of wear and tear be reached prior to reaching an output quantity individually guaranteed, or

the customer proves that the wear and tear is due to a fault of the mould or tool, or an operating error by us, then we will bear the full cost of replacement. The above provisions shall apply accordingly to the moulds and tools acquired as replacement.

Agreements to the contrary must be in writing.

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“RESEARCH? ... WHAT IT MEANS TO ME:  
**REMAINING CURIOUS**, ALWAYS ASKING NEW  
QUESTIONS ... AND WHEN IT'S A QUESTION OF  
GLASSWARE, THEN I ALWAYS ASK FOR DURAN®.”

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WHY DURAN® IS SPECIAL

Optimum physical and chemical characteristics make DURAN® glass the material of choice for laboratory applications:

- Uniform wall thickness distribution throughout the product
- Therefore better mechanical stability and greater thermal shock resistance
- Prevents material stresses, which may cause the glass to crack when subject to heating or cooling
- Benefits: improves safety for laboratory technicians, increases the lifespan of laboratory glassware and protects valuable substances

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