



Scientific Equipment di Gnudi Andrea e Antonella S.r.l. Tel. 051501153 www.favs.it • info@favs.it

• Power supply for electrophoresis cell

- 1 | Easy to use the LCD display and control panel
- 2 | Constant voltage, current or power.
- 3 | Programed value and actual value are indicated at the same time
- 4 | Programmable time function
- 5 | Memory for 10 electrophoresis programs
- 6 4 parallel outputs to connect up to 4 cuvettes at the same time.
- 7 | With automatic detection of lack of load, overload, short circuit, rapid change of resistance and system overheating



Code	ZFD014	ZFD016
Output voltaje	10-300 V	10 - 600 V
Output current	1-400 m A	1 - 500mA
Output power	1-120 W	1 - 300 W
Resolution	1V, 1 mA , 1W	
Timer	1 min to 99 h 59 min	
Backlit	LCD display	
Outputs	4 parallel outputs	
Dimensions	230x237x118mm	
Power supply	220V :50Hz	